

**Fisherbrand**<sup>®</sup>  
QUALITY. RELIABILITY. VALUE.

# The Everyday Fisherbrand Product Guide



# EXPERIENCE THE DIFFERENCE



WITH FISHER SCIENTIFIC  
LAB ESSENTIALS

**Fisherbrand®**

*Quality You Demand, Value You Appreciate*

When making buying decisions, look for the Fisherbrand name—your assurance of quality, reliability, and value. Our wide selection of products is always in stock and readily available, making it easy for you to purchase the right product every time.

At Fisher Scientific, we stand wholeheartedly behind the Fisherbrand name as we proudly extend to you the Fisherbrand Promise: If our product does not perform to published specifications, we will replace it—no questions asked.

## FISHER SCIENTIFIC LAB EQUIPMENT

*Designed and manufactured with quality and performance in mind.*

accumet • accuSpin • FisherBiotech • Isotemp  
Maxima • Micromaster • PowerGen • Stereomaster

## CHEMICALS

*Thousands of high-purity reagents and fine chemicals for research and testing.*

Fisher Chemical • Fisher BioReagents

# Table of Contents

---

Aerosol Coolant.....	1	Microplates.....	20
Apparel.....	1	Microplate Sealing Supplies.....	21
Applicators.....	1	Microscope Coverslips.....	22
Bags.....	1	Microscope Slide Holders.....	22
Beakers.....	2	Microscope Slides.....	22
Bench Protectors.....	2	Petri Dishes.....	23
Bottles.....	3	Pipet Tips.....	23
Boxes.....	10	Pipets.....	26
Brushes.....	11	Pipettors.....	28
Buffers.....	11	Pipettor and Pipet Accessories.....	28
Carboys and Tanks.....	11	Racks.....	28
Carts.....	11	Refrigerator and Freezer Accessories.....	29
Cell Scrapers.....	12	Samplers.....	29
Cell Strainers.....	12	Scissors.....	29
Chromatography Supplies.....	12	Spatulas.....	30
Closures.....	12	Stir Bars.....	30
Counting Devices.....	13	Storage Microplates.....	31
Cuvets.....	14	Tape, Labels, and Labeling Tapes.....	31
Cylinders.....	14	Thermometers - Nonmercury.....	32
Detergents and Cleaners.....	14	Timers.....	32
Eye Protection.....	15	Tubes.....	33
Film and Foil Wrapping.....	15	Tubing.....	36
Filters.....	15	Vials.....	36
Forceps.....	17	Weighing Dishes.....	40
Funnels.....	18	Wipes.....	40
Gloves.....	18		
Inoculating Loops and Needles.....	19		
Insulated Containers.....	19		
Knives and Knife Blades.....	19		
Marking Pens.....	20		
Matting, Controlled Environment.....	20		





# Aerosol Coolant

**19-098-912** - Pk. of 100, Cs. of 5  
White; Diameter: 28 in. (71cm)

## Fisherbrand Super Friendly Freeze'It



**23-022-524** - Cs. of 6  
Slide Mounting Media; Fisherbrand Super Friendly Freeze'It; Aerosol; For cooling liquid, semisolid, and solid samples; 11 oz.; Noncorrosive; Nonflammable; Odorless; Completely volatile gas; 5 in. polyethylene snorkel

# Apparel

## Beard Covers

### Fisherbrand Polypropylene Beard Cover



**19-098-917** - Pk. of 100, Cs. of 5  
Beard Cover; Fisherbrand; Polypropylene; Type: Latex-free; Size: Regular  
**19-098-918** - Pk. of 100, Cs. of 5  
X-Large

## Boot and Shoe Covers

### Fisherbrand Conductive Shoe Covers



**19-098-925** - Pk. of 50, Cs. of 3  
Covers, Shoe; Fisherbrand; Polypropylene; Conductive carbon strip; Blue; Size: Regular  
**19-098-926** - Pk. of 50, Cs. of 3  
Size: X-Large

### Fisherbrand Heavyweight Polypropylene Shoe Covers



**19-098-923** - Pk. of 50, Cs. of 3  
Shoe Covers; Polypropylene; Fisherbrand; Type: With Nonskid Surface; Size: Regular  
**19-098-924** - Pk. of 50, Cs. of 3  
Nonskid Surface, X-Large  
**19-098-921** - Pk. of 50, Cs. of 3  
Without Nonskid Surface, Regular  
**19-098-922** - Pk. of 50, Cs. of 3  
Without Nonskid Surface, X-Large

## Caps Hoods Hairnets

### Fisherbrand Latex-Free Polypropylene Bouffant Caps



**19-098-913** - Pk. of 100, Cs. of 5  
Cap, Bouffant; Fisherbrand; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic nonlatex band stitch in with soft floss for comfortable and secure fit; 21 in.; White; Flat packed Double polybag  
**19-098-914** - Pk. of 100, Cs. of 5  
24 in.; White  
**19-098-915** - Pk. of 100, Cs. of 5  
Blue; Diameter: 21 in. (53cm)  
**19-098-916** - Pk. of 100, Cs. of 5  
Blue; Diameter: 24 in. (61cm)  
**19-098-910** - Pk. of 100, Cs. of 5  
White; Diameter: 18 in. (46cm)

# Facemasks

## Fisherbrand Facemasks



**18-960B** - Cs. of 500  
Facemasks, With Earloops; Fisherbrand; Thin middle layer of meltblown polypropylene filter media lies between spunbonded polypropylene inner and blue layers ultrasonically welded to reduce particle generation; Filtration efficiency is 99%  
**18-960A** - Cs. of 600  
With Earloops and Clip  
**18-960C** - Cs. of 300  
With Ties

## Lab Coats and Jackets

### Fisherbrand Men's Long Lab Coats



**01-369A** - Ea.  
Lab Coat, Fisherbrand; Men; Blue Poplin; Three Pockets; Single-Breasted Five-Button Front; Vented Back; Side Slits 43 in. (109cm) Long; Small Chest Size 34 in.  
**01-369B** - Ea.  
Chest Size: 36; Blue  
**01-369D** - Ea.  
Chest Size: 40; Blue

# Applicators

## Fisherbrand Sterile Plastic Applicators with Cotton Tips

**14-959-96A** - Cs. of 1000  
Applicator Stick, Fisherbrand; Plastic shaft (polystyrene); Sterile; Cotton tip; 6 in. (15cm); 1/10 in. dia.; One per pack  
**14-959-96C** - Cs. of 2000  
Two per pack

# Bags

## Autoclave Bags

### Fisherbrand Polyethylene Biohazard Autoclave Bags



**01-832A** - Cs. of 100  
Bag, Biohazard Autoclave; Fisherbrand; Polyethylene; Plain; Clear; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); Thickness: 4.0mil; Size: 24 x 36 in.  
**01-827A** - Pk. of 200  
Clear polyethylene; Printed; Size: 8 x 12 in.; 1.5mil  
**01-827B** - Pk. of 200  
Clear polyethylene; Printed; Size: 12 x 24 in.; 1.5mil  
**01-827C** - Pk. of 200  
Clear polyethylene; Printed; Size: 24 x 30 in.; 1.5mil  
**01-827D** - Pk. of 200  
Clear polyethylene; Printed; Size: 24 x 36 in.; 1.5mil  
**01-814-1** - Pk. of 200  
Clear; Plain; Polyethylene; Thickness: 1.5mil; Size: 12 W x 24 in. L  
**01-814-2** - Pk. of 200  
Clear; Plain; Polyethylene; Thickness: 1.5mil; Size: 24 W x 30 in. L  
**01-814-3** - Pk. of 200  
Clear; Plain; Polyethylene; Thickness: 1.5mil; Size: 24 W x 36 in. L

- 01-829A** - Pk. of 200  
Red polyethylene; Printed; Size: 8 x 12 in.; 1.5mil
- 01-829B** - Pk. of 200  
Red polyethylene; Printed; Size: 14 x 19 in.; 1.5mil
- 01-829C** - Pk. of 200  
Red polyethylene; Printed; Size: 19 x 23 in.; 1.5mil
- 01-831A** - Cs. of 100  
Red polyethylene; Printed; Size: 24 x 36 in.; 4mil
- 01-829D** - Pk. of 200  
Red polyethylene; Printed; Size: 25 x 35 in.; 1.5mil
- 01-829F** - Pk. of 200  
Red polyethylene; Printed; Size: 33 x 30 in.; 1.5mil
- 01-829E** - Pk. of 100  
Red polyethylene; Printed; Size: 38 x 48 in.; 1.5mil
- 01-831B** - Cs. of 50  
Red polyethylene; Printed; Size: 38 x 48 in.; 4mil

**Fisherbrand Polyethylene Biohazard Autoclave Bags w/out Indicators**

- 22-044-560** - Cs. of 100  
Bags, Biohazard; Autoclave; Fisherbrand; Polyethylene; Without indicators; Clear; 1.5mil thick; Size: 24 x 36 in.

**Fisherbrand Polypropylene Biohazard Autoclave Bags**



- 01-826-1** - Pk. of 200  
Bag, Biohazard Autoclave; Fisherbrand; Polypropylene; Plain; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 8 x 12 in.
- 01-826-2** - Pk. of 200  
Clear polypropylene; Plain; Size: 10 x 10 in.; 2mil
- 01-826-3** - Pk. of 200  
Clear polypropylene; Plain; Size: 10 x 12 in.; 2mil
- 01-826-4** - Pk. of 200  
Clear polypropylene; Plain; Size: 10 x 15 in.; 2mil
- 01-826-5** - Pk. of 200  
Clear polypropylene; Plain; Size: 12 x 24 in.; 2mil
- 01-826-6** - Pk. of 200  
Clear polypropylene; Plain; Size: 24 x 30 in.; 2mil
- 01-826-7** - Pk. of 200  
Clear polypropylene; Plain; Size: 24 x 36 in.; 2mil
- 01-826A** - Pk. of 200  
Clear polypropylene; Printed; Size: 12 x 24 in.; 2mil
- 01-826B** - Pk. of 200  
Clear polypropylene; Printed; Size: 24 x 30 in.; 2mil
- 01-826C** - Pk. of 200  
Clear polypropylene; Printed; Size: 24 x 36 in.; 2mil
- 01-826D** - Pk. of 200  
Clear polypropylene; Printed; Size: 31 x 38 in.; 2mil
- 01-826E** - Pk. of 100  
Clear polypropylene; Printed; Size: 37 x 48 in.; 2mil
- 01-830A** - Pk. of 200  
Red polypropylene; Plain; Size: 8 x 12 in.; 2mil
- 01-830B** - Pk. of 200  
Red polypropylene; Plain; Size: 14 x 19 in.; 2mil
- 01-830C** - Pk. of 200  
Red polypropylene; Plain; Size: 19 x 23 in.; 2mil
- 01-830D** - Pk. of 200  
Red polypropylene; Plain; Size: 25 x 35 in.; 2mil
- 01-830E** - Pk. of 100  
Red polypropylene; Plain; Size: 37 x 48 in.; 2mil
- 01-828B** - Pk. of 200  
Red polypropylene; Printed; Size: 14 x 19 in.; 2mil
- 01-828C** - Pk. of 200  
Red polypropylene; Printed; Size: 19 x 23 in.; 2mil
- 01-828D** - Pk. of 200  
Red polypropylene; Printed; Size: 25 x 35 in.; 2mil
- 01-828E** - Pk. of 100  
Red polypropylene; Printed; Size: 37 x 48 in.; 2mil
- 01-828A** - Pk. of 200  
Red polypropylene; Printed; Size: 8 x 12 in.; 2mil

## Sample Bags

**Fisherbrand Zipper Seal Sample Bags**



- 01-816-1E** - Cs. of 1000  
Bags, Sample; Fisherbrand; Zipper seal; Transparent polyethylene; Clear, print-free surface; 2mil thick; W x L: 9 x 12 in. (23 x 30cm)
- 01-816D** - Pk. of 50, Cs. of 10  
4mil thick; 5 W x 8 in. L
- 01-816E** - Pk. of 50, Cs. of 5  
4mil thick; 9 W x 13 in. L
- 01-816-1A** - Cs. of 1000  
Bag Zip Seal Pe 2mil 2x3 in. M/cs
- 01-816-1B** - Cs. of 1000  
Bag Zip Seal Pe 2mil 3x5 in. M/cs
- 01-816-1C** - Cs. of 1000  
Bag Zip Seal Pe 2mil 4x6 in. M/cs
- 01-816-1D** - Cs. of 1000  
Bag Zip Seal Pe 2mil 5x8 in. M/cs
- 01-816A** - Pk. of 50, Cs. of 10  
Bg Zip Seal Pe 4mil 4x6 in. 50/pk
- 01-816B** - Pk. of 50, Cs. of 5  
Bg Zipseal Pe 4mil 6x13 in. 50/pk
- 01-816C** - Pk. of 50, Cs. of 5  
Bg Zipseal Pe 4mil 9x18 in. 50/pk
- 01-816F** - Pk. of 50, Cs. of 4  
Bg Zipseal Pe 4mil 9x26 in. 50/pk

## Beakers

**Fisherbrand Polypropylene Beakers**

- S30688A** - Pk. of 1, Cs. of 4  
Beaker, Polypropylene; Autoclavable; Chemical resistant; W/pouring spout and permanent graduation scales; Capacity: 2000mL

## Bench Protectors

**Fisherbrand Absorbent Surface Liners**



- 14-127-46** - Cs. of 6  
Bench Liner; Fisherbrand; Absorbent multilayers; Waterproof polypropylene backing; Size: 20 in. x 50 ft.; Dispenser box
- 14-127-47** - Cs. of 1  
Size: 20 in. x 300 ft. roll

**Fisherbrand Absorbent Underpads**



- 14-206-62** - Pk. of 50, Cs. of 5  
Underpad, Absorbent; Fisherbrand; Cushioned surface; Chemical-resistant polyethylene backing; Spunbond top layer; Absorb up to 750mL/m<sup>2</sup>; Mat; L x W: 20 x 24 in. (51 x 91cm)
- 14-206-63** - Pk. of 25, Cs. of 5  
Mat; L x W: 20 x 36 in. (51 x 91cm)
- 14-206-64** - Ea., Cs. of 4  
Roll; 20 in. x 100 ft.; Dispenser pack
- 14-206-65** - Cs. of 1  
Roll; 20 in. x 250 ft.; Dispenser pack

**Fisherbrand Polyethylene Foam Surface Liners**



- 14-127-48** - Cs. of 2  
Bench Protector; Fisherbrand; Polyethylene foam; For surface lining; Chemical resistant; Non-absorbing, non-wetting; Thickness: 1/8 in.; Size: 12 in.W x 50 ft.L roll
- 14-127-49** - Cs. of 2  
Thickness: 1/8 in.; Size: 18 in.W x 50 ft.L roll
- 14-127-50** - Cs. of 2  
Thickness: 1/8 in.; Size: 30 in.W x 50 ft.L roll
- 14-127-52** - Cs. of 24  
Thickness: 1/4 in.; Size: 10-1/2 x 12 in. sheets
- 14-127-53** - Cs. of 12  
Thickness: 1/4 in.; Size: 16 x 21 in. sheets

# Bottles

## Glass Jars

**Fisherbrand Amber Straight-Sided Jars**



- 02-912-298** - Cs. of 24  
Jar, Straight-Sided; Fisherbrand; Amber glass; 4 oz. (125mL); Screw cap size: 58mm-400
- 02-912-045** - Cs. of 24  
4 oz.; With closure
- 02-911-859** - Cs. of 12  
8 oz.; With closure

**Fisherbrand Certified Clean Clear Glass Straight-Sided Jars**



- 02-912-309** - Cs. of 24  
Jar; Fisherbrand; Clear Glass; Straight-Sided; Certified clean; With attached Faced PE-lined closures; 2 oz. (60mL), 53mm-400 screw cap size, Comes with Certificate of Analysis
- 02-912-310** - Cs. of 24  
4 oz.
- 02-912-311** - Cs. of 24  
8 oz.
- 02-912-312** - Cs. of 12  
17 oz.
- 02-912-313** - Cs. of 12  
32 oz.

**Fisherbrand Clear Glass Testing Jars**



- 02-912-116** - Cs. of 180  
Jars and jugs; Fisherbrand; Jars, clear glass testing; Jar only; Capacity: 2 oz. (60mL); 43mm-400 screw cap size
- 02-911-883** - Cs. of 180  
With polyvinyl-lined cap; 2 oz. (60mL)
- 02-911-964** - Cs. of 180  
With Teflon Faced PE-lined caps; 2 oz. (60mL)

**Fisherbrand Clear Straight-Sided Jars with Polyvinyl-lined Closures**



- 02-911-920** - Cs. of 48  
Jars, Straight-Sided clear glass; Fisherbrand; with attached Polyvinyl-lined closures, 1 oz. (30mL), 43mm-400 screw cap size
- 02-911-773** - Cs. of 24  
2 oz.
- 02-911-775** - Cs. of 24  
4 oz.
- 02-911-776** - Cs. of 12  
8 oz.
- 02-911-772** - Cs. of 12  
16 oz.

- 02-911-774** - Cs. of 12  
32 oz.

**Fisherbrand Clear Straight-Sided Jars with PTFE Faced PE-lined Closures**




- 02-911-912** - Cs. of 48  
Jar; Clear; Straight-Sided; Fisherbrand; with attached PTFE Faced PE-lined closures, 1 oz. (30mL), 43mm-400 screw cap size; Cs. of 48
- 02-912-306** - Cs. of 384  
1 oz. (30mL); Cs. of 384
- 02-911-763** - Cs. of 24  
2 oz. (60mL); Cs. of 24
- 02-911-791** - Cs. of 144  
2 oz. (60mL); Cs. of 144
- 02-911-764** - Cs. of 24  
4 oz. (125mL)
- 02-911-765** - Cs. of 12  
8 oz. (250mL)
- 02-911-761** - Cs. of 12  
16 oz. (500mL)
- 02-911-762** - Cs. of 12  
32 oz. (1000mL)

**Fisherbrand Clear Straight-Sided Jars with Rubber-lined Closures**



- 02-912-345** - Cs. of 24  
Jar; Clear; Straight-Sided; Fisherbrand; with attached Rubber-lined closures, 4 oz. (125mL), 58mm-400 screw cap size

**Fisherbrand Clear Straight-Sided Jars without Closures**



- 02-912-304** - Cs. of 384  
Jar; Clear; Straight-Sided; Fisherbrand; without closures; 1 oz. (30mL); Screw cap size: 43mm-400
- 02-911-996** - Cs. of 144  
2 oz. (60mL)
- 02-911-997** - Cs. of 24  
4 oz. (125mL)
- 02-911-876** - Cs. of 12  
8 oz. (250mL)
- 02-912-118** - Cs. of 12  
16 oz. (500mL)
- 02-912-305** - Cs. of 12  
32 oz. (1000mL)

**Fisherbrand Safety-Coated Clear Wide-Mouth Jars**



- 02-911-825** - Cs. of 12  
Jar, clear glass; Fisherbrand; Wide mouth; Safety coated, With attached Polyvinyl-lined closures, Capacity: 8 oz. (250mL), 70mm-400 screw cap size, Wide mouth
- 02-911-824** - Cs. of 12  
With polyvinyl-lined cap; 16 oz.
- 02-911-867** - Cs. of 12  
With polyvinyl-lined cap; 32 oz.
- 02-911-827** - Cs. of 12  
With PTFE-faced PE-lined cap; 8 oz.
- 02-911-868** - Cs. of 12  
With PTFE-faced PE-lined cap; 16 oz.
- 02-911-721** - Cs. of 12  
With PTFE-faced PE-lined cap; 32 oz.



## Glass Jugs

### Fisherbrand Amber Glass Jugs



- 02-912-271** - Cs. of 6  
Jug; Amber Glass; Fisherbrand; With attached Teflon Faced PE-lined closures, 80 oz. (2365mL); 38mm-400 screw cap size
- 02-911-862** - Cs. of 4  
130 oz. (3840mL); with Poly-Seal-lined closures
- 02-912-277** - Cs. of 4  
130 oz. (3840mL); with Polyvinyl cap lining
- 02-912-270** - Cs. of 4  
130 oz. (3840mL); with Rubber-lined closures
- 02-912-276** - Cs. of 4  
130 oz. (3840mL); with Teflon Faced PE-lined closures

### Fisherbrand Clear Glass Jugs



- 02-912-307** - Cs. of 6  
Jug, Clear, Glass, Fisherbrand; with attached Poly-Seal-lined closures, 64 oz. (2000mL), 38mm-400 screw cap size
- 02-912-279** - Cs. of 6  
64 oz. (2000mL); with polyvinyl closures
- 02-912-275** - Cs. of 6  
64 oz. (2000mL); with Rubber-lined closures
- 02-912-280** - Cs. of 6  
64 oz. (2000mL); with Teflon Faced PE-lined closures
- 02-912-272** - Cs. of 4  
130 oz. (3840mL); with Poly-Seal-lined closures
- 02-912-273** - Cs. of 4  
130 oz. (3840mL); with polyvinyl closures
- 02-912-274** - Cs. of 4  
130 oz. (3840mL); with Rubber-lined closures
- 02-912-278** - Cs. of 4  
130 oz. (3840mL); with Teflon Faced PE-lined closures

## Petroleum Bottles

### Fisherbrand Clear Oil Sample Bottle



- 02-911-472** - Cs. of 144  
Bottle, Sample; Fisherbrand; Clear oil; Clear glass; Capacity: 4 oz. (125mL); With cork plug; Case of 144

## Plastic Wash Bottles

### Fisherbrand Labeled Wash Bottles



- 03-409-23Q** - Pk. of 6  
Bottle, Wash; Fisherbrand; Polyethylene; Wide Mouth; Color-coded label; 500mL; for Acetone
- 03-409-23A** - Pk. of 6  
Acetone; Capacity: 500mL; Pack of 6
- 03-409-23L** - Pk. of 6  
Bleach; Capacity: 500mL
- 03-409-23M** - Pk. of 6  
Deionized Water; Capacity: 500mL
- 03-409-23K** - Pk. of 6  
Distilled Water; Capacity: 500mL
- 03-409-23B** - Pk. of 6  
Ethanol; Capacity: 500mL
- 03-409-23S** - Pk. of 6  
Ethyl Ketone, 500mL capacity
- 03-409-23C** - Pk. of 6  
Isopropanol; Capacity: 500mL
- 03-409-23R** - Pk. of 6  
Machine Oil; Capacity: 500mL
- 03-409-23D** - Pk. of 6  
Methanol; Capacity: 500mL
- 03-409-23J** - Pk. of 6  
Saline; Capacity: 500mL
- 03-409-23P** - Pk. of 6  
Soap; Capacity: 500mL
- 03-409-23E** - Pk. of 6  
Toluene; Capacity: 500mL
- 03-409-23N** - Pk. of 6  
Universal; Capacity: 500mL
- 03-409-23H** - Pk. of 6  
Variety Pack; Capacity: 500mL
- 03-409-23G** - Pk. of 6  
Water; Capacity: 500mL

## Round Glass Narrow Mouth Bottles

### Fisherbrand Amber Boston Rounds with Poly-Seal-lined Closures



- 02-911-861** - Cs. of 720  
Bottle; Amber Boston Round; Fisherbrand; With attached Poly-Seal-lined closures; Narrow mouth; 1/2 oz. (15mL), 18mm-400 screw cap size; Cs. of 720
- 02-912-339** - Cs. of 48  
1 oz. (30mL); Cs. of 48
- 02-911-734** - Cs. of 24  
2 oz. (60mL); Cs. of 24
- 02-911-898** - Cs. of 288  
2 oz. (60mL); Cs. of 288
- 02-911-735** - Cs. of 24  
4 oz. (125mL); Cs. of 24
- 02-911-895** - Cs. of 160  
4 oz. (125mL); Cs. of 160
- 02-911-736** - Cs. of 12  
8 oz. (250mL); Cs. of 12
- 02-911-900** - Cs. of 108  
8 oz. (250mL); Cs. of 108
- 02-911-743** - Cs. of 12  
16 oz. (500mL); Cs. of 12
- 02-911-878** - Cs. of 60  
16 oz. (500mL); Cs. of 60
- 02-911-744** - Cs. of 12  
32 oz. (1000mL); Cs. of 12
- 02-911-899** - Cs. of 30  
32 oz. (1000mL); Cs. of 30



### Fisherbrand Amber Boston Rounds with Polyvinyl-lined Closures



**02-911-793** - Cs. of 720  
Bottles, Amber Boston Rounds; Fisherbrand; Narrow mouth; With attached Polyvinyl-lined closures; 1/2 oz. (15mL); 18mm-400 screw cap size; Cs. of 720

**02-912-337** - Cs. of 48  
1 oz. (30mL); Cs. of 48

**02-911-746** - Cs. of 24  
2 oz. (60mL); Cs. of 24

**02-911-742** - Cs. of 288  
2 oz. (60mL); Cs. of 288

**02-911-748** - Cs. of 24  
4 oz. (125mL); Cs. of 24

**02-911-723** - Cs. of 160  
4 oz. (125mL); Cs. of 160

**02-911-749** - Cs. of 12  
8 oz. (250mL); Cs. of 12

**02-911-724** - Cs. of 108  
8 oz. (250mL); Cs. of 108

**02-911-745** - Cs. of 12  
16 oz. (500mL); Cs. of 12

**02-911-725** - Cs. of 60  
16 oz. (500mL); Cs. of 60

**02-911-747** - Cs. of 12  
32 oz. (1000mL); Cs. of 12

**02-911-729** - Cs. of 30  
32 oz. (1000mL); Cs. of 30

### Fisherbrand Amber Boston Rounds with PTFE-Faced PE-lined Closures



**02-912-338** - Cs. of 48  
Bottles; Boston Rounds; Fisherbrand; Amber; With attached PTFE Faced PE-lined closures; 1 oz. (30mL), 20mm-400 screw cap size, Narrow mouth

**02-911-737** - Cs. of 24  
2 oz. (60mL)

**02-911-740** - Cs. of 24  
4 oz. (125mL)

**02-911-741** - Cs. of 12  
8 oz. (250mL)

**02-911-738** - Cs. of 12  
16 oz. (500mL)

**02-911-739** - Cs. of 12  
32 oz. (1000mL)

### Fisherbrand Amber Boston Rounds with Rubber-lined Closures



**02-912-340** - Cs. of 24  
Bottles, Amber Boston Rounds; Fisherbrand; Narrow mouth; With attached rubber-lined closures; 2 oz. (60mL), 20mm-400 screw cap size

**02-912-341** - Cs. of 24  
4 oz. (125mL)

**02-912-342** - Cs. of 12  
8 oz. (250mL)

**02-912-343** - Cs. of 12  
16 oz. (500mL)

**02-912-344** - Cs. of 12  
32 oz. (1000mL)

### Fisherbrand Amber Boston Rounds without Closures



**02-911-794** - Cs. of 720  
Bottle; Amber Boston Round, Fisherbrand; Narrow mouth; Without closures; 0.5 oz. (15mL); 18mm-400 screw cap size

**02-911-908** - Cs. of 288  
2 oz. (60mL)

**02-911-904** - Cs. of 160  
4 oz. (125mL)

**02-911-905** - Cs. of 108  
8 oz. (250mL)

**02-911-902** - Cs. of 60  
16 oz. (500mL)

**02-911-903** - Cs. of 30  
32 oz. (1000mL)

### Fisherbrand Certified Clean Amber Boston Rounds



**02-911-928** - Cs. of 12  
Bottles; Fisherbrand; Amber boston round; Certified clean; With attached faced PE-lined closures; 8 oz. (250mL); 24mm-400 screw cap size; Comes with certificate of analysis

**02-911-833** - Cs. of 12  
34 oz. (1000mL)

### Fisherbrand Certified Low-Particulate Clear Boston Rounds



**02-911-889** - Cs. of 24  
Bottles; Fisherbrand; Clear boston rounds; Certified low-particulate; With attached Polyvinyl-lined closures; Narrow-mouth capacity: 4 oz. (125mL); 22mm-400 screw cap size

**02-911-892** - Cs. of 12  
8 oz. (250mL)

**02-911-879** - Cs. of 12  
16 oz. (500mL)

**02-911-914** - Cs. of 12  
32 oz. (1000mL)

### Fisherbrand Clear Boston Rounds with Poly-Seal-lined Caps



**02-911-782** - Cs. of 48  
Bottle; Clear Boston Round; Fisherbrand; with attached Poly-Seal-lined closures, 1 oz. (30mL), 20mm-400 screw cap size, Narrow mouth

**02-911-783** - Cs. of 24  
2 oz. (60mL)

**02-911-784** - Cs. of 24  
4 oz. (125mL)

**02-911-785** - Cs. of 12  
8 oz. (250mL)

**02-911-750** - Cs. of 12  
16 oz. (500mL)

**02-911-751** - Cs. of 12  
32 oz. (1000mL)

### Fisherbrand Clear Boston Rounds with Polyvinyl-lined Caps



**02-911-779** - Cs. of 48  
Bottle; Clear Boston Round; Fisherbrand; with attached Polyvinyl-lined closures, 1 oz. (30mL), 20mm-400 screw cap size, Narrow mouth; Cs. of 48

**02-911-752** - Cs. of 24  
2 oz. (60mL); Cs. of 24

**02-911-790** - Cs. of 288  
2 oz. (60mL); Cs. of 288

**02-911-780** - Cs. of 24  
4 oz. (125mL); Cs. of 24

**02-911-730** - Cs. of 160  
4 oz. (125mL); Cs. of 160

**02-911-781** - Cs. of 12  
8 oz. (250mL); Cs. of 12

**02-911-732** - Cs. of 108  
8 oz. (250mL); Cs. of 108

**02-911-777** - Cs. of 12  
16 oz. (500mL); Cs. of 12

**02-911-726** - Cs. of 60  
16 oz. (500mL); Cs. of 60

**02-911-778** - Cs. of 12  
32 oz. (1000mL); Cs. of 12

### Fisherbrand Clear Boston Rounds with Rubber-lined Caps



**02-911-756** - Cs. of 48  
Bottles; Fisherbrand; Clear boston round; With attached rubber-lined closures; Type III soda lime glass; Narrow-mouth; Clear for easy viewing of contents; capacity: 1 oz. (30mL); 48/case; 20mm-400 screw cap size

**02-911-757** - Cs. of 24  
2 oz. (60mL)

**02-911-759** - Cs. of 24  
4 oz. (125mL)

**02-911-760** - Cs. of 12  
8 oz. (250mL)

**02-911-755** - Cs. of 12  
16 oz. (500mL)

**02-911-758** - Cs. of 12  
32 oz. (1000mL)

**Fisherbrand Clear Boston Rounds without Caps**



- 02-911-940** - Cs. of 432  
Bottles; Fisherbrand; Clear boston round; Without closures; Type III soda lime glass; Narrow-mouth; Clear for easy viewing of contents; capacity: 1 oz. (30mL); 432/ case; 20mm-400 screw cap size
- 02-911-944** - Cs. of 288  
2 oz. (60mL)
- 02-911-945** - Cs. of 160  
4 oz. (125mL)
- 02-911-872** - Cs. of 108  
8 oz. (250mL)
- 02-911-939** - Cs. of 60  
16 oz. (500mL)
- 02-911-927** - Cs. of 12  
32 oz. (1000mL)

**Fisherbrand Custom Cleaned Amber Boston Rounds**



- 02-911-999** - Cs. of 24  
Bottles; Fisherbrand; Amber boston round; Custom clean; With attached PTFE faced PE-lined closures; capacity: 2 oz. (60mL); 20mm-400 screw cap size; Noncertified
- 02-912-001** - Cs. of 24  
4 oz. (125mL)
- 02-912-002** - Cs. of 12  
8 oz. (250mL)
- 02-911-998** - Cs. of 12  
16 oz. (500mL)
- 02-912-000** - Cs. of 12  
32 oz. (1000mL)

**Fisherbrand Custom Cleaned Clear Boston Rounds**




- 02-912-009** - Cs. of 48  
Bottles; Fisherbrand; Clear boston rounds; Custom clean; With attached PTFE Faced PE-lined closures; capacity: 1 oz. (30mL); 20mm-400 screw cap size; Noncertified
- 02-912-010** - Cs. of 24  
2 oz. (60mL)
- 02-912-011** - Cs. of 24  
4 oz. (125mL)
- 02-912-012** - Cs. of 12  
8 oz. (250mL)
- 02-912-013** - Cs. of 12  
16 oz. (500mL)
- 02-912-014** - Cs. of 12  
32 oz. (1000mL)

**Fisherbrand Safety-Coated Amber Boston Rounds**



- 02-911-896** - Cs. of 24  
Bottles; Fisherbrand; Amber boston round; Safety coated; With attached Poly-seal-lined closures; Narrow-mouth; capacity: 4 oz. (125mL); 22mm-400 screw cap size
- 02-911-970** - Cs. of 12  
With poly-seal-lined closures; Capacity: 8 oz. (250mL)
- 02-911-969** - Cs. of 12  
With poly-seal-lined closures; Capacity: 16 oz. (500mL)
- 02-911-901** - Cs. of 12  
With poly-seal-lined closures; Capacity: 32 oz. (1000mL)
- 02-911-897** - Cs. of 24  
With polyvinyl-lined closures; Capacity: 4 oz. (125mL)
- 02-911-813** - Cs. of 12  
With polyvinyl-lined closures; Capacity: 8 oz. (250mL)
- 02-911-812** - Cs. of 12  
With polyvinyl-lined closures; Capacity: 16 oz. (500mL)
- 02-911-816** - Cs. of 12  
With polyvinyl-lined closures; Capacity: 32 oz. (1000mL)
- 02-911-894** - Cs. of 24  
With PTFE-lined closures; Capacity: 4 oz. (125mL)
- 02-911-815** - Cs. of 12  
With PTFE-lined closures; Capacity: 8 oz. (250mL)
- 02-911-817** - Cs. of 12  
With PTFE-lined closures; Capacity: 16 oz. (500mL)
- 02-911-814** - Cs. of 12  
With PTFE-lined closures; Capacity: 32 oz. (1000mL)

**Fisherbrand Safety-Coated Clear Boston Rounds**



- 02-911-860** - Cs. of 108  
Bottles; Fisherbrand; Clear boston round; Safety coated; Without closures; Shatter-resistant; Non-slip coating; Not autoclavable; capacity: 8 oz. (250mL); 24mm-400 screw cap size
- 02-911-921** - Cs. of 60  
Bottles only; Capacity: 16 oz. (500mL)
- 02-911-922** - Cs. of 30  
Bottles only; Capacity: 32 oz. (1000mL)
- 02-911-973** - Cs. of 12  
With poly-seal-lined caps; Capacity: 8 oz. (250mL)
- 02-911-971** - Cs. of 12  
With poly-seal-lined caps; Capacity: 16 oz. (500mL)
- 02-911-972** - Cs. of 12  
With poly-seal-lined caps; Capacity: 32 oz. (1000mL)
- 02-911-818** - Cs. of 12  
With polyvinyl lined caps; Capacity: 8 oz. (250mL)
- 02-911-819** - Cs. of 12  
With polyvinyl lined caps; Capacity: 16 oz. (500mL)
- 02-911-820** - Cs. of 12  
With polyvinyl lined caps; Capacity: 32 oz. (1000mL)
- 02-911-821** - Cs. of 12  
With PTFE-faced, PE-lined caps; Capacity: 8 oz. (250mL)
- 02-911-822** - Cs. of 12  
With PTFE-faced, PE-lined caps; Capacity: 16 oz. (500mL)
- 02-911-823** - Cs. of 12  
With PTFE-faced, PE-lined caps; Capacity: 32 oz. (1000mL)
- 02-911-910** - Cs. of 12  
With Rubber-lined caps; Capacity: 8 oz. (250mL)
- 02-912-318** - Cs. of 12  
With rubber lined caps; Capacity: 16 oz. (500mL)

# Round Glass Wide Mouth Bottles

**Fisherbrand Amber Wide-Mouth Packers with Poly-Seal-lined Closures**



- 02-911-444** - Cs. of 24  
Bottle; Amber Wide-Mouth Packer; Fisherbrand; with attached Polyvinyl-lined closures, 1 oz. (30mL); 28mm-400 screw cap size; Cs. of 24
- 02-911-445** - Cs. of 24  
2 oz. (60mL); Cs. of 24
- 02-911-446** - Cs. of 24  
4 oz. (125mL); Cs. of 24

**Fisherbrand Amber Wide-Mouth Packers with Polyvinyl-lined Closures**



- 02-911-852** - Cs. of 720  
Bottle; Amber Wide-Mouth Packer; Fisherbrand; with attached Polyvinyl-lined closures, 0.20 oz. (6mL), 20mm-400 screw cap size; Cs. of 720
- 02-911-857** - Cs. of 912  
0.27 oz. (8mL); Cs. of 912
- 02-911-799** - Cs. of 624  
0.5 oz. (15mL); Cs. of 624
- 02-911-432** - Cs. of 24  
1 oz. (30mL); Cs. of 24
- 02-911-447** - Cs. of 432  
1 oz. (30mL); Cs. of 432
- 02-911-433** - Cs. of 24  
2 oz. (60mL); Cs. Of 24
- 02-911-448** - Cs. of 216  
2 oz. (60mL); Cs. of 216
- 02-911-434** - Cs. of 24  
4 oz. (125mL); Cs. of 24
- 02-911-449** - Cs. of 180  
4 oz. (125mL); Cs. of 180
- 02-911-435** - Cs. of 24  
8 oz. (250mL); Cs. of 24

**02-911-450** - Cs. of 84  
8 oz. (250mL); Cs. of 84  
**02-911-436** - Cs. of 12  
16 oz. (500mL); Cs. of 12  
**02-911-451** - Cs. of 60  
16 oz. (500mL); Cs. of 60  
**02-911-437** - Cs. of 12  
32 oz. (1000mL); Cs. of 12  
**02-911-452** - Cs. of 36  
32 oz. (1000mL); Cs. of 36

### Fisherbrand Amber Wide-Mouth Packers with PTFE-Faced PE-lined Closures



**02-911-438** - Cs. of 24  
Bottles; Amber wide-mouth packers, Fisherbrand; with attached PTFE Faced PE-lined closures; 1 oz. (30mL); 28mm-400 screw cap size  
**02-911-439** - Cs. of 24  
2 oz. (60mL)  
**02-911-440** - Cs. of 24  
4 oz. (125mL)  
**02-911-441** - Cs. of 24  
8 oz. (250mL)  
**02-911-442** - Cs. of 12  
16 oz. (500mL)  
**02-911-443** - Cs. of 12  
32 oz. (1000mL)  
**02-911-796** - Cs. of 6  
42 oz. (1250mL)  
**02-911-854** - Cs. of 12  
85 oz. (2500mL)

### Fisherbrand Amber Wide-Mouth Packers without Closures



**02-912-299** - Cs. of 720  
Bottle; Amber Wide-Mouth Packer; Fisherbrand; Without closures; 0.20 oz. (6mL), 20mm-400 screw cap size  
**02-911-798** - Cs. of 624  
0.5 oz. (15mL)  
**02-911-935** - Cs. of 432  
1 oz. (30mL)  
**02-911-936** - Cs. of 216  
2 oz. (60mL)  
**02-911-888** - Cs. of 180  
4 oz. (125mL)  
**02-911-938** - Cs. of 84  
8 oz. (250mL)  
**02-911-934** - Cs. of 60  
16 oz. (500mL)  
**02-911-937** - Cs. of 36  
32 oz. (1000mL)

### Fisherbrand Certified Clean Amber Wide-Mouth Packers



**02-912-004** - Cs. of 24  
Bottles; Fisherbrand; Amber Wide-mouth packers; Certified clean; With attached PTFE Faced PE-lined closures; Cleaned per EPA guidelines; capacity: 1 oz. (30mL); 28mm-400 screw cap size  
**02-912-005** - Cs. of 24  
2 oz. (60mL)  
**02-912-007** - Cs. of 12  
4 oz. (125mL)  
**02-912-008** - Cs. of 12  
8 oz. (250mL)  
**02-912-003** - Cs. of 12  
17 oz. (500mL)  
**02-912-006** - Cs. of 12  
34 oz. (1000mL)

### Fisherbrand Certified Low-Particulate Graduated Volumetric Bottles



**02-911-865** - Cs. of 48  
Bottles; Fisherbrand; Certified low-particulate graduated volumetric; With attached polyvinyl-lined closures; Containers and caps are cleaned with USP purified water; With PTFE caps; capacity: 1 oz. (30mL)  
**02-911-866** - Cs. of 48  
With PTFE caps; Capacity: 2 oz. (60mL)  
**02-912-039** - Cs. of 24  
With PTFE caps; Capacity: 4 oz. (125mL)

**02-912-040** - Cs. of 24  
With PTFE caps; Capacity: 8 oz. (250mL)  
**02-911-880** - Cs. of 24  
With PTFE caps; Capacity: 16 oz. (500mL)  
**02-912-314** - Cs. of 24  
With PTFE-faced PE-lined caps; Capacity: 4 oz. (125mL)  
**02-912-315** - Cs. of 24  
With PTFE-faced PE-lined caps; Capacity: 8 oz. (250mL)  
**02-912-316** - Cs. of 24  
With PTFE-faced PE-lined caps; Capacity: 16 oz. (500mL)

### Fisherbrand Clear Standard Wide Mouth Bottles with Polyvinyl-lined Closures



**02-911-731** - Cs. of 24  
Bottles; Wide Mouth; Fisherbrand; Clear; Standard; Type III soda-lime glass; Rounded Shoulders; With attached Polyvinyl-lined closure, 4 oz. (125mL), 48mm-400 screw cap size  
**02-911-916** - Cs. of 24  
8 oz. (250mL), 58mm-400 screw cap size  
**02-911-917** - Cs. of 24  
16 oz. (500mL), 63mm-400 screw cap size  
**02-911-919** - Cs. of 12  
32 oz. (1000mL); 70mm-400 screw cap size  
**02-911-891** - Cs. of 6  
64 oz. (2000mL), 83mm-400 screw cap size  
**02-911-918** - Cs. of 4  
129 oz. (3810mL); 89mm-400 screw cap size

### Fisherbrand Clear Standard Wide Mouth Bottles with PTFE-Faced PE-Lined Closures



**02-911-768** - Cs. of 24  
Bottles; Wide Mouth; Fisherbrand; Clear; Standard; With attached PTFE Faced PE-lined closure; 4 oz. (125mL); 48mm-400 screw cap size  
**02-911-770** - Cs. of 24  
8 oz. (250mL); 58mm-400 screw cap size  
**02-911-767** - Cs. of 24  
16 oz. (500mL); 63mm-400 screw cap size  
**02-911-913** - Cs. of 12  
32 oz. (1000mL); 70mm-400 screw cap size  
**02-911-769** - Cs. of 6  
64 oz. (2000mL), 83mm-400 screw cap size  
**02-911-766** - Cs. of 4  
129 oz. (3810mL), 89mm-400 screw cap size

### Fisherbrand Clear Standard Wide Mouth Bottles without closures



**02-911-942** - Cs. of 24  
Bottles; Fisherbrand; Clear standard Wide-mouth bottles; Without closures; Excellent choice for general reservoir requirements; Clear for easy viewing of contents; capacity: 4 oz. (125mL); 48mm-400 screw cap size  
**02-911-943** - Cs. of 24  
8 oz. (250mL)  
**02-911-871** - Cs. of 24  
16 oz. (500mL)  
**02-912-114** - Cs. of 12  
32 oz. (1000mL)  
**02-911-869** - Cs. of 6  
64 oz. (2000mL)  
**02-912-117** - Cs. of 4  
129 oz. (3810mL)

### Fisherbrand Custom Cleaned Amber Wide-Mouth Packers



**02-911-930** - Cs. of 24  
Bottles; Fisherbrand; Amber Wide-mouth packers; Custom cleaned; With attached PTFE Faced PE-lined closures; Containers and caps are noncertified and precleaned with USP Purified Water; capacity: 2 oz. (60mL); 33mm-400 screw cap size  
**02-911-932** - Cs. of 24  
4 oz. (125mL)  
**02-911-933** - Cs. of 24  
8 oz. (250mL)  
**02-911-929** - Cs. of 12  
16 oz. (500mL)



**02-911-931** - Cs. of 12  
32 oz. (1000mL)

**Fisherbrand Graduated Valumetric Bottles with Foil-lined Closures**



**02-911-923** - Cs. of 48  
Bottles, Graduated Valumetric glass; With attached Foil-lined closures, 1 oz. (30mL), 33mm-400 screw cap size  
**02-911-924** - Cs. of 48  
2 oz. (60mL)  
**02-911-925** - Cs. of 24  
4 oz. (125mL)  
**02-911-926** - Cs. of 24  
8 oz. (250mL)  
**02-911-911** - Cs. of 24  
16 oz. (500mL)

**Fisherbrand Graduated Valumetric Bottles with Poly-Seal-lined Closures**



**02-911-463** - Cs. of 48  
Bottles; Fisherbrand; Graduated Valumetric glass bottles; With attached Poly-seal-lined closures; Wide-mouth for easy filling; Clear for easy viewing of contents; capacity: 1 oz. (30mL); 48/case; 33mm-400 screw cap size  
**02-911-470** - Cs. of 432  
Capacity: 1 oz. (30mL); 432/case  
**02-911-464** - Cs. of 48  
Capacity: 2 oz. (60mL); 48/case  
**02-911-471** - Cs. of 288  
Capacity: 2 oz. (60mL); 288/case

**Fisherbrand Graduated Valumetric Bottles with Polyvinyl-lined Closures**



**02-911-453** - Cs. of 48  
Bottle, Graduated Valumetric; Glass; Fisherbrand; Attached Polyvinyl-lined closures; 1 oz. (30mL); 33mm-400 screw cap size; Wide mouth; Clear; Round  
**02-911-465** - Cs. of 432  
1 oz. (30mL); 432/case  
**02-911-454** - Cs. of 48  
2 oz. (60mL), 48/case  
**02-911-466** - Cs. of 288  
2 oz. (60mL); 88/case  
**02-911-455** - Cs. of 24  
4 oz. (125mL), 24/case  
**02-911-467** - Cs. of 144  
4 oz. (125mL); 144/case  
**02-911-456** - Cs. of 24  
8 oz. (250mL); 24/case  
**02-911-468** - Cs. of 96  
8 oz. (250mL); 96/case  
**02-911-457** - Cs. of 24  
16 oz. (500mL); 24/case  
**02-911-469** - Cs. of 48  
16 oz. (500mL); 48/case

**Fisherbrand Graduated Valumetric Bottles with PTFE Faced PE-lined Closures**



**02-911-458** - Cs. of 48  
Bottle, Graduated Valumetric; Glass; Fisherbrand; PTFE Faced PE-lined Closures; 1 oz. (30mL); 33mm-400 screw cap size; Wide mouth; Clear; Round  
**02-911-459** - Cs. of 48  
2 oz. (60mL)  
**02-911-460** - Cs. of 24  
4 oz. (125mL)  
**02-911-461** - Cs. of 24  
8 oz. (250mL)  
**02-911-462** - Cs. of 24  
16 oz. (500mL)

**Fisherbrand Graduated Valumetric Bottles with Rubber-lined Closures**



**02-911-473** - Cs. of 48  
Bottles; Fisherbrand; Graduated valumetric glass bottles; With attached rubber-lined closures; Wide-mouth for easy filling; Autoclavable; capacity: 1 oz. (30mL); 33mm-400 screw cap size  
**02-911-474** - Cs. of 48  
Capacity: 2 oz. (60mL)  
**02-911-475** - Cs. of 24  
Capacity: 4 oz. (125mL)  
**02-911-476** - Cs. of 24  
Capacity: 8 oz. (250mL)  
**02-912-303** - Cs. of 24  
Capacity: 16 oz. (500mL)

**Fisherbrand Graduated Valumetric Bottles without Closures**



**02-911-941** - Cs. of 432  
Bottles; Fisherbrand; Graduated valumetric glass bottles; Without closures; Wide-mouth for easy filling; capacity: 1 oz. (30mL); 33mm-400 screw cap size  
**02-912-255** - Cs. of 288  
Capacity: 2 oz. (60mL)  
**02-912-256** - Cs. of 144  
Capacity: 4 oz. (125mL)  
**02-912-257** - Cs. of 96  
Capacity: 8 oz. (250mL)  
**02-912-119** - Cs. of 48  
Capacity: 16 oz. (500mL)

**Fisherbrand Safety-Coated Amber Wide-Mouth Packers**



**02-912-317** - Cs. of 24  
Bottles; Wide-mouth packers; Safety-coated; Amber glass; Not autoclavable; Shatter-resistance; Nonslip coating; Bottle only; capacity: 8 oz. (250mL); 45mm-400 screw cap size  
**02-911-906** - Cs. of 12  
Bottle only; Capacity: 16 oz. (500mL)  
**02-911-907** - Cs. of 12  
Bottle only; Capacity: 32 oz. (1000mL)  
**02-911-832** - Cs. of 24  
With polyvinyl-lined caps; Capacity: 8 oz. (250mL)  
**02-911-828** - Cs. of 12  
With polyvinyl-lined caps; Capacity: 16 oz. (500mL)  
**02-911-829** - Cs. of 12  
With polyvinyl-lined caps; Capacity: 32 oz. (1000mL)  
**02-911-863** - Cs. of 24  
With PTFE-faced PE-lined caps; Capacity: 8 oz. (250mL)  
**02-911-830** - Cs. of 12  
With PTFE-faced PE-lined caps; Capacity: 16 oz. (500mL)  
**02-911-831** - Cs. of 12  
With PTFE-faced PE-lined caps; Capacity: 32 oz. (1000mL)



# Round Plastic Narrow Mouth Bottles

## Fisherbrand Amber Narrow-Mouth HDPE Bottles



- 03-007-46** - Pk. of 12, Cs. of 6  
Bottle; Narrow-Mouth; Fisherbrand; HDPE; Amber; Protect light-sensitive materials during sampling, shipping or storage; 125mL
- 03-007-47** - Pk. of 12, Cs. of 6  
250mL
- 03-007-48** - Pk. of 12, Cs. of 4  
500mL
- 03-007-49** - Pk. of 6, Cs. of 4  
1L

## Fisherbrand Narrow-Mouth Field Sample Bottles



- 02-895A** - Cs. of 500  
Bottle, Field-Sample; Fisherbrand; HDPE; Bulk-Packaged; Screw Cap Size: 24mm; Case of: 500; Capacity: 4 oz.
- 02-895B** - Cs. of 250  
8 oz. (250mL)
- 02-895C** - Cs. of 125  
16 oz. (500mL)
- 02-895D** - Cs. of 50  
32 oz. (1000mL)

# Round Plastic Wide Mouth Bottles

## Fisherbrand Amber Wide-Mouth HDPE Bottles



- 03-007-52** - Pk. of 12, Cs. of 6  
Bottle; Wide-Mouth; Fisherbrand; HDPE; Amber; 125mL
- 03-007-53** - Pk. of 12, Cs. of 6  
250mL
- 03-007-54** - Pk. of 12, Cs. of 4  
500mL
- 03-007-55** - Pk. of 6, Cs. of 4  
1L

## Fisherbrand HDPE Round Bottles



- 02-896-2A** - Pk. of 12, Cs. of 6  
Bottle, Round; Fisherbrand; HDPE; Leakproof; May be radiation, EtO, or chemically sterilized; Not autoclavable; With polypropylene screw cap; Capacity: 1 oz. (30mL); Screw cap size: 28mm; Translucent
- 02-896-2B** - Pk. of 12, Cs. of 6  
Capacity: 2 oz.
- 02-896-2D** - Pk. of 12, Cs. of 6  
Capacity: 8 oz.
- 02-896-2C** - Pk. of 12, Cs. of 6  
Capacity: 4 oz.
- 02-896-2E** - Pk. of 12, Cs. of 4  
Capacity: 16 oz.
- 02-896-2F** - Pk. of 6, Cs. of 4  
Capacity: 32 oz.

## Fisherbrand Wide Mouth LDPE Bottles



- 02-896-1A** - Pk. of 12, Cs. of 6  
Bottle; Fisherbrand; LDPE; Ungraduated; Wide-Mouth; Screw Cap Size: 28mm; Capacity: 1 oz.
- 02-896-1B** - Pk. of 12, Cs. of 6  
2 oz.
- 02-896-1C** - Pk. of 12, Cs. of 6  
4 oz.
- 02-896-1D** - Pk. of 12, Cs. of 6  
8 oz.
- 02-896-1E** - Pk. of 12, Cs. of 4  
16 oz.
- 02-896-1F** - Pk. of 6, Cs. of 4  
32 oz.

## Fisherbrand Wide Mouth Polypropylene Bottles



- 02-896A** - Pk. of 12, Cs. of 6  
Bottle, General-purpose; Fisherbrand; Wide-mouth; Polypropylene; With polypropylene linerless screw closures; Rigid walls; Chemical resistant; Capacity: 1 oz. (30mL); Screw cap size: 43mm
- 02-896B** - Pk. of 12, Cs. of 6  
2 oz.
- 02-896C** - Pk. of 12, Cs. of 6  
4 oz.
- 02-896D** - Pk. of 12, Cs. of 6  
8 oz.
- 02-896E** - Pk. of 12, Cs. of 4  
16 oz.
- 02-896F** - Pk. of 6, Cs. of 4  
32 oz.

## Fisherbrand Wide-Mouth Field Sample Bottles



- 02-895-1A** - Cs. of 500  
Bottles; HDPE; Field-Sample; Fisherbrand; Bulk-Packaged; Screw Cap Size: 38mm; Case of: 500; Capacity: 4 oz. (125mL)
- 02-895-1B** - Cs. of 250  
8 oz.
- 02-895-1C** - Cs. of 125  
16 oz.
- 02-895-1D** - Cs. of 50  
32 oz.

# Square Glass Wide Mouth Bottles

## Fisherbrand Certified Low-Particulate Clear French Squares



- 02-911-890** - Cs. of 24  
Bottles; Fisherbrand; Certified low-particulate; Clear French; Containers and caps are cleaned with USP purified water; PTFE-lined caps included; With Certificate of Analysis; Capacity: 4 oz. (125mL); 33mm-415 screw cap size
- 02-911-893** - Cs. of 24  
8 oz. (250mL)
- 02-911-885** - Cs. of 12  
32 oz. (1000mL)

**Fisherbrand Clear French Squares with Poly-Seal-lined Closures**



- 02-911-402** - Cs. of 48  
Bottles; Fisherbrand; Clear french squares with Poly-seal-lined; Wide-mouth; Easy to fill; Clear glass allow easy viewing of contents; capacity: 0.5 oz. (15mL); 48/case; 20mm-400 screw cap size
- 02-911-403** - Cs. of 48  
1 oz. (30mL); Cs. of 48
- 02-911-404** - Cs. of 48  
2 oz. (60mL); Cs. of 48
- 02-911-405** - Cs. of 24  
4 oz. (125mL); Cs. of 24
- 02-911-875** - Cs. of 120  
4 oz. (125mL); Cs. of 120

**Fisherbrand Clear French Squares with Polyvinyl-lined Closures**



- 02-912-302** - Cs. of 576  
Bottles; Fisherbrand; Clear french squares with polyvinyl-lined; Wide-mouth; Easy to fill; Clear glass allow easy viewing of contents; capacity: 0.5 oz. (15mL); 576/case; 20mm-400 screw cap size
- 02-911-379** - Cs. of 48  
0.5 oz. (15mL); Cs. of 48
- 02-911-380** - Cs. of 48  
1 oz. (30mL); Cs. of 48
- 02-911-397** - Cs. of 280  
1 oz. (30mL); Cs. of 280
- 02-911-381** - Cs. of 48  
2 oz. (60mL); Cs. of 48
- 02-911-398** - Cs. of 240  
2 oz. (60mL); Cs. of 240
- 02-911-382** - Cs. of 24  
4 oz. (125mL); Cs. of 24
- 02-911-399** - Cs. of 120  
4 oz. (125mL); Cs. of 120
- 02-911-383** - Cs. of 24  
8 oz. (250mL); Cs. of 24
- 02-911-400** - Cs. of 84  
8 oz. (250mL); Cs. of 84
- 02-911-384** - Cs. of 24  
16 oz. (500mL); Cs. of 24
- 02-911-401** - Cs. of 40  
16 oz. (500mL); Cs. of 40
- 02-912-017** - Cs. of 12  
32 oz. (1000mL); Cs. of 12

**Fisherbrand Clear French Squares with PTFE Faced PE-lined Closures**



- 02-911-385** - Cs. of 48  
Bottle; Clear French Square; Fisherbrand; with attached Teflon Faced PE-lined closures, 1/2 oz. (15mL), 20mm-400 screw cap size
- 02-911-386** - Cs. of 48  
1 oz. (30mL)
- 02-911-387** - Cs. of 48  
2 oz. (60mL)
- 02-911-388** - Cs. of 24  
4 oz. (125mL)
- 02-911-389** - Cs. of 24  
8 oz. (250mL)
- 02-911-390** - Cs. of 24  
16 oz. (500mL)
- 02-912-016** - Cs. of 12  
32 oz. (1000mL)

**Fisherbrand Clear French Squares with Rubber-lined Closures**



- 02-911-391** - Cs. of 48  
Bottle; Clear French Square; Fisherbrand; with attached Rubber-lined closures, 0.5 oz. (15mL); 20mm-400 screw cap size
- 02-911-392** - Cs. of 48  
1 oz. (30mL)
- 02-911-393** - Cs. of 48  
2 oz. (60mL)
- 02-911-394** - Cs. of 24  
4 oz. (125mL)
- 02-911-395** - Cs. of 24  
8 oz. (250mL)
- 02-911-396** - Cs. of 24  
16 oz. (500mL)
- 02-912-300** - Cs. of 12  
32 oz. (1000mL)

**Fisherbrand Clear French Squares without Closures**



- 02-912-301** - Cs. of 576  
Bottle; Clear French Square; Fisherbrand; Without closures, 1/2 oz. (15mL), 20mm-400 screw cap size
- 02-911-946** - Cs. of 280  
1 oz. (30mL)
- 02-911-947** - Cs. of 240  
2 oz. (60mL)
- 02-911-887** - Cs. of 120  
4 oz. (125mL)
- 02-912-297** - Cs. of 84  
8 oz. (250mL)
- 02-911-877** - Cs. of 40  
16 oz. (500mL)
- 02-911-884** - Cs. of 24  
32 oz. (1000mL)

**Fisherbrand Clear Glass Tablet Squares**



- 02-911-881** - Cs. of 288  
Bottles; Fisherbrand; Clear glass tablet squares; Packed with caps attached; With Polyvinyl lined; capacity: 1 oz. (30mL); 33mm-400 screw cap size
- 02-911-993** - Cs. of 288  
Capacity: 1 oz. (30mL); With Poly-seal-lined caps
- 02-911-994** - Cs. of 288  
Capacity: 2 oz. (60mL); With PTFE-faced, PE-lined caps; 288/Cs.
- 02-911-995** - Cs. of 144  
Capacity: 2 oz. (60mL); With PTFE-faced, PE-lined caps; 144/Cs.
- 02-912-296** - Cs. of 144  
Capacity: 2 oz. (60mL); Bottles only
- 02-912-308** - Cs. of 144  
Capacity: 2 oz. (60mL); With Poly-seal-lined caps

**Boxes**

**Fisherbrand Cryo Storage Boxes**



- 03-395-402** - Pk. of 5, Cs. of 10  
Cryogenic Storage Box; Polypropylene; 81 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 12.7L x 12.7W x 5.1cm H; 5/Pk.; Blue
- 14-826-340** - Pk. of 5, Cs. of 10  
81-Place; Green
- 03-395-400** - Pk. of 5, Cs. of 10  
81-Place; Natural
- 03-395-404** - Pk. of 5, Cs. of 10  
81-Place; Orange
- 03-395-403** - Pk. of 5, Cs. of 10  
81-Place; Purple

- 03-395-401** - Pk. of 5, Cs. of 10  
81-Place; Assorted
- 03-391-201** - Pk. of 5, Cs. of 10  
100-Place; Blue
- 03-391-202** - Pk. of 5, Cs. of 10  
100-Place; Green
- 03-391-200** - Pk. of 5, Cs. of 10  
100-Place; Natural
- 03-391-204** - Pk. of 5, Cs. of 10  
100-Place; Orange
- 03-391-203** - Pk. of 5, Cs. of 10  
100-Place; Purple
- 03-391-205** - Pk. of 5, Cs. of 10  
100-Place; Assorted

### Fisherbrand Cryo/Freezer Boxes



- 03-395-100** - Ea.  
Fiberbox, Cryo; Argos Technologies; White; 3 x 3 x 2 in.
- 03-395-112** - Ea.  
3 x 3 x 2 in.; with Dividers
- 03-395-102** - Ea., Cs. of 200  
5 x 5 x 2 in.
- 03-395-113** - Ea., Cs. of 200  
5 x 5 x 2 in.; with 81 Dividers
- 03-395-114** - Ea., Cs. of 200  
5 x 5 x 2 in.; with Dividers
- 03-395-103** - Ea.  
5 x 5 x 3 in.
- 03-395-101** - Ea.  
5 x 5 x 4 in.

### Fisherbrand Dividers for Cryo/Freezer Boxes

- 03-395-108** - Ea., Dozen  
Cell Dividers, Cryo, Argos Technologies, 16 Cell
- 03-395-110** - Ea., Dozen  
25 Cell
- 03-395-109** - Ea., Dozen  
36 Cell
- 03-395-107** - Ea., Dozen  
49 Cell
- 03-395-106** - Ea., Dozen  
64 Cell
- 03-395-105** - Ea., Dozen  
81 Cell
- 03-395-104** - Ea., Dozen  
100 Cell
- 03-395-111** - Ea., Dozen  
144 Cell

## Brushes

### Fisherbrand Adaptable-To-Container Brush



- 03-565** - Ea., Pk. of 12  
Brush, Adaptable-to-Container; Fisherbrand; Single-row bristles; Heavy wire handle with swivel joint; Overall Length: 12 in.; Bristle dia. 1-1/2 x 4-1/4 in. L

### Fisherbrand Stainless-Steel Handled Brushes



- 03-578** - Pk. of 12, Cs. of 12  
Brush, Test tube; Fisherbrand; Stiff, white nylon bristles; Radial tufted end; Stainless-steel handle is rustproof, corrosion resistant; Bristles: 3/4 x 3-1/2 in.; Length: 9 in.
- 03-572-5** - Pk. of 12, Cs. of 12  
Overall Length: 11 in.; Bristle L x D: 4 x 1-3/8 in.

## Buffers

### Buffer-Pac Color-Coded Solutions (Certified), Fisher Chemical

- SB105** - Pk. of 3, Cs. of 6  
Buffer-Pac Color-coded Solutions; Certified; Fisher Chemical; Poly Bottles; 1 each: SB101-500, pH 4.00 (Red), SB107-500, pH 7.00 (Yellow); SB115-500, pH 10.00 (Blue)

## Carboys and Tanks

### Fisherbrand Polypropylene Wide-Mouth Carboys with Spigots



- 12-565-638** - Ea., Cs. of 6  
Carboy; Fisherbrand; Polypropylene; Wide-Mouth; with spigot; 100mm closure. Capacity: 1-1/4 gal. (5L)
- 12-565-639** - Ea., Cs. of 6  
Capacity: 2-5/8 gal. (10L)
- 12-565-640** - Ea., Cs. of 4  
Capacity: 5-1/4 gal. (20L)

### Fisherbrand Polypropylene Wide-Mouth Carboys without Spigots



- 12-565-633** - Ea., Cs. of 6  
Carboys, Wide-mouth; Without spigot; Polypropylene; Fisherbrand; With polypropylene screw cap; Graduated TO DELIVER in liters; Chemical resistant; Molded handgrips; Leakproof; Capacity: 2-5/8 gal. (10L); Closure size: 100mm
- 12-565-634** - Ea., Cs. of 4  
Capacity: 5 1/4 gal. (20L)

## Carts

### Fisher Scientific Laboratory Carts



- 22-363-250** - Ea.  
Cart, Laboratory; Rugged, welded 1 in. square, epoxy-coated steel frame with off-white powder-coated finish; 400 lb. capacity; Chemical Lab Cart
- 22-363-251** - Ea.  
Glassware Cart with 2 Baskets
- 22-363-253** - Ea.  
Glassware Cart with 4 Baskets
- 22-363-255** - Ea.  
Laboratory Cart with 2 Removable SS Trays



22-363-252 - Ea.  
Pan Cart  
22-363-254 - Ea.  
Portable Table

# Cell Scrapers

## Fisherbrand Cell Lifters



**08-773-1** - Cs. of 100  
Cell Lifter; Disposable; Polyethylene; Fisherbrand; Blade Length: 2cm; Handle Length: 18cm; Sterile; Individually wrapped

## Fisherbrand Cell Scrapers



**08-773-2** - Cs. of 100  
Cell Scraper; Fisherbrand; Handle: Polystyrene; Blade: Polyethylene; Blade Length: 1.8cm; Handle Length: 25cm; Sterile; Individually wrapped  
**08-773-3** - Cs. of 100  
O.D. x L: 16 x 100mm; 250/tray


## Fisherbrand Cell Spreaders



**08-100-10** - Pk. of 25, Cs. of 4  
Cell Spreader; Fisherbrand; Width: 30mm; Autoclavable; Nonsterile  
**08-100-12** - Pk. of 25, Cs. of 4  
Width: 30mm; Sterile  
**08-100-11** - Pk. of 25, Cs. of 4  
Width: 60mm; Nonsterile  
**08-100-13** - Pk. of 25, Cs. of 4  
Width: 60mm; Sterile

# Cell Strainers


## Fisherbrand Cell Strainers



**22-363-547** - Cs. of 50  
Cell Strainer; Fisherbrand; 40um; Color: Blue; Polypropylene frame with nylon mesh; Sterile; Individually packed; Universal fit for all 50mL conical tubes  
**22-363-548** - Cs. of 50  
70um  
**22-363-549** - Cs. of 50  
100um

# Chromatography Supplies

## Fisherbrand Pure Cellulose Chromatography Paper



**05-714-1** - Pk. of 100  
Paper, Chromatography; Fisherbrand; Thickness: 0.19mm (thin); Pure Cellulose; 46 x 57cm (18 x 22 in.) sheets  
**05-714-4** - Pk. of 100  
Thickness: 0.35mm

# Closures

## Caps

### Fisherbrand EnviroWare Borosilicate Glass EPA Vials Replacement Caps

**03-340-7L** - Cs. of 144  
Caps; Kimble; Replacement; For Enviroware Borosilicate Glass EPA Vials; Style: White polypropylene screw thread closures; Cap Type: With open tops; GPI Thread: 24-400

### Fisherbrand Screw Caps for Disposable Glass Tubes



**14-957-82E** - Cs. of 1000  
Cap, Phenolic Screw Cap with Rubber Liner. 18-415GPI Thread., Fisherbrand  
**14-959-36B** - Cs. of 1000  
GPI Thread: 15-415; Polypropylene  
**14-957-82A** - Cs. of 1000  
Phenolic, 13-425 Thread  
**14-959-36A** - Cs. of 1000  
Phenolic, 15-415 Thread

## Corks

### Fisherbrand Cork Stoppers - Laboratory Grade



**07-782C** - Pk. of 500  
Stoppers, Cork; Fisherbrand; Size: 0; Laboratory Grade  
**07-782B** - Pk. of 500  
Size: 00; Laboratory Grade  
**07-782D** - Pk. of 500  
Size: 1; Laboratory Grade  
**07-782E** - Pk. of 500  
Size: 2; Laboratory Grade  
**07-782A** - Pk. of 500  
Size: 3/0; Laboratory Grade  
**07-782F** - Pk. of 500  
Size: 3; Laboratory Grade  
**07-782G** - Pk. of 500  
Size: 4; Laboratory Grade  
**07-782H** - Pk. of 500  
Size: 5; Laboratory Grade  
**07-782J** - Pk. of 500  
Size: 6; Laboratory Grade  
**07-782K** - Pk. of 500  
Size: 7; Laboratory Grade  
**07-781Y** - Pk. of 100  
Size: 20; Laboratory Grade

### Fisherbrand Cork Stoppers - Superior Grade



**07-791A** - Pk. of 100  
Stoppers, Cork; Fisherbrand; Size: 0; Superior Grade  
**07-791B** - Pk. of 100  
Size: 1; Superior Grade  
**07-791C** - Pk. of 100  
Size: 2; Superior Grade  
**07-791D** - Pk. of 100  
Size: 3; Superior Grade  
**07-791E** - Pk. of 100  
Size: 4; Superior Grade  
**07-791F** - Pk. of 100  
Size: 5; Superior Grade  
**07-791G** - Pk. of 100  
Size: 6; Superior Grade  
**07-791H** - Pk. of 100  
Size: 7; Superior Grade  
**07-791J** - Pk. of 100  
Size: 8; Superior Grade  
**07-791K** - Pk. of 100  
Size: 9; Superior Grade  
**07-791L** - Pk. of 100  
Size: 10; Superior Grade  
**07-791M** - Pk. of 100  
Size: 11; Superior Grade



## Rubber Stoppers

### Fisherbrand One-Hole Rubber Stoppers



**14-135B** - Pk. of 80, Cs. of 5  
Stopper, Rubber; Fisherbrand; Size No: 0; Hole Size:  
5mm

**14-135A** - Pk. of 109, Cs. of 5  
Size No: 00; Hole Size: 3mm

**14-135C** - Pk. of 59, Cs. of 5  
Size No: 1; Hole Size: 4mm

**14-135D** - Pk. of 54, Cs. of 5  
Size No: 2; Hole Size: 5mm

**14-135E** - Pk. of 41, Cs. of 5  
Size No: 3; Hole Size: 5mm

**14-135F** - Pk. of 35, Cs. of 5  
Size No: 4; Hole Size: 5mm

**14-135G** - Pk. of 29, Cs. of 5  
Size No: 5; Hole Size: 5mm

**14-135H** - Pk. of 27, Cs. of 5  
Size No: 5-1/2; Hole Size: 5mm

**14-135J** - Pk. of 21, Cs. of 5  
Size No: 6; Hole Size: 5mm

**14-135K** - Pk. of 19, Cs. of 5  
Size No: 6-1/2; Hole Size: 5mm

**14-135L** - Pk. of 15, Cs. of 5  
Size No: 7; Hole Size: 5mm

**14-135M** - Pk. of 12, Cs. of 5  
Size No: 8; Hole Size: 5mm

**14-135N** - Pk. of 10, Cs. of 5  
Size No: 9; Hole Size: 5mm

**14-135P** - Pk. of 8, Cs. of 5  
Size No: 10; Hole Size: 5mm

**14-135Q** - Pk. of 6, Cs. of 5  
Size No: 11; Hole Size: 5mm

**14-135R** - Pk. of 5, Cs. of 5  
Size No: 12; Hole Size: 5mm

**14-135S** - Pk. of 4, Cs. of 5  
Size No: 13; Hole Size: 5mm

**14-135T** - Ea., Cs. of 8  
Size No: 14; Hole Size: 5mm

**14-135U** - Ea., Cs. of 7

Size No: 15; Hole Size: 5mm

**14-141T** - Ea., Cs. of 13

Size No.: 13-1/2

**14-141V** - Ea., Cs. of 8

Size No.: 14

**14-141W** - Cs. of 7

Size No.: 15

**14-141X** - Cs. of 4

Size No.: 16

### Fisherbrand Two-Hole Rubber Stoppers

**14-140A** - Pk. of 113, Cs. of 5

Stoppers, Rubber Two Hole; Fisherbrand; Size No.: 00; Hole Size: 3mm

**14-140B** - Pk. of 82, Cs. of 5

Size No.: 0; Hole Size: 3mm

**14-140C** - Pk. of 63, Cs. of 5

Size No.: 1; Hole Size: 4mm

**14-140D** - Pk. of 57, Cs. of 5

Size No.: 2; Hole Size: 5mm

**14-140E** - Pk. of 43, Cs. of 5

Size No.: 3; Hole Size: 5mm

**14-140F** - Pk. of 37, Cs. of 5

Size No.: 4; Hole Size: 5mm

**14-140G** - Pk. of 30, Cs. of 5

Size No.: 5; Hole Size: 5mm

**14-140H** - Pk. of 22, Cs. of 5

Size No.: 6; Hole Size: 5mm

**14-140J** - Pk. of 16, Cs. of 5

Size No.: 7; Hole Size: 5mm

**14-140K** - Pk. of 13, Cs. of 5

Size No.: 8; Hole Size: 5mm

**14-140L** - Pk. of 11, Cs. of 5

Size No.: 9; Hole Size: 5mm

**14-140M** - Pk. of 8, Cs. of 5

Size No.: 10; Hole Size: 5mm

**14-140N** - Pk. of 6, Cs. of 5

Size No.: 11; Hole Size: 5mm

**14-140P** - Pk. of 5, Cs. of 5

Size No.: 12; Hole Size: 5mm

**14-140Q** - Pk. of 4, Cs. of 5

Size No.: 13; Hole Size: 5mm

**14-140R** - Ea., Cs. of 8

Size No.: 14; Hole Size: 5mm

**14-140S** - Ea., Cs. of 7

Size No.: 15; Hole Size: 5mm

### Fisherbrand Solid Neoprene Rubber Stoppers



**14-141A** - Pk. of 68, Cs. of 5  
Stopper, Rubber; Fisherbrand; Size No.: 0; Synthetic  
neoprene resists petroleum products, fats, oils, and  
organic chemicals

**14-141B** - Pk. of 48, Cs. of 5  
Size No.: 1

**14-141C** - Pk. of 44, Cs. of 5  
Size No.: 2

**14-141D** - Pk. of 32, Cs. of 5  
Size No.: 3

**14-141E** - Pk. of 27, Cs. of 5  
Size No.: 4

**14-141F** - Pk. of 23, Cs. of 5  
Size No.: 5

**14-141G** - Pk. of 17, Cs. of 5  
Size No.: 6

**14-141H** - Pk. of 15, Cs. of 5  
Size No.: 6-1/2

**14-141J** - Pk. of 12, Cs. of 5  
Size No.: 7

**14-141K** - Pk. of 10, Cs. of 5  
Size No.: 8

**14-141L** - Pk. of 8, Cs. of 5  
Size No.: 9

**14-141M** - Pk. of 8, Cs. of 5  
Size No.: 10

**14-141N** - Pk. of 7, Cs. of 5  
Size No.: 10-1/2

**14-141P** - Pk. of 6, Cs. of 5  
Size No.: 11

**14-141Q** - Pk. of 4, Cs. of 5  
Size No.: 11-1/2

**14-141R** - Pk. of 5, Cs. of 5  
Size No.: 12

**14-141S** - Pk. of 4, Cs. of 5  
Size No.: 13

## Counting Devices

### Fisherbrand Digital Counter



**07-905-2** - Ea.

Counter, Digital; Fisherbrand Reads to 99,999; Water-  
and shock-resistant; Easy-to-read 1/4 in. LCD; Solid-state  
construction; Battery power for two years' continuous  
service

### Fisherbrand Hand Tally Counter



**07-905-6** - Ea.

Counter, Hand Tally; Fisherbrand; Tallies to 9999; Resets  
to 0; Pocket size; Chrome finish

## Cuvets

### Fisherbrand Electroporation Cuvettes Plus



- FB101** - Pk. of 50  
Cuvettes, Electroporation; Fisherbrand Cuvettes Plus; 1mm gap with imbedded aluminum electrodes, Capacity: 90uL; Field strengths up to 25 kV/cm, Gamma irradiated for guaranteed sterility, Disposable
- FB102** - Pk. of 50  
2mm gap with imbedded aluminum electrodes, Capacity: 400uL
- FB104** - Pk. of 50  
4mm gap with imbedded aluminum electrodes, Capacity: 800uL

## Cylinders

### Glass Graduated Cylinders Rated To Contain

#### Fisherbrand Class B Standard Cylinders, TC

- 08-557-1A** - Pk. of 2, Cs. of 12  
Cylinder, Graduated; Fisherbrand; Class B Glass; Calibrated To Contain; White-ceramic graduations; Handy pouring spouts; Capacity: 10mL; Subdivisions: 0.2mL
- 08-557-1B** - Ea., Cs. of 18  
Capacity: 25mL; Subdivisions: 0.2mL
- 08-557-1C** - Ea., Cs. of 18  
Capacity: 50mL; Subdivisions: 1mL
- 08-557-1D** - Ea., Cs. of 12  
Capacity: 100mL; Subdivisions: 1mL
- 08-557-1E** - Ea., Cs. of 12  
Capacity: 250mL; Subdivisions: 2mL
- 08-557-1F** - Ea., Cs. of 8  
Capacity: 500mL; Subdivisions: 5mL
- 08-557-1G** - Ea., Cs. of 4  
Capacity: 1000mL; Subdivisions: 10mL
- 08-557-1H** - Ea., Cs. of 4  
Capacity: 2000mL; Subdivisions: 20mL

### Glass Graduated Cylinders Rated To Deliver

#### Fisherbrand Cylinders with Double Metric Scales, TD



- 08-555A** - Ea., Cs. of 36  
Cylinders, Graduated; Fisherbrand; Class B; Glass; With double ceramic metric scales; Calibrated To Deliver; Heavy base; Pouring lip; Capacity: 10mL; Subdivisions: 0.2mL; Tolerance: +/-0.1mL
- 08-555B** - Ea., Cs. of 36  
Capacity: 25mL; Subdivisions: 0.5mL; Tolerance: +/-0.20mL
- 08-555C** - Ea., Cs. of 36  
Capacity: 50mL; Subdivisions: 1mL; Tolerance: +/-0.34mL
- 08-555D** - Ea., Cs. of 36  
Capacity: 100mL; Subdivisions: 1mL; Tolerance: +/-1.0mL
- 08-555E** - Ea., Cs. of 18  
Capacity: 250mL; Subdivisions: 2mL; Tolerance: +/-2.0mL
- 08-555F** - Ea., Cs. of 12  
Capacity: 500mL; Subdivisions: 5mL; Tolerance: +/-4.0mL
- 08-555G** - Ea., Cs. of 8  
Capacity: 1000mL; Subdivisions: 10mL; Tolerance: +/-6.0mL


## Plastic Graduated Cylinders

### Fisherbrand Polymethylpentene Graduated Cylinders



- 03-007-30** - Ea., Cs. of 24  
Cylinder, Graduated; Fisherbrand; Polymethylpentene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 10mL; Subdivisions: 0.2mL
- 03-007-31** - Ea., Cs. of 18  
Capacity: 25mL; Subdivisions: 0.5mL
- 03-007-32** - Ea., Cs. of 18  
Capacity: 50mL; Subdivisions: 1mL
- 03-007-33** - Ea., Cs. of 12  
Capacity: 100mL; Subdivisions: 1mL
- 03-007-34** - Ea., Cs. of 12  
Capacity: 250mL; Subdivisions: 2mL
- 03-007-35** - Ea., Cs. of 8  
Capacity: 500mL; Subdivisions: 5mL
- 03-007-36** - Ea., Cs. of 6  
Capacity: 1000mL; Subdivisions: 10mL
- 03-007-37** - Ea., Cs. of 4  
Capacity: 2000mL; Subdivisions: 20mL

### Fisherbrand Polypropylene Graduated Cylinders



- 03-007-38** - Ea., Cs. of 12  
Cylinder, Graduated; Fisherbrand; Polypropylene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 10mL; Subdivisions: 0.2mL
- 03-007-39** - Ea., Cs. of 12  
Capacity: 25mL; Subdivisions: 0.5mL
- 03-007-40** - Ea., Cs. of 12  
Capacity: 50mL; Subdivisions: 1mL
- 03-007-41** - Ea., Cs. of 8  
Capacity: 100mL; Subdivisions: 1mL
- 03-007-42** - Ea., Cs. of 8  
Capacity: 250mL; Subdivisions: 2mL
- 03-007-43** - Ea., Cs. of 6  
Capacity: 500mL; Subdivisions: 5mL
- 03-007-44** - Ea., Cs. of 4  
Capacity: 1000mL; Subdivisions: 10mL
- 03-007-45** - Ea., Cs. of 2  
Capacity: 2000mL; Subdivisions: 20mL

## Detergents and Cleaners

### Fisherbrand Sparkleen 1 Detergent



- 04-320-5** - Ea.  
Detergent; Fisherbrand Sparkleen 1; 100 lb. (45kg); For manual washing; Removes organic and inorganic deposits; Contains phosphates

## Fisherbrand Versa-Clean Concentrate



**04-342** - Ea., Cs. of 12  
Cleaner; Fisherbrand Versa-Clean; 33-7/8 oz.; For use in ultrasonic cleaners, high-pressure sprayers, scrubbers/soak tanks; For use on plastic, glass, rubber, vinyl, most metals; Biodegradeable; Phosphate free; Nontoxic; Noncaustic  
**04-343** - Ea., Cs. of 4  
1-1/16 gal. (4L)

## Eye Protection

### Fisherbrand Impact and Chemical Splash Goggles



**19-181-503** - Ea., Pk. of 12  
Goggles, Impact; Fisherbrand; Soft PVC frame; Adjustable elastic strap; Direct vent; Polycarbonate lens has antifog, antiscratch coating  
**19-181-504** - Ea., Pk. of 12  
Chemical Splash

### Fisherbrand 200 Series Spectacles



**19-130-2088** - Pk. of 12, Cs. of 10  
Spectacles, Safety; Fisherbrand 200 Series; One-piece lens styling; Molded-in sideshields and browguards; Adjustable spatula temples; Black frame; Clear Hardcoat lens; Absorbs >99.9% of UV light; Meets ANSI Z87.11989 and CSA Z94.3  
**19-130-2089** - Pk. of 12, Cs. of 10  
Blue frame; Clear Hardcoat Lens

## Film and Foil Wrapping

### Fisherbrand Aluminum Foil



**01-213** - Ea., Cs. of 12  
Foil, Aluminum; Fisherbrand; 0.00045 x 3 in. (0.001 x 7.6cm); 1/2 lb. Roll; Width: 7.6cm; Length: 96m; Thickness: 0.011mm  
**01-213-1** - Pk. of 1000, Cs. of 10  
Sheets; L x W: 6 x 6 in  
**01-213-100** - Ea., Cs. of 24  
Standard-gauge roll; 12 in. W x 25 ft. L  
**01-213-101** - Ea., Cs. of 16  
Standard-gauge roll; 12 in. W x 200 ft. L  
**01-213-102** - Ea., Cs. of 9  
Standard-gauge roll; 12 in. W x 500 ft. L  
**01-213-103** - Ea., Cs. of 12  
Heavy-gauge roll; 18 in. W x 25 ft. L  
**01-213-104** - Ea., Cs. of 12  
Heavy-gauge roll; 18 in. W x 50 ft. L  
**01-213-105** - Ea.  
Heavy-gauge roll; 18 in. W x 500 ft. L

### Fisherbrand Clear Plastic Lab Wrap



**22-305-654** - Ea.  
Wrap, Plastic; Fisherbrand; High-cling and low-tangle properties; 100% polyethylene; Ultraclean; Carton with metal tear-off bar; Clear; Roll W x L: 12 in. x 100 ft. (30.5cm x 30.5m)

## Filters

### Centrifugal Filter Devices

#### Fisherbrand Single Cytology Funnels



**10-358** - Pk. of 50, Cs. of 10  
Funnel, Cytology; Single; Fisherbrand; With Brown Filter Card; volumes up to 0.25mL or less; sealing caps included; two deposition areas on each slide; For the Shandon CytoSpin Centrifuge; Case  
**10-354** - Pk. of 50, Cs. of 10  
With white filter card; 50/Pk.  
**10-355** - Cs. of 500  
With white filter card; 500/Cs.

#### Fisherbrand Double Cytology Funnel with White Filter Card



**10-356** - Pk. of 25, Cs. of 10  
Funnel, Cytology; Double; Fisherbrand; With White Filter Card; volumes up to 0.25mL or less; sealing caps included; two deposition areas on each slide; For the Shandon CytoSpin Centrifuge; Pack  
**10-353** - Cs. of 500  
Double; Bulk

## Glass Fiber Filters and Prefilters

### Fisherbrand Glass Fiber Circles



**09-804-21B** - Pk. of 100  
Filters; Glass Fiber Circles; Fisherbrand; Ideal as prefilters and for gravimetric analyses; Fast flow; High retention; Diameter: 21mm; Grade G2  
**09-804-21C** - Pk. of 100  
Diameter: 21mm; Grade G4  
**09-804-21A** - Pk. of 100  
Diameter: 21mm; Grade G6  
**09-804-21D** - Pk. of 100  
Diameter: 21mm; Grade G8  
**09-804-24C** - Pk. of 100  
Diameter: 24mm; Grade G4  
**09-804-24A** - Pk. of 100  
Diameter: 24mm; Grade G6  
**09-804-25D** - Pk. of 100  
Diameter: 25mm; Grade G8  
**09-804-42B** - Pk. of 100  
Diameter: 42.5mm; Grade G2  
**09-804-42C** - Pk. of 100  
Diameter: 42.5mm; Grade G4  
**09-804-42A** - Pk. of 100  
Diameter: 42.5mm; Grade G6  
**09-804-42D** - Pk. of 100  
Diameter: 42.5mm; Grade G8  
**09-804-55B** - Pk. of 100  
Diameter: 55mm; Grade G2  
**09-804-55C** - Pk. of 100  
Diameter: 55mm; Grade G4  
**09-804-55A** - Pk. of 100  
Diameter: 55mm; Grade G6  
**09-804-55D** - Pk. of 100  
Diameter: 55mm; Grade G8

- 09-804-70C** - Pk. of 100  
Diameter: 70mm; Grade G4
- 09-804-70A** - Pk. of 100  
Diameter: 70mm; Grade G6
- 09-804-70D** - Pk. of 100  
Diameter: 70mm; Grade G8
- 09-804-90C** - Pk. of 100  
Diameter: 90mm; Grade G4
- 09-804-90A** - Pk. of 100  
Diameter: 90mm; Grade G6
- 09-804-110A** - Pk. of 100  
Diameter: 110mm; Grade G6
- 09-804-150A** - Pk. of 100  
Diameter: 150mm; Grade G6

## Qualitative Filter Paper

### Fisherbrand Fluted Qualitative Circles



- 09-790-14D** - Pk. of 100  
Filter Paper; Fluted; Qualitative Circles; Fisherbrand;  
Diameter: 12.5cm; Grade: P8--Fluted
- 09-790-14E** - Pk. of 100  
Diameter: 15cm
- 09-790-14F** - Pk. of 100  
Diameter: 18.5cm
- 09-790-14G** - Pk. of 100  
Diameter: 24cm

### Fisherbrand Qualitative Grade Plain Circles and Sheets - P2 Grade



- 09-803-5A** - Pk. of 100  
Filter Paper, Qualitative; Fisherbrand; Grade P2;  
Cellulose fiber; Particle Retention: 1-5um; Porosity: Fine;  
Filtration Speed: Herzberg, 1400 seconds; Flowrate:  
Slow-5mL/min.; Circles; 5.5cm dia.
- 09-803-5B** - Pk. of 100  
Circles; 7cm dia.
- 09-803-5C** - Pk. of 100  
Circles; 9cm dia.
- 09-803-5D** - Pk. of 100  
Circles; 11cm dia.
- 09-803-5E** - Pk. of 100  
Circles; 12.5cm dia.
- 09-803-5F** - Pk. of 100  
Circles; 15cm dia.

### Fisherbrand Qualitative Grade Plain Circles and Sheets - P4 Grade



- 09-803-6A** - Pk. of 100  
Filter Paper, Qualitative; Fisherbrand; Grade P4;  
Cellulose fiber; Particle Retention: 4-8um; Porosity:  
Fine; Filtration Speed: Herzberg, 730 seconds; Flowrate:  
Slow-25mL/min.; Circles; 4.25cm dia.
- 09-803-6B** - Pk. of 100  
Circles; 5.5cm dia.
- 09-803-6C** - Pk. of 100  
Circles; 7cm dia.
- 09-803-6D** - Pk. of 100  
Circles; 9cm dia.
- 09-803-6E** - Pk. of 100  
Circles; 11cm dia.
- 09-803-6F** - Pk. of 100  
Circles; 12.5cm dia.
- 09-803-6G** - Pk. of 100  
Circles; 15cm dia.
- 09-803-6H** - Pk. of 100  
Circles; 18.5cm dia.
- 09-803-6J** - Pk. of 100  
Circles; 20.5cm dia.
- 09-803-6K** - Pk. of 100  
Circles; 24cm dia.

### Fisherbrand Qualitative Grade Plain Circles and Sheets - P5 Grade



- 09-801A** - Pk. of 100  
Filter Paper, Qualitative; Fisherbrand; Grade P5;  
Cellulose fiber; Particle Retention: 5-10um; Porosity:  
Medium; Filtration Speed: Herzberg, 230 seconds;  
Flowrate: Slow-60mL/min.; Circles; 7cm dia.
- 09-801AA** - Pk. of 100  
Circles; 2.5cm dia.
- 09-801B** - Pk. of 100  
Circles; 9cm dia.
- 09-801BB** - Pk. of 100  
Circles; 3.5cm dia.
- 09-801C** - Pk. of 100  
Circles; 11cm dia.
- 09-801D** - Pk. of 100  
P5 Grade; Dia: 12.5cm
- 09-801E** - Pk. of 100  
Circles; 15cm dia.
- 09-801F** - Pk. of 100  
Circles; 19cm dia.
- 09-801G** - Pk. of 100  
Circles; 24cm dia.
- 09-801H** - Pk. of 100  
Circles; 33cm dia.
- 09-801J** - Pk. of 100  
Circles; 4.25cm dia.
- 09-801K** - Pk. of 100  
Circles; 5.5cm dia.
- 09-802** - Pk. of 500  
Sheets; 50 x 50cm
- 09-802-1A** - Pk. of 100  
Sheets; 46 x 57cm

### Fisherbrand Qualitative Grade Plain Circles and Sheets - P8 Grade



- 09-795A** - Pk. of 100  
Filter Paper, Qualitative; Fisherbrand; Grade P8;  
Cellulose fiber; Particle Retention: 20-25um; Porosity:  
Coarse; Filtration Speed: Herzberg, 55 seconds;  
Flowrate: Fast-160mL/min.; Circles; 5.5cm dia.
- 09-795AA** - Pk. of 100  
Circles; 2.5cm dia.
- 09-795B** - Pk. of 100  
Circles; 7cm dia.
- 09-795BB** - Pk. of 100  
Circles; 3.5cm dia.
- 09-795C** - Pk. of 100  
Circles; 9cm dia.
- 09-795D** - Pk. of 100  
Circles; 11cm dia.
- 09-795E** - Pk. of 100  
P8 Grade; Dia: 12.5cm
- 09-795F** - Pk. of 100  
Circles; 15cm dia.
- 09-795G** - Pk. of 100  
Circles; 18.5cm dia.
- 09-804** - Pk. of 100  
Circles; 19cm dia.
- 09-795H** - Pk. of 100  
Circles; 20.5cm dia.
- 09-795J** - Pk. of 100  
Circles; 24cm dia.
- 09-795K** - Pk. of 100  
Circles; 33cm dia.
- 09-802-1B** - Pk. of 100  
Sheets; 46 x 57cm
- 09-800** - Pk. of 500  
Sheets; 50 x 50cm

### Fisherbrand Wet-Strengthened Qualitative Circles



- 09-790-12C** - Pk. of 100  
Filter Paper; Qualitative; Fisherbrand; Wet-Strengthened;  
Circles; Coarse cellulose fiber; Creped surface; Fast  
flowrate; Particle Retention >25um; Diameter: 9cm
- 09-790-12D** - Pk. of 100  
Diameter: 11cm
- 09-790-12E** - Pk. of 100  
Diameter: 12.5cm
- 09-790-12F** - Pk. of 100  
Diameter: 15cm



**09-790-12G** - Pk. of 100  
Diameter: 18.5cm

# Quantitative Filter Paper

## Fisherbrand Quantitative Grade Circles



- 09-790-2A** - Pk. of 100  
Filter Paper, Quantitative; Fisherbrand; Circles; Grade: Q2; Particle retention: 1-5um; Porosity: Fine; Filtration speed: Slow; Surface: Smooth, dense texture; Diameter: 5.5cm
- 09-790-2B** - Pk. of 100  
Grade Q2; Diameter: 7cm
- 09-790-2C** - Pk. of 100  
Grade Q2; Diameter: 9cm
- 09-790-2D** - Pk. of 100  
Grade Q2; Diameter: 11cm
- 09-790-2E** - Pk. of 100  
Grade Q2; Diameter: 12.5cm
- 09-790-2F** - Pk. of 100  
Grade Q2; Diameter: 15cm
- 09-790-2G** - Pk. of 100  
Grade Q2; Diameter: 18.5cm
- 09-790-4C** - Pk. of 100  
Grade Q5; Diameter: 9cm
- 09-790-4D** - Pk. of 100  
Grade Q5; Diameter: 11cm
- 09-790-4E** - Pk. of 100  
Grade Q5; Diameter: 12.5cm
- 09-790-4F** - Pk. of 100  
Grade Q5; Diameter: 15cm
- 09-790A** - Pk. of 100  
Grade Q8; Diameter: 5.5cm
- 09-790B** - Pk. of 100  
Grade Q8; Diameter: 7cm
- 09-790C** - Pk. of 100  
Grade Q8; Diameter: 9cm
- 09-790D** - Pk. of 100  
Grade Q8; Diameter: 11cm
- 09-790E** - Pk. of 100  
Grade Q8; Diameter: 12.5cm
- 09-790F** - Pk. of 100  
Grade Q8; Diameter: 15cm
- 09-790G** - Pk. of 100  
Grade Q8; Diameter: 18.5cm

# Syringe Filters

## Fisherbrand Syringe Filters



- 09-719A** - Pk. of 50  
Filters, Syringe; Fisherbrand; Low hold-up volume; Membrane: MCE; Pore Size: 0.22um; Diameter: 25mm; Sterile
- 09-719B** - Pk. of 50  
Membrane: MCE; Pore Size: 0.45um; Diameter: 25mm; Sterile
- 09-720-5** - Pk. of 100  
Membrane: Nylon; Pore Size: 0.22um; Diameter: 13mm; Nonsterile
- 09-719E** - Pk. of 50  
Membrane: Nylon; Pore Size: 0.22um; Diameter: 25mm; Nonsterile
- 09-719C** - Pk. of 50  
Membrane: Nylon; Pore Size: 0.22um; Diameter: 25mm; Sterile
- 09-720-6** - Pk. of 100  
Membrane: Nylon; Pore Size: 0.45um; Diameter: 13mm; Nonsterile
- 09-719-5** - Pk. of 1000  
Membrane: Nylon; Pore Size: 0.45um; Diameter: 25mm; Nonsterile; Pk/1000
- 09-719F** - Pk. of 50  
Membrane: Nylon; Pore Size: 0.45um; Diameter: 25mm; Nonsterile; Pk/50
- 09-719D** - Pk. of 50  
Membrane: Nylon; Pore Size: 0.45um; Diameter: 25mm; Sterile

- 09-720-7** - Pk. of 100  
Membrane: PTFE; Pore Size: 0.22um; Diameter: 13mm; Nonsterile
- 09-719G** - Pk. of 50  
Membrane: PTFE; Pore Size: 0.22um; Diameter: 25mm; Nonsterile
- 09-719-7** - Pk. of 1000  
Membrane: PTFE; Pore size: 0.45microm; Diameter: 25mm; Nonsterile; Pack of 1000
- 09-720-8** - Pk. of 100  
Membrane: PTFE; Pore Size: 0.45um; Diameter: 13mm; Nonsterile
- 09-719H** - Pk. of 50  
Membrane: PTFE; Pore Size: 0.45um; Diameter: 25mm; Nonsterile
- 09-720-3** - Pk. of 100  
Membrane: PVDF; Pore Size: 0.22um; Diameter: 13mm; Sterile
- 09-720-4** - Pk. of 100  
Membrane: PVDF; Pore Size: 0.45um; Diameter: 13mm; Sterile

## Fisherbrand Wall-Mount Dispenser with 25mm Syringe Filters



- 09-730-16** - Pk. of 150  
Filters; Syringe; in Wall-Mount Dispenser; Fisherbrand; Wall-mount; With 25mm syringe filters; Velcro fastening strips; Reclosable tip-out drawer; Membrane Type: Nylon; Pore Size: 0.2um
- 09-730-17** - Pk. of 150  
Membrane Type: Nylon; Pore Size: 0.45um
- 09-730-19** - Pk. of 150  
Membrane Type: PTFE; Pore Size: 0.2um
- 09-730-21** - Pk. of 150  
Membrane Type: PTFE; Pore Size: 0.45um
- 09-730-25** - Pk. of 150  
Membrane Type: PVDF; Pore Size: 0.22um
- 09-730-27** - Pk. of 150  
Membrane Type: PVDF; Pore Size: 0.45um

# Forceps

## Fisherbrand Dissecting Blunt-Pointed Forceps



- 08-895** - Ea., Cs. of 10  
Forceps, Dissecting Blunt-Pointed; Fisherbrand; Length: 4-1/8 in.; Style: Straight; Case of: 10
- 08-887** - Ea., Cs. of 12  
Length: 5 in.; Style: curved; Case of: 10
- 08-890** - Ea., Cs. of 10  
Length: 5 in.; Style: Straight; Case of: 10

## Fisherbrand General-Purpose Broad-Tipped Forceps



- 10-300** - Ea., Cs. of 10  
Forceps, General-purpose; Fisherbrand; Broad-tipped; Highly polished, ASTM standard stainless steel; Grooved handles; Latex-free; Autoclavable; Passivated; Length: 4-1/2 in. (11.4cm)

## Fisherbrand General-Purpose Curved Forceps



- 10-275** - Ea., Cs. of 12  
Forceps; Fisherbrand; Stainless Steel; Curved; Autoclavable; Length: 5-7/8 in.

## Fisherbrand General-Purpose Extra-Long Forceps



- 10-316A** - Ea., Cs. of 12  
Forceps, General-Purpose Extra-Long; Fisherbrand; Chromium-plated steel; Autoclavable; Length: 7-7/8 in.
- 10-316B** - Ea., Cs. of 12  
Length: 9-7/8 in.
- 10-316C** - Ea., Cs. of 12  
Length: 11-13/16 in.

**Fisherbrand General-Purpose Offset-Tipped Forceps**



**10-295** - Ea., Cs. of 12  
 Forceps, General-purpose; Fisherbrand; Offset-tipped; FDA certified stainless steel; Grooved handles; Latex-free; Length: 4-3/8 in. (11.1cm)

**Fisherbrand General-Purpose Pinning Forceps**



**10-270** - Ea., Cs. of 12  
 Forceps, Pinning; Fisherbrand; Stainless Steel; Extra-Passivated; Style: Curved, Serrated; 4-1/2 in. (114mm)

**Fisherbrand General-Purpose PTFE Forceps**



**02-617-161** - Ea.  
 Forceps, Fisherbrand, General-purpose PTFE, Style: With sharp tips; Length: 3-15/16 in.  
**02-617-162** - Ea.  
 Sharp Tips; Length: 5-7/8 in.  
**02-617-163** - Ea.  
 Sharp Tips; Length: 7-7/8 in.  
**02-617-158** - Ea.  
 Square Tips; Length: 3-15/16 in.  
**02-617-159** - Ea.  
 Square Tips; Length: 5-7/8 in.  
**02-617-160** - Ea.  
 Square Tips; Length: 7-7/8 in.

**Fisherbrand General-Purpose Student Forceps**



**10-280** - Pk. of 10  
 Forceps, Fisherbrand, General-purpose for students, SS, straight, Length: 5 in.

**Fisherbrand Hemostats Halstead Mosquito Forceps**



**13-812-10** - Ea., Cs. of 4  
 Forceps, Mosquito; Fisherbrand; Passivated; Serrated Jaws; Box-lock joint; Curved; Length: 5-1/8 in.  
**13-812-8** - Ea.  
 Straight; Length: 5 in.

**Fisherbrand Hemostats Kelly Forceps**




**13-812-14** - Ea.  
 Forceps, Kelly; Fisherbrand; Serrated jaws; Box-lock joint; Length: 5-1/2 in.; Straight  
**08-907** - Ea.  
 Straight

**Fisherbrand Hemostats Rankin Forceps**



**13-812-45** - Ea., Cs. of 4  
 Forceps, Rankin; Fisherbrand; Stainless Steel; Length: 6-1/4 in.

**Fisherbrand Hemostats Rochester-Ochsner Forceps**



**13-812-16** - Ea.  
 Forceps, Rochester-Ochsner; Fisherbrand; Serrated jaws; Box-lock joint; Length 6-1/4 in.  
**13-812-18** - Ea.  
 Length: 7 1/4 in.

**Fisherbrand Hemostats Rochester-Pean Serrated Forceps**



**13-812-27** - Ea., Cs. of 4  
 Forceps, Rochester-Pean; Fisherbrand; Serrated jaws; Autoclaveable; Curved; Length: 7-1/4 in.  
**13-812-24** - Ea.  
 Straight; Length: 5-1/2 in.  
**13-812-26** - Ea., Cs. of 4  
 Straight; Length: 6-1/4 in.

**Fisherbrand Hemostats Self-Closing Artery Forceps**



**08-905** - Ea., Cs. of 10  
 Forcep, Self-closing; Fisherbrand; Serrated tips; Blunt points; Length: 2-1/4 in. Curved

**13-812-4** - Ea.  
 Straight

# Funnels

## Filling Funnels

**Fisherbrand Powder Funnels**



**10-500** - Pk. of 12, Cs. of 3  
 Funnel; Powder; Fisherbrand; Polypropylene; External ribbing prevents vapor lock; Parallel stem; Autoclavable; Top dia. x H: 65 x 66mm; Capacity: 50mL  
**10-500-1** - Pk. of 12, Cs. of 3  
 Capacity: 100mL  
**10-500-2** - Pk. of 6, Cs. of 4  
 Capacity: 225mL  
**10-500-3** - Pk. of 4, Cs. of 6  
 Capacity: 750mL

## Filtering Funnels

**Fisherbrand Heavy-Duty/Utility Funnels**



**10-500-5** - Pk. of 12, Cs. of 6  
 Funnels; Fisherbrand; Heavy-Duty utility; Polypropylene; External ridge/ribbing prevent vapor lock; Autoclavable; Top dia. x H: 65 x 90mm; Capacity: 40mL  
**10-500-6** - Pk. of 12, Cs. of 6  
 Capacity: 100mL  
**10-500-7** - Pk. of 12, Cs. of 6  
 Capacity: 140mL  
**10-500-8** - Pk. of 12, Cs. of 6  
 Capacity: 210mL  
**10-500-9** - Pk. of 3, Cs. of 8  
 Capacity: 410mL  
**10-500-10** - Pk. of 3, Cs. of 4  
 Capacity: 500mL

**Fisherbrand Long-Stem Analytical Funnels**



**10-500-11** - Pk. of 12, Cs. of 3  
 Funnels, Filtering; Fisherbrand; Polypropylene; External ridge and ribs prevent vapor lock; Ribbed inside slopes 58 deg. ; Autoclavable; Uses standard filter paper; Top dia. x stem L: 34 x 52mm; Height: 75mm; Capacity: 12mL  
**10-500-12** - Pk. of 12, Cs. of 3  
 Capacity: 23mL  
**10-500-13** - Pk. of 12, Cs. of 3  
 Capacity: 37mL  
**10-500-14** - Pk. of 12, Cs. of 3  
 Capacity: 50mL  
**10-500-15** - Pk. of 6, Cs. of 6  
 Capacity: 100mL  
**10-500-16** - Pk. of 6, Cs. of 4  
 Capacity: 150mL  
**10-500-17** - Pk. of 4, Cs. of 6  
 Capacity: 225mL

# Gloves

**Fisherbrand Powdered Vinyl Exam Gloves**



**19-041-189A** - Pk. of 100, Cs. of 10  
 Gloves, Exam; Fisherbrand; Vinyl; Powdered; Nonsterile; Smooth finish; Ambidextrous; Size: X-Small  
**19-041-189B** - Pk. of 100, Cs. of 10  
 Size: Small  
**19-041-189C** - Pk. of 100, Cs. of 10  
 Size: Medium  
**19-041-189D** - Pk. of 100, Cs. of 10  
 Size: Large  
**19-041-189E** - Pk. of 100, Cs. of 10  
 Size: X-Large

### Fisherbrand Powder-Free Latex Gloves with Aloe



**19-050-548A** - Pk. of 100, Cs. of 10  
Glove, Exam; Latex, Powder-Free, Fisherbrand, Small, 9-1/2 in. long, 5-6mil thick, Ambidextrous, Non-Sterile, Natural color, Beaded cuff, Medically approved 510K and meets ASTM standards, Aloe coating inside glove  
**19-050-548B** - Pk. of 100, Cs. of 10  
Medium  
**19-050-548C** - Pk. of 100, Cs. of 10  
Large  
**19-050-548D** - Pk. of 100, Cs. of 10  
X-Large

### Fisherbrand Powder-Free Nitrile Exam Gloves




**19-130-1597B** - Pk. of 100, Cs. of 10  
Gloves, Exam; Nitrile; Powder-free; Fisherbrand; Odor-free; beaded cuff for easier donning; Latex-free; nonsterile; ambidextrous; Length: 9-1/2 in. (240mm); 6mil thickness at fingertip; Small  
**19-130-1597C** - Pk. of 100, Cs. of 10  
Medium  
**19-130-1597D** - Pk. of 100, Cs. of 10  
Large  
**19-130-1597F** - Pk. of 100, Cs. of 10  
XX-Large

### Fisherbrand Powder-Free Nitrile Exam Gloves with Aloe



**19-050-550A** - Pk. of 100, Cs. of 10  
Glove, Exam; Nitrile, Powder-Free, Fisherbrand, Small, 9-1/2 in. long, 5-6mil thick, Ambidextrous, Non-Sterile, Mint Green color, Beaded cuff, Medically approved 510K and meets ASTM standards, Aloe coating inside glove  
**19-050-550B** - Pk. of 100, Cs. of 10  
Medium  
**19-050-550C** - Pk. of 100, Cs. of 10  
Large  
**19-050-550D** - Pk. of 100, Cs. of 10  
X-Large

### Fisherbrand Powder-Free Vinyl Exam Gloves



**19-041-190A** - Pk. of 100, Cs. of 10  
Gloves, Vinyl Exam; Powder-Free; Fisherbrand; Nonsterile; Medical-grade; Smooth finish; Ambidextrous; Exceed FDA minimum requirements; Size: X-Small  
**19-041-190B** - Pk. of 100, Cs. of 10  
Small  
**19-041-190C** - Pk. of 100, Cs. of 10  
Medium  
**19-041-190D** - Pk. of 100, Cs. of 10  
Large  
**19-041-190E** - Pk. of 100, Cs. of 10  
X-Large

## Inoculating Loops and Needles

### Fisherbrand Disposable Needles



**13-070-2** - Pk. of 960, Cs. of 10  
Inoculating Loop; Disposable; Fisherbrand; Needle; Color: Natural; 10 pack

### Fisherbrand Loops



**22-363-608** - Cs. of 250  
Loop, Flexible; Fisherbrand; Inoculating; Disposable; Volume: 1uL; Color: Light Blue; Individual wrap

### Fisherbrand Loops and Needles



**22-363-595** - Pk. of 960, Cs. of 10  
Loop, Flexible; Fisherbrand; Inoculating; Disposable; Volume: 1uL; Color: Light Blue  
**22-363-596** - Pk. of 960, Cs. of 10  
Needle; Color: White  
**22-363-597** - Pk. of 960, Cs. of 10  
Flexible Loop; Volume: 10uL; Color: Yellow  
**22-363-598** - Cs. of 1000  
Rigid Loop; Volume: 1uL; Color: Clear  
**22-363-599** - Cs. of 1000  
Rigid Loop; Volume: 10uL, Light Blue  
**22-363-600** - Cs. of 250  
Flexible Loop; Volume: 10uL; Color: Yellow; Individual wrapped  
**22-363-601** - Cs. of 250  
CombiLoop; Flexible, Volume: 10/1uL; Color: White; Individually wrapped

## Insulated Containers

### Fisherbrand Polyurethane Ice Buckets



**02-591-44** - Ea.  
Bucket, Ice; Fisherbrand; Polyurethane; For use with liquid nitrogen, dry ice, water ice and salt solutions; Color: Black; 4.5L  
**02-591-45** - Ea.  
Color: Blue  
**02-591-46** - Ea.  
Color: Red

## Knives and Knife Blades

### Fisherbrand Stainless-Steel Blades for Dissecting Knife Handles



**08-916-5A** - Pk. of 100, Cs. of 100  
Blades, Knife; Fisherbrand; Surgical quality; Sterile; Individually foil-wrapped; Exceptionally strong and uniformly sharp; Blades for sizes 3, 5, 7 handles; Size: 10  
**08-916-5B** - Pk. of 100, Cs. of 100  
Size: 11  
**08-916-5C** - Pk. of 100, Cs. of 100  
Size: 12  
**08-916-5D** - Pk. of 100, Cs. of 100  
Size: 15

## Handles

### Fisherbrand Disposable-Blade Dissecting Knife Handles



**08-914-5** - Ea.  
Knife Handles, Dissecting; Fisherbrand; Disposable-blade; Re-usable Handles; Slender; Contoured; Length: 6-5/16 in.  
**08-917-5** - Ea.  
Flat with corrugated grip; Length: 5-9/32 in.



## Scalpels

### Fisherbrand Metal Scalpels



**08-920A** - Ea., Cs. of 12  
Scalpel, Dissecting; Fisherbrand; Easy to clean; Retains cutting edge; Cutting Edge: 1 in.  
**08-920B** - Ea., Cs. of 12  
Cutting Edge: 1-1/2 in.

## Marking Pens

### Fisherbrand Marking Pens



**13-379** - Pk. of 10  
Pen, Marking; Fisherbrand; Quick-drying, no unpleasant odor; Writes on cold and damp surfaces; Permanent on most surfaces; Resistant to abrasion, water, and alcohol; Extra Fine tip; Ink Color: Blue; 10/pk  
**13-379-1** - Pk. of 10  
Extra-Fine Tip; Black  
**13-379-2** - Pk. of 10  
Extra-Fine Tip; Red  
**13-379-3** - Pk. of 10  
Extra-Fine Tip; Green  
**13-379-4** - Pk. of 10  
Fine Tip; Black  
**13-379-5** - Pk. of 10  
Fine Tip; Red  
**13-379-6** - Pk. of 10  
Fine Tip; Blue  
**13-379-7** - Pk. of 10  
Fine Tip; Green

### Fisherbrand Securiline Marker II Pens



**14-905-30** - Pk. of 10  
Pens, Marking; Fisherbrand Securiline Marker II; Line width: 0.015 in.; Ink is resistant to many normal laboratory solutions; Black  
**14-905-35** - Pk. of 10  
Red

## Matting, Controlled Environment

### Fisherbrand Hard Surface Adhesive Mats




**19-181-516** - Cs. of 4  
Mat; Hard Surface Adhesive; Fisherbrand; Reorder Notification on Eighteenth Sheet; Blue Sheets on Blue Mat; High Tack; Nominal L x W: 36 x 18 in. (91 x 46cm); High Tack  
**19-181-517** - Cs. of 4  
High Tack; Nominal L x W: 46 x 18 in. (117 x 46cm)  
**19-181-518** - Cs. of 4  
Moderate Tack; Nominal L x W: 46 x 18 in. (117 x 46cm)  
**19-181-519** - Cs. of 4  
Moderate Tack; L x W: 36 x 18 in. (91 x 46cm)  
**19-181-520** - Cs. of 4  
Moderate Tack; L x W: 46 x 24 in. (117 x 61cm)

## Microplates


### Assay Microplates

#### Fisherbrand 96-Well Plates



**21-377-203** - Cs. of 50  
Plate, 96-Well; Fisherbrand; Microwell; Nonsterile; 0.45mL; Certified; Clear; polystyrene; flat bottom; bulk packs of 50 plates  
**21-377-204** - Cs. of 50  
Nonsterile; Not Certified  
**21-377-205** - Cs. of 40  
Sterile; Certified

#### Fisherbrand Clear Polystyrene 96-Well Plates



**12-565-500** - Cs. of 80  
Microplate, Assay, 96-well; Fisherbrand; Clear; Round well; 250uL well volume; Untreated; Nonsterile; Polystyrene; Without lid  
**12-565-501** - Cs. of 80  
Flat well

#### Fisherbrand Clear Polystyrene 384-Well Plates



**12-565-506** - Cs. of 80  
Plate, Polystyrene; Fisherbrand; 384-well; Clear; square; 10-100uL working volume; nonsterile; untreated; Ideal for low-volume assays such as colorimetric, enzymatic, cloning, replication, array, and storage  
**12-565-511** - Cs. of 180  
Lid

## Genomic Microplates

#### Fisherbrand 96-Well Low-Profile, Skirted PCR Plate




**14-230-237** - Cs. of 25  
Microplate, PCR, 96-Well, Low-Profile; Fisherbrand; Well Capacity: 200uL; Raised rim; Virgin polypropylene; Skirted; Color: Natural; Cut-away corner; Alphanumeric grid; Cleanroom produced  
**14-230-238** - Cs. of 25  
Red  
**14-230-239** - Cs. of 25  
Blue  
**14-230-240** - Cs. of 25  
Green  
**14-230-241** - Cs. of 25  
Yellow

#### Fisherbrand 96-Well Nonskirted PCR Plates



**14-230-232** - Cs. of 25  
Microplate, PCR, 96-well; Fisherbrand; Well Capacity: 300uL; Raised rim; Virgin polypropylene; Nonskirted; Color: Natural; Cut-away corner; Alphanumeric grid; Cleanroom produced  
**14-230-233** - Cs. of 25  
Red  
**14-230-234** - Cs. of 25  
Blue  
**14-230-235** - Cs. of 25  
Green  
**14-230-236** - Cs. of 25  
Yellow

#### Fisherbrand 96-Well Semi-Skirted PCR Plates



**14-230-244** - Cs. of 25  
Microplate, PCR, 96-Well; Fisherbrand; Well Capacity: 300uL; Raised rim; Virgin polypropylene; Semi-skirted; Color: Natural; Cut-away corner; Alphanumeric grid; Cleanroom produced



- 14-230-245** - Cs. of 25  
Red
- 14-230-246** - Cs. of 25  
Blue
- 14-230-247** - Cs. of 25  
Green
- 14-230-248** - Cs. of 25  
Yellow

### Fisherbrand 96-Well Tall-Chimney PCR Plate



**14-230-242** - Cs. of 50  
Microplate, PCR, 96-well, Tall-Chimney; Fisherbrand; Well Capacity: 340uL; Raised well rim; Virgin polypropylene; Nonskirted; Color: Natural; Cut-away corner; Alphanumeric grid; Cleanroom produced

### Fisherbrand 384-Well Skirted PCR Plate



**14-230-243** - Cs. of 50  
Microplate, PCR, 384-Well; Fisherbrand; Max. Well Capacity: 40uL; Raised rim; Virgin polypropylene; Skirted; Color: Natural; Cut-away corners; Alphanumeric grid; Cleanroom produced

## PCR Tubes

### Fisherbrand 0.2mL Flat- and Domed-Cap PCR Tubes



- 14-230-205** - Cs. of 1000  
Tube, PCR; Fisherbrand; Domed cap; Natural; 0.2mL; Virgin polypropylene; Thin-wall design; Cleanroom produced
- 14-230-206** - Cs. of 1000  
Domed cap; Red
- 14-230-207** - Cs. of 1000  
Domed cap; Blue
- 14-230-208** - Cs. of 1000  
Domed cap; Green
- 14-230-209** - Cs. of 1000  
Domed cap; Yellow
- 14-230-225** - Cs. of 1000  
Flat cap; Natural
- 14-230-226** - Cs. of 1000  
Flat cap; Red
- 14-230-227** - Cs. of 1000  
Flat cap; Blue
- 14-230-228** - Cs. of 1000  
Flat cap; Green
- 14-230-229** - Cs. of 1000  
Flat cap; Yellow

### Fisherbrand 0.2mL PCR Tube Strips



- 14-230-210** - Cs. of 250  
Tube Strips, PCR; Fisherbrand; 8-tube strips; Domed cap strips; Natural; 0.2mL; Virgin polypropylene; Thin-wall design; Cleanroom produced
- 14-230-211** - Cs. of 250  
Domed cap strips; Red
- 14-230-212** - Cs. of 250  
Domed cap strips; Blue
- 14-230-213** - Cs. of 250  
Domed cap strips; Green
- 14-230-214** - Cs. of 250  
Domed cap strips; Yellow
- 14-230-215** - Cs. of 250  
Flat cap strips; Natural
- 14-230-216** - Cs. of 250  
Flat cap strips; Red
- 14-230-217** - Cs. of 250  
Flat cap strips; Blue
- 14-230-218** - Cs. of 250  
Flat cap strips; Green
- 14-230-219** - Cs. of 250  
Flat cap strips; Yellow
- 14-230-220** - Cs. of 250  
Without cap strips; Natural
- 14-230-221** - Cs. of 250  
Without cap strips; Red

- 14-230-222** - Cs. of 250  
Without cap strips; Blue
- 14-230-223** - Cs. of 250  
Without cap strips; Green
- 14-230-224** - Cs. of 250  
Without cap strips; Yellow


### Fisherbrand 0.5mL Flat-Cap PCR Tubes



- 14-230-200** - Cs. of 1000  
Tube, PCR; Fisherbrand; Flat cap; Natural; 0.5mL; Virgin polypropylene; Thin-wall design; Cleanroom produced
- 14-230-201** - Cs. of 1000  
Red
- 14-230-202** - Cs. of 1000  
Blue
- 14-230-203** - Cs. of 1000  
Green
- 14-230-204** - Cs. of 1000  
Yellow

## PCR Caps

### Fisherbrand Flat and Domed 8-Cap Strips



- 14-230-230** - Cs. of 250  
Cap Strip; Fisherbrand; For use with 96-well PCR plates or PCR tube strips; Flat caps; 8 caps per strip; Virgin polypropylene; Cleanroom produced; Large end tabs for easy removal; Labeled A-H for strip orientation
- 14-230-231** - Cs. of 250  
Domed Cap Strips

## Microplate Sealing Supplies

### Fisherbrand MicroMat Cap Mat



- 08-408-219** - Pk. of 5  
Microplate Sealing Supplies; Cap mat; MicroMat; For Use with: 96-well plates; Fisherbrand; For pressure sealing thermal cyclers

### Fisherbrand Silicone Plate Seals



- 03-396-49** - Pk. of 10  
Seal, Microplate; Fisherbrand; Silicone; Clear; Lightweight; For Costar 96-Well Polypropylene Plates and Assay Blocks, Whatman 96-Well Multi-Chem Microplates; Nonslit
- 03-396-50** - Pk. of 10  
For Costar 96-Well Polypropylene Plates and Assay Blocks, Whatman 96-Well Multi-Chem Microplates; Preslit
- 03-396-51** - Pk. of 10  
For Fisherbrand Round-Bottom 96-Well Polypropylene Plates and Assay Blocks, Nunc 96-Well Storage Plates; Nonslit
- 03-396-52** - Pk. of 10  
For Fisherbrand Round-Bottom 96-Well Polypropylene Plates and Assay Blocks, Nunc 96-Well Storage Plates; Preslit
- 03-396-53** - Pk. of 10  
For Fisherbrand 96-Well 2.2mL Assay Blocks, Whatman 96-Well Polypropylene Storage Plates (Square-Well Plate); Nonslit
- 03-396-54** - Pk. of 10  
For Fisherbrand 96-Well 2.2mL Assay Blocks, Whatman 96-Well Polypropylene Storage Plates (Square-Well Plate); Preslit
- 03-396-55** - Pk. of 10  
For Fisherbrand 384-Well Polypropylene Plates, Costar Brand 384-Well Polypropylene Storage Block (Square-Well Plate); Nonslit

**03-396-56** - Pk. of 10  
 For Fisherbrand 384-Well Polypropylene Plates, Costar Brand 384-Well Polypropylene Storage Block (Square-Well Plate); Preslit

# Microscope Coverslips

## Fisherbrand Cover Glass

- 22-266-858** - Pk. of 1, Cs. of 10  
 Fisherbrand Coverglass; No. 1.5; 24 x 30mm; 1 oz. pack
- 22-266-882** - Pk. of 1, Cs. of 10  
 No. 1.5; 24 x 60mm
- 22-266-908** - Pk. of 1, Cs. of 10  
 No. 2; 22 x 30mm
- 22-266-940** - Pk. of 1, Cs. of 10  
 No. 2; 24 x 30mm
- 22-266-973** - Pk. of 1, Cs. of 10  
 No. 2; 24 x 60mm
- 22-031-145** - Pk. of 152, Cs. of 10  
 No. 1; 15mm dia. circle
- 22-293-232** - Pk. of 1, Cs. of 10  
 No. 2; 12mm dia. circle

## Fisherbrand Coverslips

- 22-037-169** - Pk. of 1, Cs. of 10  
 Fisherbrand Coverslip; 24 x 40mm; 1mm thickness; Pack of 1
- 22-037-170** - Pk. of 1, Cs. of 10  
 24 x 50mm
- 22-037-171** - Pk. of 1, Cs. of 10  
 24 x 60mm

## Fisherfinest Superslip Cover Glass

- 12-545-89** - Pk. of 1, Cs. of 10  
 Coverslip; Fisherfinest Superslip; 24 x 60mm; 1 oz. pack

# Microscope Slide Holders

## Fisherbrand Color-Coded ABS Plastic Boxes for 25 or 100 Slides



- 03-448-6** - Ea., Cs. of 72  
 Box, Microscope Slide; Fisherbrand; Color-coded; ABS plastic; Hold 3 x 1 in. and 75 x 25mm slides; With hinged lid; Numbered slide partitions; Stackable; Capacity: Hold 25 Slides; Color: Blue
- 03-448-7** - Ea., Cs. of 72  
 Capacity: 25 Slides; Color: Green
- 03-448-8** - Ea., Cs. of 72  
 Capacity: 25 Slides; Color: Red
- 03-448-9** - Ea., Cs. of 72  
 Capacity: 25 Slides; Color: Yellow
- 03-448-10** - Ea., Cs. of 72  
 Capacity: 25 Slides; Color: White
- 03-448-1** - Ea., Cs. of 72  
 Capacity: 100 Slides; Color: Blue
- 03-448-2** - Ea., Cs. of 72  
 Capacity: 100 Slides; Color: Green
- 03-448-3** - Ea., Cs. of 72  
 Capacity: 100 Slides; Color: Red
- 03-448-4** - Ea., Cs. of 72  
 Capacity: 100 Slides; Color: Yellow
- 03-448-5** - Ea., Cs. of 72  
 Capacity: 100 Slides; Color: White

# Microscope Slides

## Fisherbrand ColorFrost Microscope Slides with Clipped Corners

- 22-037-160** - Pk. of 144, Cs. of 10  
 Slides, Microscope, ColorFrost; Fisherbrand; 45-degree clipped corners; Color-coded frosted end; Chemical resistant; Marks easily; 25 x 75mm; Thickness: 1mm; Blue
- 22-037-161** - Pk. of 144, Cs. of 10  
 Green
- 22-037-162** - Pk. of 144, Cs. of 10  
 Pink
- 22-037-163** - Pk. of 144, Cs. of 10  
 Yellow

## Fisherbrand Extra-thick Microscope Slides

- 22-267-062** - Pk. of 144, Cs. of 10  
 Slides, Microscope; Fisherbrand; Extra-thick; 1.2mm, clear, corrosion-resistant glass; Precleaned; Flat surfaces and ground edges; Size: 3 x 1 in.; Sides frosted: Double
- 22-267-096** - Pk. of 144, Cs. of 10  
 Size: 3 x 1 in.; Type: plain glass
- 22-267-005** - Pk. of 144, Cs. of 5  
 Size: 3 x 1.5 in.; Type: plain glass
- 22-267-013** - Pk. of 144, Cs. of 5  
 Size: 3 x 2 in.; Type: plain glass

## Fisherbrand Microscope Slides - Frosted

- 22-265-446** - Pk. of 144, Cs. of 25  
 Fisherbrand Microscope Slide; 3 x 1 in.; Frosted; Pack of 144
- 12-550-33** - Pk. of 144, Cs. of 10  
 3 x 1 in.; Frosted area
- 12-550-343** - Pk. of 144, Cs. of 10  
 3 x 1 in.; Frosted at one End

## Fisherbrand Microscope Slides - Plain

- 12-550-A3** - Pk. of 144, Cs. of 10  
 Slides, Microscope; Fisherbrand; Plain; 75x25mm; Pack of 144

## Fisherbrand Superfrost Excell Microscope Slides



- 22-034-985** - Pk. of 144, Cs. of 10  
 Slides; Microscopy; Fisherbrand; Superfrost; Excell; Superior adhesion and wettability; For immunochemistry laboratory applications; Adhesion surface for heat induced epitope retrieval; Uniform hydrophilic surface

## Fisherbrand Superfrost Plus and Colorfrost Plus Microscope Slides



- 22-230-890** - Pk. of 144, Cs. of 10  
 Slides; Microscopy; Fisherbrand; ColorFrost; Special treatment to minimize tissue loss while staining; Ideal for frozen sections, cyto centrifuge preparations, standard Papanicolaou smears; Color: aqua; Size: 75 x 25mm
- 12-550-17** - Gross, Cs. of 10  
 Colorfrost Plus; Blue
- 22-230-892** - Pk. of 144, Cs. of 10  
 Colorfrost Plus; Burnt Orange
- 22-230-891** - Pk. of 144, Cs. of 10  
 Colorfrost Plus; Gold
- 12-550-18** - Gross, Cs. of 10  
 Colorfrost Plus; Green
- 12-550-16** - Pk. of 144, Cs. of 10  
 Colorfrost Plus; Lilac
- 12-550-19** - Gross, Cs. of 10  
 Colorfrost Plus; Pink
- 12-550-16A** - Pk. of 144, Cs. of 10  
 Colorfrost Plus; Tan
- 12-550-20** - Gross, Cs. of 10  
 Colorfrost Plus; Yellow
- 12-550-15** - Gross, Cs. of 10  
 Superfrost Plus; White

## Fisherfinest Premium Plain Glass Microscope Slides




**12-544-1** - Gross, Cs. of 10  
Slides, Microscopy; Premium; Fisherfinest; Plain glass;  
Clear; Premium quality, water-white glass; Precleaned;  
Size: 3 x 1 in.

**12-544-4** - Gross, Cs. of 10  
Size: 75 x 25mm

# Petri Dishes

## Fisherbrand Compartmentalized Disposable Petri Dishes



**08-757-150** - Cs. of 500  
Petri Dish, Disposable; Sectioned plates; W/numbered  
sections; Radiation sterilized; Polystyrene; 25/  
polyethylene sleeve; Size: 100 dia. x 15mm H; Type: I (2  
sections)

**08-757-90** - Cs. of 500  
Beveled Ridge; 20 per sleeve

**08-757-91** - Cs. of 500  
Raised Ridge; 20 per sleeve

**08-757-9A** - Cs. of 500  
Beveled Ridge

**08-758-2** - Cs. of 500  
Size: 100 dia. x 15mm H; Type: X (4 sections)

**08-758-3** - Cs. of 500  
Size: 100 dia. x 15mm H; Type: Y (3 sections)

## Fisherbrand Disposable Contact Plate



**08-757-11C** - Cs. of 500  
Plates, Contact; Disposable; Fisherbrand; 65 x 15mm  
Sterile; Sixteen 10 x 10mm square grid; Marked  
numerically and alphabetically

## Fisherbrand Petri Dish

**08-757-9B** - Cs. of 500  
Dish, Petri; Fisherbrand; Disposable; Stackable; Vented; Radiation sterilized;  
Raised Ridge; 100 x 15mm

## Fisherbrand Petri Dishes with Clear Lids



**08-757-12** - Cs. of 500  
Dish, Petri; Fisherbrand; Clear Lid; Radiation-sterilized;  
Polystyrene; Raised Ridge; 100 Dia. x 15mmH; Edges  
are square; Three venting ribs in lid; Stackable;  
Disposable

**08-757-13** - Cs. of 500  
Beveled Ridge; 100 O.D. x 15mm H

**08-757-14** - Cs. of 100  
Beveled Ridge; 150 Dia. x 15mmH

# Pipet Tips

## Specialty Pipetter Tips

### Fisherbrand Dispenser Pipet Tips



**13-683-61H** - Pk. of 100  
Pipet Tips, Dispenser; Fisherbrand; Polypropylene; inner  
plunger of stainless steel; designed to fit Eppendorf  
Repeater Pipetter; Dispense from 8 to 48 times without  
refilling; +/-2.5% accuracy; long, thin; Nonsterile; 0.05mL

**13-683-61A** - Pk. of 100  
Nonsterile; 0.5mL

**13-683-61B** - Pk. of 100  
Nonsterile; 1.25mL

**13-683-61C** - Pk. of 100  
Nonsterile; 2.5mL

**13-683-61D** - Pk. of 100  
Nonsterile; 5.0mL

**13-683-61E** - Pk. of 100  
Nonsterile; 12.5mL

**13-683-61F** - Pk. of 100  
Nonsterile; 25mL

**13-683-61G** - Pk. of 100  
Nonsterile; 50mL

**13-683-62A** - Pk. of 100  
Sterile; 0.5mL

**13-683-62C** - Pk. of 100  
Sterile; 2.5mL

**13-683-62D** - Pk. of 100  
Sterile; 5.0mL

**13-683-62F** - Pk. of 100  
Sterile; 50mL

**13-683-62E** - Pk. of 100  
Sterile; 12.5mL

**13-683-62G** - Pk. of 100  
Sterile; 25mL

## Fisherbrand Gel-Loading Tips



**02-707-88** - Pk. of 1020, Cs. of 10  
Tips, Pipet; Gel-loading; Fisherbrand; ideal for DNA  
sequencing or loading electrophoresis gels; Round tip;  
Natural color; Volume, length: 0.5-10uL, 2.3 in. ; O.D.:  
0.25mm; Packaging: LR, 5 x 204

**02-707-89** - Pk. of 1020, Cs. of 10  
Volume, length: 0.5-10uL, 2.3 in. ; O.D.: 0.25mm;  
Packaging: LRS, 5 x 204

**02-707-86** - Pk. of 1632, Cs. of 3  
Volume, length: 0.5-10uL, 2.3 in. ; O.D.: 0.6mm;  
Packaging: LR, 8 x 204

**02-707-182** - Pk. of 1020, Cs. of 5  
Volume, length: 1-200uL, 2.75 in. ; O.D.: 0.6mm ;  
Packaging: LR, 5 x 204

**02-707-81** - Pk. of 1020, Cs. of 5  
Volume, length: 1-200uL, 2.75 in. ; O.D.: 0.6mm ;  
Packaging: LRS, 5 x 204

**02-707-181** - Pk. of 1000, Cs. of 10  
Volume, length: 1-200uL, 2.75 in. ; O.D.: 0.6mm;  
Packaging: BP

**02-707-138** - Pk. of 960, Cs. of 5  
Volume, length: 1-200uL, 2.75 in. ; O.D.: 0.6mm;  
Packaging: HR, 10 x 96

**02-707-139** - Pk. of 960, Cs. of 5  
Volume, length: 1-200uL, 2.75 in. ; O.D.: 0.6mm;  
Packaging: HRS, 10 x 96

**02-707-82** - Pk. of 1000, Cs. of 10  
Volume, length: 1uL, 2.25 in. ; O.D.: 0.6mm; Packaging:  
BP

**02-707-83** - Pk. of 1020, Cs. of 5  
Volume, length: 1uL, 2.25 in. ; O.D.: 0.6mm; Packaging:  
LR, 5 x 204

**02-707-84** - Pk. of 1020, Cs. of 5  
Volume, length: 1uL, 2.25 in. ; O.D.: 0.6mm; Packaging:  
LRS, 5 x 204

## Fisherbrand Micro Pipet Tips



**02-707-146** - Pk. of 1000, Cs. of 10  
Pipet Tips; Fisherbrand Micro, Reach; Volume: 0.1-10uL,  
1-1/4 in. Nano; BP; Sterile; 1000/Pk.

**02-707-147** - Pk. of 960, Cs. of 5  
Volume: 0.1-10uL, 1-1/4 in. Nano; LR, 10 x 96; Sterile;  
960/Pk.

**21-277-2A** - Pk. of 1000, Cs. of 10  
Volume: 0.1-10uL, 1-1/4 in.; Sterile; 1000/Pk.

**21-277-2B** - Pk. of 960, Cs. of 5  
Volume: 0.1-10uL, 1-1/4 in.; HR, 10 x 96; Sterile; 1000/  
Pk.

**02-681-440** - Pk. of 960, Cs. of 5  
Volume: 0.1-10uL, 1-1/4 in.; HRS, 10 x 96; Sterile; 960/  
Pk.

**02-707-135** - Pk. of 480, Cs. of 10  
Volume: 0.1-10uL, 1-1/4 in.; RI, 5 x 96; Sterile; 480/Pk.

**02-707-149** - Pk. of 1000, Cs. of 10  
Volume: 0.1-10uL, 1-11/16 in. Reach; BP; Sterile; 1000/  
Pk.

**02-707-136** - Pk. of 960, Cs. of 5  
Volume: 0.1-10uL, 1-11/16 in. Reach; HR, 10 x 96;  
Sterile; 960/Pk.

**02-681-437** - Pk. of 960, Cs. of 5  
Volume: 0.1-10uL, 1-11/16 in. Reach; HRS, 10 x 96;  
Sterile; 960/Pk.

**02-707-137** - Pk. of 960, Cs. of 5  
Volume: 0.1-10uL; Sterile; Hinged Rack; 10 x 96



**21-197-2F** - Pk. of 1000, Cs. of 10  
Volume: 0.5-10uL, 1-13/16 in.; BP; Sterile; 1000/Pk.  
**21-197-2E** - Pk. of 960, Cs. of 5  
Volume: 0.5-10uL, 1-13/16 in.; HR, 10 x 96; Sterile; 960/  
Pk.

**Fisherbrand Premium Plus Gel Loading Flex Tips**



**05-408-152** - Pk. of 200, Cs. of 2  
Pipetter tips; Specialty; Gel-loading; Fisherbrand;  
Premium Plus; MiniFlex; Tip style: flat; Outer diameter:  
0.17mm; Packaging: 200/pk.; Fit ultra-micro pipettors;  
45mm long; 2uL graduation; 0.1 to 10uL volume  
**05-408-154** - Pk. of 1000, Cs. of 10  
MiniFlex; Tip style: round; Outer diameter: 0.57mm; 200/  
bag  
**05-408-153** - Pk. of 200, Cs. of 2  
MiniFlex; Tip style: round; Outer diameter: 0.57mm; 200/  
pk.  
**05-408-150** - Pk. of 200, Cs. of 2  
Multiflex; Tip style: flat; Outer diameter: 0.37mm; 200/pk  
**05-408-151** - Pk. of 200, Cs. of 2  
Multiflex; Tip style: round; Outer diameter: 0.57mm; 200/  
pk.

**Fisherbrand SureOne Pipet Tips: Nonfilter**



**02-707-449** - Pk. of 1000, Cs. of 10  
Tips, Pipet; Fisherbrand; SureOne; Nonfilter; Beveled Tip;  
Graduated at 10, 50, 100uL; Volume: 1-200uL Length: 2  
in.; Natural; BP  
**02-707-415** - Pk. of 960, Cs. of 5  
Beveled Tip; Graduated at 10, 50, 100uL; Volume:  
1-200uL Length: 2 in.; Natural; HR, 10 x 96  
**02-707-416** - Pk. of 960, Cs. of 5  
Beveled Tip; Graduated at 10, 50, 100uL; Volume:  
1-200uL Length: 2 in.; Natural; HRS, 10 x 96  
**02-707-417** - Pk. of 960, Cs. of 5  
Beveled Tip; Graduated at 10, 50, 100uL; Volume:  
1-200uL Length: 2 in.; Natural; RD, 10 x 96  
**02-707-448** - Pk. of 1000, Cs. of 10  
Beveled Tip; Graduated at 10, 50, 100uL; Volume:  
1-200uL Length: 2 in.; Yellow; BP  
**02-707-412** - Pk. of 960, Cs. of 5  
Beveled Tip; Graduated at 10, 50, 100uL; Volume:  
1-200uL Length: 2 in.; Yellow; HR, 10 x 96  
**02-707-413** - Pk. of 960, Cs. of 5  
Beveled Tip; Graduated at 10, 50, 100uL; Volume:  
1-200uL Length: 2 in.; Yellow; HRS, 10 x 96  
**02-707-414** - Pk. of 960, Cs. of 5  
Beveled Tip; Graduated at 10, 50, 100uL; Volume:  
1-200uL Length: 2 in.; Yellow; RD, 10 x 96  
**02-707-450** - Pk. of 1000, Cs. of 10  
Beveled Tip; Volume: 1-200uL; Length: 1-1/4 in.; Natural;  
BP  
**02-707-451** - Pk. of 1000, Cs. of 10  
Beveled Tip; Volume: 1-200uL; Length: 1-1/4 in.; Yellow;  
BP  
**02-707-418** - Pk. of 960, Cs. of 5  
Beveled Tip; Volume: 1-200uL; Length: 2 in.; Natural; HR,  
10 x 96  
**02-707-419** - Pk. of 960, Cs. of 5  
Beveled Tip; Volume: 1-200uL; Length: 2 in.; Natural;  
HRS, 10 x 96  
**02-707-420** - Pk. of 960, Cs. of 5  
Beveled Tip; Volume: 1-200uL; Length: 2 in.; Natural; RD,  
10 x 96  
**02-707-421** - Pk. of 960, Cs. of 5  
Beveled Tip; Volume: 1-200uL; Length: 2 in.; Yellow; HR,  
10 x 96  
**02-707-422** - Pk. of 960, Cs. of 5  
Beveled Tip; Volume: 1-200uL; Length: 2 in.; Yellow;  
HRS, 10 x 96  
**02-707-423** - Pk. of 960, Cs. of 5  
Beveled Tip; Volume: 1-200uL; Length: 2 in.; Yellow; RD,  
10 x 96  
**02-707-447** - Pk. of 1000, Cs. of 10  
Beveled Tip; Volume: 5-300uL; Length: 2 in.; Natural; BP  
**02-707-409** - Pk. of 960, Cs. of 5  
Beveled Tip; Volume: 5-300uL; Length: 2 in.; Natural; HR,  
10 x 96  
**02-707-410** - Pk. of 960, Cs. of 5  
Beveled Tip; Volume: 5-300uL; Length: 2 in.; Natural;  
HRS, 10 x 96

**02-707-411** - Pk. of 960, Cs. of 5  
Beveled Tip; Volume: 5-300uL; Length: 2 in.; Natural; RD,  
10 x 96  
**02-707-453** - Pk. of 1000, Cs. of 10  
Micropoint Tip; Eppendorf Style; Volume: 0.5-10uL;  
Length: 1-1/4 in.; Natural; BP  
**02-707-433** - Pk. of 960, Cs. of 5  
Micropoint Tip; Eppendorf Style; Volume: 0.5-10uL;  
Length: 1-1/4 in.; Natural; HR, 10 x 96  
**02-707-434** - Pk. of 960, Cs. of 5  
Micropoint Tip; Eppendorf Style; Volume: 0.5-10uL;  
Length: 1-1/4 in.; Natural; HRS, 10 x 96  
**02-707-405** - Pk. of 1000, Cs. of 10  
Micropoint Tip; Graduated at 100, 200, 500, 1000uL;  
Volume: 100-1250uL; Length: 2-3/4 in.; Blue; BP  
**02-707-406** - Pk. of 960, Cs. of 5  
Micropoint Tip; Graduated at 100, 200, 500, 1000uL;  
Volume: 100-1250uL; Length: 2-3/4 in.; Blue; HR, 10 x 96  
**02-707-407** - Pk. of 960, Cs. of 5  
Micropoint Tip; Graduated at 100, 200, 500, 1000uL;  
Volume: 100-1250uL; Length: 2-3/4 in.; Blue; HRS, 10 x  
96  
**02-707-408** - Pk. of 960, Cs. of 5  
Micropoint Tip; Graduated at 100, 200, 500, 1000uL;  
Volume: 100-1250uL; Length: 2-3/4 in.; Blue; RI, 10 x 96  
**02-707-400** - Pk. of 1000, Cs. of 10  
Micropoint Tip; Graduated at 100, 200, 500, 1000uL;  
Volume: 100-1250uL; Length: 2-3/4 in.; Natural; BP  
**02-707-401** - Pk. of 960, Cs. of 5  
Micropoint Tip; Graduated at 100, 200, 500, 1000uL;  
Volume: 100-1250uL; Length: 2-3/4 in.; Natural; HR, 10 x  
96  
**02-707-402** - Pk. of 960, Cs. of 5  
Micropoint Tip; Graduated at 100, 200, 500, 1000uL;  
Volume: 100-1250uL; Length: 2-3/4 in.; Natural; HRS, 10  
x 96  
**02-707-403** - Pk. of 960, Cs. of 5  
Micropoint Tip; Graduated at 100, 200, 500, 1000uL;  
Volume: 100-1250uL; Length: 2-3/4 in.; Natural; RI, 10 x  
96  
**02-707-455** - Pk. of 1000, Cs. of 10  
Micropoint Tip; Graduated at 2.5uL; Extended Length;  
Volume: 0.1-10uL; Length: 1-1/4 in.; Natural; BP  
**02-707-440** - Pk. of 960, Cs. of 5  
Micropoint Tip; Graduated at 2.5uL; Extended Length;  
Volume: 0.1-10uL; Length: 1-1/4 in.; Natural; HR, 10 x 96  
**02-707-441** - Pk. of 960, Cs. of 5  
Micropoint Tip; Graduated at 2.5uL; Extended Length;  
Volume: 0.1-10uL; Length: 1-1/4 in.; Natural; HRS, 10 x  
96  
**02-707-436** - Pk. of 960, Cs. of 5  
Micropoint Tip; Graduated at 2uL; Volume: 0.1-10uL;  
Length: 1-1/4 in.; Natural; HR, 10 x 96  
**02-707-437** - Pk. of 960, Cs. of 5  
Micropoint Tip; Graduated at 2uL; Volume: 0.1-10uL;  
Length: 1-1/4 in.; Natural; HRS, 10 x 96  
**02-707-438** - Pk. of 960, Cs. of 5  
Micropoint Tip; Graduated at 2uL; Volume: 0.1-10uL;  
Length: 1-1/4 in.; Natural; RD, 10 x 96  
**02-707-454** - Pk. of 1000, Cs. of 10  
Micropoint Tip; Graduated at 2uL; Volume: 0.1-10uL;  
Length: 2 in.; Natural; BP  
**02-707-446** - Pk. of 1000, Cs. of 10  
Thin Wall Micropoint Tip; Volume: 1-200uL; Length: 2 in.;  
Natural; BP  
**02-707-424** - Pk. of 960, Cs. of 5  
Thin Wall Micropoint Tip; Volume: 1-200uL; Length: 2 in.;  
Natural; HR, 10 x 96  
**02-707-425** - Pk. of 960, Cs. of 5  
Thin Wall Micropoint Tip; Volume: 1-200uL; Length: 2 in.;  
Natural; HRS, 10 x 96  
**02-707-426** - Pk. of 960, Cs. of 5  
Thin Wall Micropoint Tip; Volume: 1-200uL; Length: 2 in.;  
Natural; RD, 10 x 96  
**02-707-452** - Pk. of 1000, Cs. of 5  
Thin Wall Micropoint Tip; Volume: 1-200uL; Length: 2 in.;  
Yellow; BP  
**02-707-427** - Pk. of 960, Cs. of 5  
Thin Wall Micropoint Tip; Volume: 1-200uL; Length: 2 in.;  
Yellow; HR, 10 x 96  
**02-707-428** - Pk. of 960, Cs. of 5  
Thin Wall Micropoint Tip; Volume: 1-200uL; Length: 2 in.;  
Yellow; HRS, 10 x 96

**02-707-429** - Pk. of 960, Cs. of 5  
Thin Wall Micropoint Tip; Volume: 1-200uL; Length: 2 in.;  
Yellow; RD, 10 x 96

## Universal Pipetter Tips

### Fisherbrand Redi-Tip Pipet Tips



**02-707-502** - Pk. of 960, Cs. of 5  
Pipet Tips; Redi-Tip; General Purpose; Volume: 1-200uL;  
Yellow; Sterile; Length 2 in.; Hinged rack; 10 x 96

**02-707-505** - Pk. of 960, Cs. of 5  
1-200uL; Nonsterile

**02-707-509** - Pk. of 1000, Cs. of 5  
201-1000uL; Sterile

**02-707-510** - Pk. of 1000, Cs. of 5  
201-1000uL; Nonsterile

**02-681-436** - Cs. of 96  
Empty Rack for 200uL tips

**02-681-438** - Ea., Cs. of 50  
Empty Rack for 1000uL tips

**02-707-500** - Pk. of 1000, Cs. of 10  
General Purpose; 1-200uL; Yellow; Bulk Pack

**02-707-119** - Pk. of 960, Cs. of 5  
General Purpose; Volume: 1-300uL; Lift-off rack, 10 x 96;  
Natural; Sterile

**02-707-508** - Pk. of 1000, Cs. of 5  
General-Purpose; Volume, length: 101-1000uL, 2-3/4 in.;  
Color: blue; Packaging: HR, 10 x 100; 1000/PK

**02-707-507** - Pk. of 1000, Cs. of 10  
General-Purpose; Volume/Length: 101-1000uL, 2-3/4 in.;  
Color: blue; Packaging: Bulk

**02-681-157** - Pk. of 1000, Cs. of 10  
Reference; 0 to 300uL; Natural; BP

**02-681-158** - Pk. of 960, Cs. of 5  
Reference; 0-300uL; Natural; LR, 10 x 96

**02-681-454** - Pk. of 1000, Cs. of 10  
Reference; 1 to 300uL; Natural; Bulk pack

**02-681-455** - Pk. of 960, Cs. of 5  
Reference; 1 to 300uL; Natural; Lift-off rack, 10 x 96

**02-681-164** - Pk. of 1000, Cs. of 5  
Reference; 101-1000uL; Blue; Lift-off rack

**02-681-166** - Pk. of 1000, Cs. of 5  
Reference; 101-1000uL; Blue; Stack rack, 2 x 500

**02-681-169** - Pk. of 1000, Cs. of 5  
Reference; 101-1000uL; Natural; Lift-off rack; LR, 10 x  
100

**02-681-170** - Pk. of 1000, Cs. of 5  
Reference; 101-1000uL; Natural; Lift-off rack; LRS, 10 x  
100

**02-681-171** - Pk. of 1000, Cs. of 5  
Reference; 101-1000uL; Natural; Stack rack, 2 x 500

**02-681-146** - Pk. of 960, Cs. of 5  
Reference; Graduated; Volume, length: 1-2000uL, 2 in.;  
Natural; Packaging: LR, 10 x 96

**02-681-147** - Pk. of 960, Cs. of 5  
Reference; Graduated; Volume, length: 1-2000uL, 2 in.;  
Natural; Packaging: LRS, 10 x 96

**02-681-150** - Pk. of 1000, Cs. of 10  
Reference; Graduated; Volume, length: 1-2000uL, 2 in.;  
Natural; Packaging: SR, 5 x 200

**02-681-151** - Pk. of 1000, Cs. of 10  
Reference; Graduated; Volume, length: 1-200uL, 2 in.;  
Yellow; Packaging: BP

**02-681-154** - Pk. of 960, Cs. of 5  
Reference; Graduated; Volume, length: 1-200uL, 2 in.;  
Yellow; Packaging: HR, 10 x 96

**02-681-155** - Pk. of 960, Cs. of 5  
Reference; Graduated; Volume, length: 1-200uL, 2 in.;  
Yellow; Packaging: HRS, 10 x 96

**02-681-152** - Pk. of 960, Cs. of 5  
Reference; Graduated; Volume, length: 1-200uL, 2 in.;  
Yellow; Packaging: LR, 10 x 96

**02-681-156** - Pk. of 1000, Cs. of 10  
Reference; Graduated; Volume, length: 1-200uL, 2 in.;  
Yellow; Packaging: SR, 5 x 200

**02-681-153** - Pk. of 960, Cs. of 5  
Reference; Graduated; Volume/length: 1-200uL, 2 in.;  
Yellow; Packaging: LRS, 10 x 96

**02-681-163** - Pk. of 1000, Cs. of 10  
Reference; Volume, length: 101-1000uL, 2-3/4 in.; Blue;  
Packaging: BP

**02-681-165** - Pk. of 1000, Cs. of 5  
Reference; Volume, length: 101-1000uL, 2-3/4 in.; Blue;  
Packaging: Lift-Off Rack; Sterile; 10 x 96

**02-681-168** - Pk. of 1000, Cs. of 10  
Reference; Volume, length: 101-1000uL, 2-3/4 in.;  
Natural; Packaging: Bulk Pack

**02-681-457** - Pk. of 960, Cs. of 5  
Reference; Volume, Length: 1-200uL, 2 in.; Natural;  
Packaging: HRS, 10 x 96

**02-707-503** - Pk. of 480, Cs. of 10  
Type: General-Purpose; Nonsterile; Volume: 1-200uL  
volume, 2 in. long; Yellow; Packaging: Rack Inserts: 5 x  
96

**02-707-504** - Pk. of 1000, Cs. of 10  
Type: Reference Tips; Volume: 1-200uL volume, 2 in.  
long; Natural; Packaging: Bulk Pack

**02-707-506** - Pk. of 480, Cs. of 10  
Type: Reference Tips; Volume: 1-200uL volume, 2 in.  
long; Natural; Packaging: Rack Inserts, 5 x 96

**02-707-511** - Pk. of 1000, Cs. of 5  
Universal; 101-1000uL; Nonsterile

**02-681159** - Pk. of 960, Cs. of 5  
Volume: 0-300uL; Packaging Lift off Rack Sterile

**02-707-501** - Pk. of 960, Cs. of 5  
Volume: 1-200uL; Yellow; Nonsterile

### Fisherbrand Stack-Rack Space-Saver Tips



**02-707-114** - Pk. of 480, Cs. of 10  
Pipet Tips; Stack-Rack Space-Saver; Tip: Beveled;  
Volume: 1-200uL; Length: 2 in., Graduated; 5 racks of 96;  
For Pipettors: Pipetman, Eppendorf, Oxford BenchMate,  
Costar One, 8 and 12 Pette, Finnpiette, Socorex,  
Titertek Digital, MLA

**02-707-105** - Pk. of 1000, Cs. of 10  
Beveled; Natural; Volume: 1-200uL; Length: 2 in.; 5 racks  
of 200

**02-707-112** - Pk. of 480, Cs. of 10  
Beveled; Natural; Volume: 1-200uL; Length: 2 in.; 5 racks  
of 96

**02-707-115** - Pk. of 480, Cs. of 10  
Beveled; Yellow; Volume: 1-200uL; Length: 2 in.,  
Graduated; 5 racks of 96

**02-707-106** - Pk. of 1000, Cs. of 10  
Beveled; Yellow; Volume: 1-200uL; Length: 2 in.; 5 racks  
of 200

**02-707-113** - Pk. of 480, Cs. of 10  
Beveled; Yellow; Volume: 1-200uL; Length: 2 in.; 5 racks  
of 96

**02-707-108** - Pk. of 1000, Cs. of 5  
Standard; Blue; Volume: 101-1000uL; Length: 2-3/4 in.; 2  
stacks/5 racks of 100

**02-707-107** - Pk. of 1000, Cs. of 5  
Standard; Natural; Volume: 101-1000uL; Length: 2-3/4  
in.; 2 stacks/5 racks of 100

**02-707-100** - Pk. of 1000, Cs. of 10  
Standard; Natural; Volume: 1-200uL; Length: 2 in.; 5  
racks of 200

**02-707-110** - Pk. of 480, Cs. of 10  
Standard; Natural; Volume: 1-200uL; Length: 2 in.; 5  
racks of 96

**02-707-101** - Pk. of 1000, Cs. of 10  
Standard; Yellow; Volume: 1-200uL; Length: 2 in.; 5 racks  
of 200

**02-707-111** - Pk. of 480, Cs. of 10  
Standard; Yellow; Volume: 1-200uL; Length: 2 in.; 5 racks  
of 96

### Fisherbrand Standardization Tips



**02-681-142** - Pk. of 960, Cs. of 5  
Pipet Tips; Universal; Fisherbrand; Standardization;  
Pliable top opening for better fit on different types  
of pipettors; Volume/Length: 0-200uL, 2 in.; Yellow;  
Packaging: LRS, 10 x 96

**02-681-145** - Pk. of 1000, Cs. of 10  
Volume/Length: 1-200uL, 2 in.; Natural; Packaging: BP

# Pipets

## Disposable Pipets

### Fisherbrand Borosilicate Glass Disposable Serological Pipets with Regular Tip, Short Length



- 13-678-36A** - Cs. of 300  
Pipet, Serological; Fisherbrand; Borosilicate Glass; Disposable; with regular tip; Short Length; Capacity: 1mL; Subdivisions: 1/100mL; Sterile and Plugged
- 13-678-36B** - Cs. of 200  
Capacity: 5mL; Subdivisions: 1/10mL; Indiv. Wrapped
- 13-678-36C** - Cs. of 200  
Capacity: 10mL; Subdivisions: 2/10mL; Indiv. Wrapped
- 13-678-36D** - Cs. of 200  
Capacity: 25mL; Subdivisions: 2/10mL; Indiv. Wrapped
- 13-678-36E** - Cs. of 50  
Capacity: 50mL; Subdivisions: 5/10mL; Indiv. Wrapped

### Fisherbrand Borosilicate Glass Disposable Serological Pipets with Regular Tip, Standard Length



- 13-678-25A** - Cs. of 1000  
Pipet, Serological; Fisherbrand; Borosilicate Glass; Disposable; with regular tip; Standard Length; Capacity: 1mL; Subdivisions: 1/10mL; Sterile and Plugged; Packaging: Multipack
- 13-678-25B** - Cs. of 1000  
Capacity: 1mL; Subdivisions: 1/100mL; Multipack
- 13-678-25C** - Cs. of 500  
Capacity: 2mL; Subdivisions: 1/100mL; Multipack
- 13-678-25D** - Cs. of 500  
Capacity: 5mL; Subdivisions: 1/10mL; Multipack
- 13-678-25E** - Cs. of 500  
Capacity: 10mL; Subdivisions: 1/10mL; Multipack
- 13-678-26A** - Cs. of 500  
Capacity: 1mL; Subdivisions: 1/10mL; Packaging: Canisters
- 13-678-27B** - Cs. of 1000  
Capacity: 1mL; Subdivisions: 1/10mL; Individually Wrapped
- 13-678-27C** - Cs. of 1000  
Capacity: 1mL; Subdivisions: 1/100mL; Indiv. Wrapped
- 13-678-27D** - Cs. of 500  
Capacity: 2mL; Subdivisions: 1/100mL; Indiv. Wrapped
- 13-678-27E** - Cs. of 500  
Capacity: 5mL; Subdivisions: 1/10mL; Indiv. Wrapped
- 13-678-27F** - Cs. of 500  
Capacity: 10mL; Subdivisions: 1/10mL; Indiv. Wrapped

### Fisherbrand Borosilicate Glass Disposable Serological Pipets with Regular Tip, Standard Length, Nonsterile



- 13-678-31E** - Cs. of 1000  
Pipet, Serological; Fisherbrand; Borosilicate Glass; Disposable; Capacity: 1mL; Subdivisions: 1/10mL; w/ Regular Tip, Standard Length; Nonsterile; Nonplugged
- 13-678-31F** - Cs. of 1000  
Capacity: 1mL; Subdivisions: 1/100mL
- 13-678-31H** - Cs. of 500  
Capacity: 5mL; Subdivisions: 1/10mL; Case of 500
- 13-678-31J** - Cs. of 500  
Capacity: 10mL; Subdivisions: 1/10mL; 500/Cs

### Fisherbrand Disposable Borosilicate-Glass Pasteur Pipets




- 13-678-20A** - Cs. of 720  
Pipets, Pasteur; Fisherbrand; Disposable; Borosilicate glass; Tips held to close diameter tolerances; Ends are square cut; Designed from ASTM E 732 specifications; Plugs not included; Approx. Length: 5-3/4 in. (146mm); 720/Cs.
- 13-678-20B** - Cs. of 1440  
Approx. Length: 5-3/4 in. (146mm); 1440/Cs.
- 13-678-20C** - Cs. of 720  
Approx. Length: 9 in. (229mm); 720/Cs.
- 13-678-20D** - Cs. of 1440  
Approx. Length: 9 in. (229mm); 1440/Cs.

### Fisherbrand Disposable Flint-Glass Pasteur Pipets in Stand-Up Case



- 22-230-482** - Cs. of 1000  
Pipet, Pasteur; Fisherbrand; Disposable; Flint-Glass; Approx. Length: 5-3/4 in.; in Stand-Up Case


### Fisherbrand Polystyrene Disposable Serological Pipets with Magnifier Stripe



- 13-676-10G** - Cs. of 1000  
Pipet, Serological; Fisherbrand; Black Magnifier Stripe; Disposable; Polystyrene; Sterile; nonwetttable; Individual; All-plastic; Yellow; 1mL; +/-0.02mL tolerance; 1/100mL Subdiv.; 2/10mL Neg. Grads.; Ideal for tissue culture work
- 13-676-10B** - Cs. of 1000  
Capacity 1mL; Multipack
- 13-678-11B** - Cs. of 1000  
Capacity 1mL; Subdivision 1/100mL; Neg Grads 2/10mL; Tolerance +/-0.02mL; Paper/Plastic; Yellow
- 13-675-3C** - Cs. of 500  
Capacity 2mL; Individually wrapped, All-plastic wrap
- 13-675-2C** - Cs. of 500  
Capacity 2mL; Subdivision 1/100mL; Neg Grads 2/10mL; Tolerance +/-0.04mL; Multipack; Green
- 13-678-11C** - Cs. of 500  
Capacity 2mL; Subdivision 1/100mL; Neg Grads 2/10mL; Tolerance +/-0.04mL; Paper/Plastic; Green
- 13-678-11D** - Cs. of 200  
Capacity 5mL; +/-0.1mL tolerance; 1/10mL Subdiv.; 2.5mL Neg. Grads; Blue
- 13-676-10H** - Cs. of 200  
Capacity 5mL; Individually wrapped, All-plastic wrap
- 13-676-10C** - Cs. of 500  
Capacity 5mL; Subdivision; 1/10mL; Neg Grads 2.5mL; Tolerance +/-0.1mL; Multipack; Blue
- 13-676-10J** - Cs. of 200  
Capacity 10mL; Individually wrapped; All-plastic wrap
- 13-676-10F** - Cs. of 500  
Capacity 10mL; Multipack
- 13-678-11E** - Cs. of 200  
Capacity 10mL; Subdivision 1/10mL; Neg Grads 3mL; Tolerance +/-0.2mL; Paper/Plastic; Orange
- 13-676-10K** - Cs. of 200  
Capacity 25mL; Individually wrapped, All-plastic wrap
- 13-676-10M** - Cs. of 200  
Capacity 25mL; Subdivision 2/10mL; Neg Grads 10mL; Tolerance +/-0.5mL; Multipack; Red
- 13-678-11** - Cs. of 200  
Capacity 25mL; Subdivision 2/10mL; Neg Grads 10mL; Tolerance +/-0.5mL; Paper/Plastic; Red
- 13-676-10R** - Cs. of 100  
Capacity 50mL; +/-1mL tolerance; 5/10mL Subdiv.; Purple
- 13-676-10Q** - Cs. of 100  
Capacity 50mL; Subdivision 5/10mL; Tolerance +/-1mL; Multipack; Purple
- 13-678-11F** - Cs. of 100  
Capacity 50mL; Subdivision 5/10mL; Tolerance +/-1mL; Paper/Plastic; Purple



## Fisherbrand Polystyrene Disposable Standard Serological Pipets



- 13-678-12B** - Cs. of 1000  
Pipet, Serological; Fisherbrand; Polystyrene; Disposable;  
Subdivisions: 1/100mL; Capacity: 1mL
- 13-678-12C** - Cs. of 500  
Capacity: 2mL
- 13-678-12D** - Cs. of 500  
Capacity: 5mL
- 13-678-12E** - Cs. of 500  
Capacity: 10mL
- 13-678-14B** - Cs. of 200  
Capacity: 25mL
- 13-678-14C** - Cs. of 100  
Capacity: 50mL

- 13-671-108D** - Cs. of 12  
Capacity: 10 in 1/10mL
- 13-671-108E** - Cs. of 12  
Capacity: 25 in 1/10mL

## Fisherbrand Long Tip Serological Pipets



- 13-674-26H** - Cs. of 12  
Pipet, Serological; Fisherbrand; Long Tip; Capacity: 2 in  
1/10mL; Accuracy: +/-0.03mL
- 13-674-26M** - Cs. of 12  
Capacity: 10 in 1/10mL; Accuracy: +/-0.06mL

## Fisherbrand Standard-Top Serological Pipets



- 13-675G** - Cs. of 12  
Pipet, Serological; Fisherbrand; Standard-top;  
Borosilicate glass; Permanent amber markings; Tuf-Temp  
treated; Calibrated To Deliver/Blow Out; Capacity: 1 in  
1/10mL; Accuracy: +/-0.02mL; Red
- 13-675F** - Cs. of 12  
Capacity: 1 in 1/100mL; Accuracy: +/-0.02mL; Yellow
- 13-675D** - Cs. of 12  
Capacity: 5/10 in 1/100mL; Accuracy: +/-0.01mL; Yellow

## Fisherbrand Standard Disposable Transfer Pipets



- 13-711-7M** - Pk. of 500, Cs. of 10  
Pipets, Transfer; Fisherbrand; Disposable; Standard;  
One-piece; Single squeeze draws up to 3.2mL; 5-7/8 in  
(15cm) long; Capacity: 7.7mL; Nongraduated
- 13-711-9D** - Pk. of 500, Cs. of 10  
Graduated; Capacity: 7.5mL; Length: 6 in.

## Volumetric Pipets

### Fisherbrand Reusable Class A Borosilicate Glass Volumetric Pipets



- 13-650-2A** - Cs. of 12  
Pipets, Volumetric; Fisherbrand; Reusable; Class A;  
Borosilicate Glass; Type: Standard; Capacity: 0.5mL
- 13-650-2B** - Cs. of 12  
Type: Standard; Capacity: 1mL
- 13-650-2C** - Cs. of 12  
Type: Standard; Capacity: 2mL
- 13-650-2D** - Cs. of 12  
Type: Standard; Capacity: 3mL
- 13-650-2E** - Cs. of 12  
Type: Standard; Capacity: 4mL
- 13-650-2F** - Cs. of 12  
Type: Standard; Capacity: 5mL
- 13-650-2G** - Cs. of 12  
Type: Standard; Capacity: 6mL
- 13-650-2H** - Cs. of 12  
Type: Standard; Capacity: 7mL
- 13-650-2J** - Cs. of 12  
Type: Standard; Capacity: 8mL
- 13-650-2K** - Cs. of 12  
Type: Standard; Capacity: 9mL
- 13-650-2L** - Cs. of 12  
Type: Standard; Capacity: 10mL
- 13-650-2M** - Ea., Cs. of 12  
Type: Standard; Capacity: 15mL
- 13-650-2N** - Ea., Cs. of 12  
Type: Standard; Capacity: 20mL
- 13-650-2P** - Ea., Cs. of 12  
Type: Standard; Capacity: 25mL
- 13-650-2Q** - Ea., Cs. of 12  
Type: Standard; Capacity: 30mL
- 13-650-2R** - Ea., Cs. of 12  
Type: Standard; Capacity: 40mL
- 13-650-2S** - Ea., Cs. of 12  
Type: Standard; Capacity: 50mL
- 13-650-2T** - Ea., Cs. of 12  
Type: Standard; Capacity: 75mL
- 13-650-2U** - Ea., Cs. of 12  
Type: Standard; Capacity: 100mL
- 13-650-2W** - Ea., Cs. of 4  
Type: Standard; Capacity: 200mL
- 13-650-2X** - Cs. of 12  
Type: Standard; Capacity: 12mL
- 13-650-2Y** - Cs. of 12  
Type: Standard; Capacity: 2.5mL
- 13-650-3A** - Cs. of 6  
Type: Serialized; Capacity: 0.5mL
- 13-650-3B** - Cs. of 6  
Type: Serialized; Capacity: 1mL
- 13-650-3C** - Cs. of 6  
Type: Serialized; Capacity: 2mL
- 13-650-3D** - Cs. of 6  
Type: Serialized; Capacity: 3mL

## Measuring Pipets

### Fisherbrand Measuring (Mohr) Standard-Taper Pipets



- 13-665A** - Cs. of 12  
Pipets, Reusable Measuring; Fisherbrand; Standard-  
Taper; Borosilicate glass; Permanent amber markings;  
Capacity: 1/10 in 1/100mL; Accuracy: +/-0.005mL; Color  
Code: White
- 13-665F** - Cs. of 12  
Capacity: 1 in 1/100mL; Accuracy: +/-0.02mL
- 13-665G** - Cs. of 12  
Capacity: 1 in 1/10mL; Accuracy: +/-0.02mL
- 13-665H** - Cs. of 12  
Capacity: 2 in 1/10mL; Accuracy: +/-0.02mL
- 13-665K** - Cs. of 12  
Capacity: 5 in 1/10mL; Accuracy: +/- 0.04mL
- 13-665M** - Cs. of 12  
Capacity: 10 in 1/10mL; Accuracy: +/- 0.06mL
- 13-665N** - Cs. of 12  
Capacity: 25 in 1/10mL; Accuracy: +/-0.10mL

## Reusable Pipets

### Fisherbrand Large Tip Opening Serological Pipets



- 13-671-108A** - Cs. of 12  
Pipet, Serological; Fisherbrand; Large Tip Opening;  
Capacity: 1 in 1/10mL; Accuracy: +/-0.05mL; Color Code:  
Red
- 13-671-108B** - Cs. of 12  
Capacity: 2 in 1/10mL
- 13-671-108C** - Cs. of 12  
Capacity: 5 in 1/10mL

- 13-650-3E** - Cs. of 6  
Type: Serialized; Capacity: 4mL
- 13-650-3F** - Cs. of 6  
Type: Serialized; Capacity: 5mL
- 13-650-3G** - Cs. of 6  
Type: Serialized; Capacity: 6mL
- 13-650-3H** - Cs. of 6  
Type: Serialized; Capacity: 7mL
- 13-650-3J** - Cs. of 6  
Type: Serialized; Capacity: 8mL
- 13-650-3K** - Cs. of 6  
Type: Serialized; Capacity: 9mL
- 13-650-3L** - Cs. of 6  
Type: Serialized; Capacity: 10mL
- 13-650-3M** - Ea., Cs. of 6  
Type: Serialized; Capacity: 15mL
- 13-650-3N** - Ea., Cs. of 6  
Type: Serialized; Capacity: 20mL
- 13-650-3P** - Ea., Cs. of 6  
Type: Serialized; Capacity: 25mL
- 13-650-3Q** - Ea., Cs. of 6  
Type: Serialized; Capacity: 30mL
- 13-650-3R** - Ea., Cs. of 6  
Type: Serialized; Capacity: 40mL
- 13-650-3S** - Ea., Cs. of 6  
Type: Serialized; Capacity: 50mL
- 13-650-3T** - Ea., Cs. of 6  
Type: Serialized; Capacity: 75mL
- 13-650-3V** - Ea., Cs. of 6  
Type: Serialized; Capacity: 100mL
- 13-650-3W** - Ea.  
Type: Serialized; Capacity: 200mL
- 13-650-3X** - Cs. of 6  
Type: Serialized; Capacity: 12mL
- 13-650-3Y** - Cs. of 6  
Type: Serialized; Capacity: 2.5mL

- 13-681-503** - Cs. of 20  
For 50mL Basins; APET; Sterile, Individually Bagged; L x W x H: 5-3/8 x 2-1/2 x 1/2 in.; Pack of: 20
- 13-681-504** - Cs. of 200  
100mL; Polystyrene; Sterile; Five per bag; L x W x H: 5-3/4 x 3-1/2 x 1-1/4 in.; Pack of: 200

**Fisherbrand Pipetter Stand Workstation**



- 13-681-24** - Ea.  
Pipet Stand Workstation; Fisherbrand; L x W x H: 6-5/8 x 16-3/8 x 9-1/2 in.; With two bins: L x W x H: 3 x 8 x 5 in.; White Acrylic with Transparent Bins

# Racks

## Microcentrifuge Tube Racks

**Fisherbrand 96-Well Racks for PCR Tubes**



- 03-448-20** - Pk. of 5, Cs. of 40  
Tube Rack, Fisherbrand; Slip-Fit Top; 96-Well; 5 x 4 x 1-1/2 in.; Reversible; Std. 8x12 Array; Snap-In-Place Lid; Easy-To-Grasp; Grid Lines on Wells; For 1.5 to 2.0mL Tubes; Assorted Colors; Molded/Raised Rungs

**Fisherbrand Interlocking Four-Way Tube Racks**



- 03-448-12** - Pk. of 5  
Rack, Tube; Fisherbrand; Interlocking Four-Way; 2 L x -1/2 W x 7 in. H; Blue; Accommodate Tubes From 0.5 to 50mL; Autoclavable; Polypropylene
- 03-448-13** - Pk. of 5  
Color: Fluorescent green
- 03-448-16** - Pk. of 5  
Color: Fluorescent Orange
- 03-448-14** - Pk. of 5  
Color: Fluorescent Pink
- 03-448-15** - Pk. of 5  
Color: Fluorescent Yellow
- 03-448-11** - Pk. of 5  
Color: Natural
- 03-448-17** - Pk. of 5  
Color: Assorted Colors

**Fisherbrand Microcentrifuge Tube Racks**



- 05-541-37** - Pk. of 5, Cs. of 4  
Rack, Microcentrifuge Tube; Reversible; Polypropylene; Capacity: 96 0.5mL tubes or 96 1.5 to 2.0mL microcentrifuge tubes; Color: Black
- 05-541-45** - Pk. of 5, Cs. of 4  
Color: Blue
- 05-541-39** - Pk. of 5, Cs. of 4  
Color: Fluorescent Blue
- 05-541-40** - Pk. of 5, Cs. of 4  
Color: Fluorescent Green
- 05-541-41** - Pk. of 5, Cs. of 4  
Color: Fluorescent Orange
- 05-541-43** - Pk. of 5, Cs. of 4  
Color: Fluorescent Yellow
- 05-541-46** - Pk. of 5, Cs. of 4  
Color: Green
- 05-541-36** - Pk. of 5, Cs. of 4  
Color: Lavender
- 05-541-29** - Pk. of 5, Cs. of 4  
Color: Natural
- 05-541-47** - Pk. of 5, Cs. of 4  
Color: Orange
- 05-541-48** - Pk. of 5, Cs. of 4  
Color: Pink

# Pipettors

**Fisherbrand Rubber Bulb-Type Safety Pipet Fillers**



- 13-681-50** - Ea., Cs. of 72  
Dispenser, Pipet Filler; Fisherbrand; Bulb-type; Safety; Rubber; With two stainless steel ball valves; One hand use; Fit pipets with O.D. from 3/16 to 3/8 in.; Color: Black
- 13-681-51** - Ea., Cs. of 72  
Red

# Pipetter and Pipet Accessories

**Fisherbrand Autoclavable Pipet Basins**



- 13-681-800** - Cs. of 20  
Pipet Basin; Fisherbrand; Autoclavable, 12-Channel with Lid; Polypropylene; Sterile; Individually Bagged; L x W x H: 5 x 2-1/4 x 1 in.
- 13-681-801** - Cs. of 50  
12-Channel with Lid; Nonsterile; L x W x H: 5 x 2-1/4 x 1 in.
- 13-681-802** - Cs. of 10  
Multichannel; Polypropylene; Nonsterile; L x W x H: 4-3/4 x 4-1/2 x 1 in.

**Fisherbrand Disposable Pipet Basins**



- 13-681-500** - Cs. of 100  
Basin, Pipet; Fisherbrand; Disposable; 50mL, PVC, Nonsterile, Bulk Packaging, L x W x H: 5-3/4 x 2-3/4 x 1 in.; Pack of: 100
- 13-681-501** - Cs. of 80  
50mL; Polystyrene; Sterile; L x W x H: 5 3/4 x 2 3/4 x 1 in.; Pack of: 80
- 13-681-502** - Cs. of 200  
50mL; Polystyrene; Sterile; Five per bag; L x W x H: 5-3/4 x 2-3/4 x 1 in.; Pack of: 200

**05-541-49** - Pk. of 5, Cs. of 4  
Color: Yellow  
**05-541-38** - Pk. of 5, Cs. of 4  
Color: Fluorescent Assorted  
**05-541-44** - Pk. of 5, Cs. of 4  
Color: Nonfluorescent Assorted

## Special Purpose Racks

### Fisherbrand Cryogenic Vial Storage Boxes



**10-500-28** - Pk. of 4, Cs. of 6  
Dispenser for Pipet Tips; Fisherbrand; Polypropylene; inner plunger of polyethylene; designed to fit Eppendorf Repeater Pipetter; Dispense from 8 to 48 times without refilling; Nonsterile; 1.25mL; +/- 0.8% accuracy; long, thin  
**10-500-29** - Pk. of 4, Cs. of 6  
Holds up to 5.0mL vials

## Refrigerator and Freezer Accessories

### Fisherbrand Stainless-Steel Freezer Racks



**03-395-632** - Ea.  
Rack, Freezer; Fisherbrand; Sturdy, rust-resistant stainless steel; Chest freezer rack; 2 in. boxes; 4 shelves  
**03-395-634** - Ea.  
Chest freezer rack; 2 in. boxes; 6 shelves  
**03-395-631** - Ea.  
Chest freezer rack; 2 in. boxes; 9 shelves  
**03-395-613** - Ea.  
Chest freezer rack; 2 in. boxes; 10 shelves  
**03-395-614** - Ea.  
Chest freezer rack; 2 in. boxes; 11 shelves  
**03-395-615** - Ea.  
Chest freezer rack; 2 in. boxes; 12 shelves  
**03-395-616** - Ea.  
Chest freezer rack; 2 in. boxes; 16 shelves  
**03-395-628** - Ea.  
Chest freezer rack; For 15mL and 50mL tube boxes; 5 shelves  
**03-395-629** - Ea.  
Chest freezer rack; For 15mL and 50mL tube boxes; 6 shelves  
**03-395-630** - Ea.  
Mini-vertical freezer rack; For mini 2 in. boxes; 9 shelves  
**03-395-638** - Ea.  
Mini-vertical freezer rack; For mini 2 in. boxes; 12 shelves  
**03-395-622** - Ea.  
Upright freezer drawer rack; 15 and 50mL Centrifuge Tube Boxes; Holds 78-50mL tubes  
**03-395-621** - Ea.  
Upright freezer drawer rack; 15 and 50mL Centrifuge Tube Boxes; Holds 120-15mL tubes  
**03-395-624** - Ea.  
Upright freezer drawer rack; 15 and 50mL Centrifuge Tube Boxes; Holds 120-15mL tubes, With handle  
**03-395-623** - Ea.  
Upright freezer drawer rack; 15 and 50mL Centrifuge Tube Boxes; Holds 160-15mL tubes  
**03-395-608** - Ea.  
Upright freezer drawer rack; 2 in. boxes; 9 shelves (3 x 3 array)  
**03-395-610** - Ea.  
Upright freezer drawer rack; 2 in. boxes; 12 shelves (4 x 3 array, 16x9x5)  
**03-395-609** - Ea.  
Upright freezer drawer rack; 2 in. boxes; 12 shelves (4 x 3 array, 22x9x5)  
**03-395-611** - Ea.  
Upright freezer drawer rack; 2 in. boxes; 16 shelves (4 x 4 array)  
**03-395-612** - Ea.  
Upright freezer drawer rack; 2 in. boxes; 20 shelves (4 x 5 array)

**03-395-618** - Ea.  
Upright freezer drawer rack; 3 in. boxes; 12 shelves (4 x 3 array)  
**03-395-626** - Ea.  
Upright freezer rack; 15mL and 50mL Tube Boxes; 6 shelves (3 x 2 array)  
**03-395-627** - Ea.  
Upright freezer rack; 100-place cell hinged-top plastic storage boxes; 8 shelves (4 x 2 array)  
**03-395-625** - Ea.  
Upright freezer rack; 100-place cell hinged-top plastic storage boxes; 12 shelves  
**03-395-600** - Ea.  
Upright freezer rack; 2 in. boxes; 9 shelves (3 x 3 array)  
**03-395-601** - Ea.  
Upright freezer rack; 2 in. boxes; 12 shelves (3 x 4 array)  
**03-395-603** - Ea.  
Upright freezer rack; 2 in. boxes; 12 shelves (4 x 3 array)  
**03-395-602** - Ea.  
Upright freezer rack; 2 in. boxes; 15 shelves (3 x 5 array)  
**03-395-604** - Ea.  
Upright freezer rack; 2 in. boxes; 16 shelves (4 x 4 array)  
**03-395-605** - Ea.  
Upright freezer rack; 2 in. boxes; 20 shelves (4 x 5 array)  
**03-395-607** - Ea.  
Upright freezer rack; 2 in. boxes; 20 shelves (5 x 4 array)  
**03-395-606** - Ea.  
Upright freezer rack; 2 in. boxes; 24 shelves (4 x 6 array)  
**03-395-617** - Ea.  
Upright freezer rack; 3 in. boxes; 12 shelves (4 x 3 array)  
**03-395-619** - Ea.  
Upright freezer rack; 96- and 384-well Microtiter Plates; Holds 84 plates w/lids, 96 w/o lids  
**03-395-620** - Ea.  
Upright freezer rack; 96- and 384-well Microtiter Plates; Holds 96 plates w/lids, 120 w/o lids

## Samplers

### Fisherbrand Reusable Polypropylene Scoops



**14-375-258** - Cs. of 12  
Scoops; Fisherbrand Reusable Clear Polypropylene; Capacity: 75cc; Length; 7-1/2 in. (19cm); Bulk packaged; Nonsterile; Steam autoclavable at 121 deg. C (250 deg. F)  
**14-375-259** - Cs. of 12  
Capacity: 100cc; Length: 7 in.  
**14-375-260** - Cs. of 12  
Capacity: 200cc; Length: 8 in.

## Scissors

### Fisherbrand Delicate Dissecting Scissors



**08-951-10** - Ea., Cs. of 10  
Scissors, Dissecting; Fisherbrand; Forged surgical-grade stainless steel; FDA certified material; Can be sterilized; Passivated; Autoclavable; Rockwell: HRC 46-50 deg. ; 4-1/2 in.; Curved  
**08-951-5** - Ea., Cs. of 10  
Straight

### Fisherbrand Doyen Abdominal Scissors



**13-804-4** - Ea., Cs. of 10  
Scissors, Abdominal; Fisherbrand; Doyen; Forged from surgical-grade stainless steel; FDA certified material; Can be sterilized; Passivated; Autoclavable; Rockwell: HRC 46-50 deg. ; Curved  
**13-804-2** - Ea., Cs. of 10  
Straight

### Fisherbrand Mayo Dissecting Scissors



**13-804-5** - Ea., Cs. of 6  
Scissors, Dissecting; Fisherbrand; Mayo style; Forged surgical-grade stainless steel; FDA certified material; Can be sterilized; Passivated; Autoclavable; Rockwell: HRC 46-50 deg. ; 5-1/2 in.; Curved



- 13-804-6 - Ea., Cs. of 6  
5 1/2 in., Straight
- 13-804-12 - Ea., Cs. of 6  
6 3/4 in., Curved
- 13-804-8 - Ea., Cs. of 6  
6 3/4 in., Straight

### Fisherbrand Sharp-Pointed Dissecting Scissors



- 08-935 - Ea., Cs. of 10  
Scissors, Dissecting; Fisherbrand; Sharp-pointed; Finely finished; Narrow blades; 1-in. edge; Sharp/sharp; FDA-certified material; Sterilizable; Passivated; Autoclavable; HRC 46-50 deg. ; Curved; Length: 4-1/8 in.
- 08-940 - Ea., Cs. of 12  
Style: Straight; Length: 4-1/2 in.

### Fisherbrand Standard Dissecting Scissors



- 08-951-25 - Ea., Cs. of 10  
Scissors, Dissecting; Fisherbrand; Forged surgical-grade stainless steel; FDA certified material; Can be sterilized; Passivated; Autoclavable; Rockwell: HRC 46-50 deg. ; Length: 5-1/8 in.; Curved
- 08-951-20 - Ea., Cs. of 10  
Style: Straight; Length: 6 in.

### Fisherbrand Straight-Blade Operating Scissors



- 13-806 - Ea., Cs. of 10  
Scissors, Operating; Fisherbrand; Straight-blade; Surgical-grade stainless steel; Sterilizable; Passivated; Autoclavable; HRC 46-50 deg. ; Blunt/blunt; Length: 4-1/2 in.
- 13-806-2 - Ea., Cs. of 10  
Blunt/sharp; Length: 4-1/2 in.
- 13-808-2 - Ea., Cs. of 10  
Blunt/sharp; Length: 5-1/2 in.
- 13-810-2 - Ea., Cs. of 10  
Blunt/sharp; Length: 6-1/2 in.
- 13-808-4 - Ea., Cs. of 10  
Sharp/sharp; Length: 5-1/2 in.
- 13-810-4 - Ea.  
Sharp/sharp; Length: 6-1/2 in.

## Spatulas

### Fisherbrand Chemi-Scraper Spatulas



- 14-373 - Ea., Pk. of 12  
Spatula, Chemi-Scraper; Fisherbrand; Type 18-8 nickel stainless steel with bright satin finish; One rounded and one square end; 7 in. L

### Fisherbrand Scoopula Spatulas



- 14-357Q - Pk. of 6  
Spatula, Scoopula; Fisherbrand; Stainless Steel; Minimizes spilling in handling small amounts of material; Will not rust; One pointed end, one blunt end; Length: 6 in.

### Fisherbrand Semimicro Spatula with One Tapered End, One Rounded End



- 14-374 - Ea., Pk. of 12  
Spatula, Semimicro; Fisherbrand; Stainless Steel; One end tapered to 0.14 in. (3.5mm) rounded tip for use in conical-bottom tubes; Rounded end has 1/4 in. (6.6mm) rounded tip for regular tubes; Length: 7 in.

### Fisherbrand Spatula/Scraper



- 14-373-25A - Ea., Pk. of 12  
Spatula/Scraper; Fisherbrand; Type 18-8 nickel stainless steel; Stem dia.: 0.156 in. (4mm); Double flat ends: one rounded, one squared; Length: 12 in.
- 14-373-25B - Ea., Pk. of 12  
Length: 18 in.

### Fisherbrand Spoonula Lab Spoons



- 14-375-10 - Pk. of 6, Cs. of 12  
Spatula; Fisherbrand; FDA-certified polished stainless steel; Spoon end is 1-1/4 x 9/16 in.; Spatula end is 2 x 5/16 in.; Length: 9 in.

## Stir Bars

### Fisherbrand 8-Agon Stir Bars



- 14-512-145 - Ea., Cs. of 50  
Stir Bar, Magnetic; Fisherbrand; 8-Agon; Octagonal; PTFE; Blue; 8 dia. x 12.7mm L
- 14-512-146 - Ea., Cs. of 50  
Blue; 8 dia. x 16mm L
- 14-512-147 - Ea., Cs. of 50  
Blue; 8 dia. x 22mm L
- 14-512-148 - Ea., Cs. of 50  
Blue; 8 dia. x 25.4mm L
- 14-512-149 - Ea., Cs. of 50  
Blue, 8 dia. x 38mm L
- 14-512-1451 - Ea., Cs. of 50  
Red; 8 dia. x 12.7mm L
- 14-512-1461 - Ea., Cs. of 50  
Red; 8 dia. x 16mm L
- 14-512-1471 - Ea., Cs. of 50  
Red; 8 dia. x 22mm L
- 14-512-1481 - Ea., Cs. of 50  
Red; 8 dia. x 25.4mm L
- 14-512-1491 - Ea., Cs. of 50  
Red; 8 dia. x 38mm L
- 14-512-1452 - Ea., Cs. of 50  
Yellow; 8 dia. x 12.7mm L
- 14-512-1462 - Ea., Cs. of 50  
Yellow; 8 dia. x 16mm L
- 14-512-1472 - Ea., Cs. of 50  
Yellow; 8 dia. x 22mm L
- 14-512-1482 - Ea., Cs. of 50  
Yellow; 8 dia. x 25.4mm L
- 14-512-1492 - Ea., Cs. of 50  
Yellow; 8 dia. x 38mm L

### Fisherbrand Blue Micro Stir Bars



- 14-512-150 - Ea., Cs. of 50  
Stir Bar, Magnetic; Fisherbrand Micro; PTFE; Blue; 1.5 dia. x 8 mm L
- 14-512-151 - Ea., Cs. of 50  
1.5 dia. x 15mm L
- 14-512-152 - Ea., Cs. of 50  
2 dia. x 2mm L
- 14-512-153 - Ea., Cs. of 50  
2 dia. x 5mm L
- 14-512-154 - Ea., Cs. of 50  
2 dia. x 7mm L
- 14-512-155 - Ea., Cs. of 50  
3 dia. x 3mm L
- 14-512-156 - Ea., Cs. of 50  
3 dia. x 10mm L

### Fisherbrand Color-Coded Dumbbell-Style Stir Bars



- 14-512-159 - Ea., Cs. of 50  
Stir Bar, Magnetic; Fisherbrand; PTFE; Dumbbell-style; Blue; 8 dia. x 32mm L
- 14-512-158 - Ea., Cs. of 50  
Red; 8 dia. x 32mm L
- 14-512-160 - Ea., Cs. of 50  
Red, 8 dia. x 52mm L

### Fisherbrand Egg-Shaped Bars



- 14-512-122 - Ea., Cs. of 50  
Stir Bar, Magnetic; Fisherbrand; Egg-shaped; PTFE; For round-bottom vessels; 6 dia. x 16mm L
- 14-513-53 - Ea., Cs. of 50  
6 dia. x 32mm L

- 14-512-121** - Ea., Cs. of 50  
9.5 dia. x 19mm L
- 14-512-120** - Ea., Cs. of 50  
12.7 dia. x 25.4mm L
- 14-512-123** - Ea., Cs. of 50  
16 dia. x 38mm L
- 14-513-54** - Ea., Cs. of 50  
19 dia. x 41mm L
- 14-513-55** - Ea., Cs. of 50  
19 dia. x 64mm L

### Fisherbrand Magnetic Micro Stirring Bars



- 14-513-65** - Ea., Cs. of 50  
Stir Bar, Magnetic; Fisherbrand; Octagonal; 5/16 in. dia.; 2 in. L; pivot rings; Alnico\_V magnet encased in PTFE; flat surface

### Fisherbrand Magnetic Stir Bar Six Packs



- 14-513-64SIX** - Pk. of 6  
Stir Bar, Magnetic; Fisherbrand Micro; PTFE; Packaged in snap-cap polyethylene vial; 1.5 dia. x 8mm L
- 14-513-63SIX** - Pk. of 6  
Micro, 3/32 in. dia.; 1/4 in. L
- 14-513-65SIX** - Pk. of 6  
Micro, 1/8 in. dia.; 3/8 in. L
- 14-513-57SIX** - Pk. of 6  
Octagonal, 1/8 in. dia.; 0.5 in. L
- 14-513-58SIX** - Pk. of 6  
Octagonal, 5/16 in. dia. x 1/2 in. L
- 14-513-59SIX** - Pk. of 6  
Octagonal, 5/16 in. dia.; 1 in. L
- 14-513-60SIX** - Pk. of 6  
Octagonal, 5/16 in. dia.; 1-1/2 in. L
- 14-513-61SIX** - Pk. of 6  
Octagonal, 5/16 in. dia.; 2 in. L
- 14-513-51SIX** - Pk. of 6  
Octagonal, 3/8 in. dia.; 1 in. L
- 14-513-66SIX** - Pk. of 6  
Octagonal, 3/8 in. dia.; 2 in. L

### Fisherbrand Magnetic Stirring Bar Assortments



- 14-512-129** - Pk. of 12, Cs. of 6  
Stir Bar, Magnetic; Fisherbrand; Assortments; PTFE; Octagon; 12 bars w/pivot ring

### Fisherbrand Octagonal Magnetic Stir Bars



- 14-513-57** - Ea., Cs. of 50  
Stir Bar, Magnetic; Fisherbrand; Octagonal; 1/8 in. dia.; 1/2 in. L; pivot rings; Alnico\_V magnet encased in PTFE; flat surface
- 14-513-58** - Ea., Cs. of 50  
5/16 in. dia.; 1/2 L
- 14-513-62** - Ea., Cs. of 50  
5/16 in. dia.; 5/8 in. L
- 14-513-59** - Ea., Cs. of 50  
5/16 in. dia.; 1 in. L
- 14-513-60** - Ea., Cs. of 50  
5/16 in. dia.; 1-1/2 in. L
- 14-513-51** - Ea., Cs. of 50  
3/8 in. dia.; 1 in. L
- 14-513-52** - Ea., Cs. of 50  
3/8 in. dia.; 1-1/2 in. L
- 14-513-66** - Ea., Cs. of 50  
3/8 in. dia.; 2 in. L
- 14-513-67** - Ea., Cs. of 50  
3/8 in. dia.; 2-1/2 in. L
- 14-513-56** - Ea., Cs. of 50  
1/2 in. dia.; 1-1/2 in. L
- 14-513-68** - Ea., Cs. of 50  
1/2 in. dia.; 3 in. L;

## Storage Microplates

### Fisherbrand 96-Well DeepWell Polypropylene Microplates


- 12-566-120** - Cs. of 50  
Microplate; Cell Culture; Fisherbrand; Deepwell; Polypropylene; 96-well; Color: Natural; Nonsterile; Well Capacity: 1mL; Working Vol.: 50 to 1uL; Autoclavable
- 12-566-121** - Cs. of 60  
Well Capacity: 2mL

### Fisherbrand 96-Well Round Bottom Polypropylene Plates



- 12-565-505** - Cs. of 80  
Plate, Polypropylene; Fisherbrand; 96-well; round bottom; Blue; withstand temps from -80 to 121 deg. C; nonsterile; 50uL well volume
- 12-565-502** - Cs. of 80  
Natural
- 12-565-503** - Cs. of 80  
Red
- 12-565-504** - Cs. of 80  
Yellow

### Fisherbrand 384-Well Polypropylene Plates



- 12-565-507** - Cs. of 120  
Plate, Polypropylene; Fisherbrand; 384-well; Natural; nonsterile; chemical resistant; withstands temps from -80 to +80 deg. C; 120uL well volume; Ideal for storage of compounds and nucleic acid libraries
- 12-565-508** - Cs. of 120  
384-Well; Red
- 12-565-509** - Cs. of 120  
384-Well; Yellow

## Tape, Labels, and Labeling Tapes

### Fisherbrand White Autoclavable Tape



- 15-903** - Ea.  
Tape, Autoclave; Fisherbrand; White with white steam-sensitive stripes that turn black during sterilization; No residue; Size: 1/2 in. W x 60 yd. L
- 15-904** - Ea.  
3/4 in. x 60 yd.
- 15-905** - Ea.  
1 in. x 60 yd.

### Micryo Dots and Strips for Cryo Storage



- 03-391-188** - Pk. of 3840  
Label, Cryo Storage; Micryo Dots for Laser Printers; Use in temperatures as low as -196 deg. C; Suitable for use on microcentrifuge tubes, cryovials, racks, boxes, and microtiter plates; White; For 0.5-1.5mL tubes
- 03-391-111** - Pk. of 3840  
Micryo Dots for Laser Printers; White; For 2.0 and 5.0mL tubes
- 03-391-105** - Pk. of 2160  
Micryo Dots for Laser Printers; White; For 15 and 50mL tubes
- 03-391-102** - Pk. of 5000  
Micryo Dots; Rainbow; For 0.5-1.5mL tubes
- 03-391-187** - Pk. of 1000  
Micryo Dots; White; For 0.5-1.5mL tubes; 1000/Roll
- 03-391-101** - Pk. of 5000  
Micryo Dots; White; For 0.5-1.5mL tubes; 5000/Roll
- 03-391-189** - Pk. of 1000  
Micryo Dots; White; For 1.5 and 2.0mL tubes
- 03-391-103** - Pk. of 1000  
Micryo Dots; White; For 15 and 50mL tubes; 1000/Roll
- 03-391-104** - Pk. of 5000  
Micryo Dots; White; For 15 and 50mL tubes; 5000/Roll
- 03-391-109** - Pk. of 2975  
Micryo Strips for Laser Printers; Rainbow; For 0.5mL tubes

- 03-391-107** - Pk. of 2125  
Micryo Strips for Laser Printers; Rainbow; For 1.5-2mL tubes
- 03-391-113** - Pk. of 3900  
Micryo Strips for Laser Printers; Rainbow; For 96/384 plates
- 03-391-108** - Pk. of 2380  
Micryo Strips for Laser Printers; White; For 0.5mL tubes
- 03-391-106** - Pk. of 1700  
Micryo Strips for Laser Printers; White; For 1.5-2mL tubes
- 03-391-112** - Pk. of 3120  
Micryo Strips for Laser Printers; White; For 96/384 plates
- 03-391-110** - Pk. of 600  
Micryo Strips for Laser Printers; White; For Racks/boxes
- 03-391-115** - Pk. of 1000  
Micryo Strips; Blue; For 1.5-2mL tubes
- 03-391-117** - Pk. of 1000  
Micryo Strips; Green; For 1.5-2mL tubes
- 03-391-119** - Pk. of 1000  
Micryo Strips; Purple; For 1.5-2mL tubes
- 03-391-116** - Pk. of 1000  
Micryo Strips; Red; For 1.5-2mL tubes
- 03-391-100** - Pk. of 1000  
Micryo Strips; White; For 1.5-2mL tubes
- 03-391-114** - Pk. of 1000  
Micryo Strips; White; For 2mL tubes
- 03-391-118** - Pk. of 1000  
Micryo Strips; Yellow; For 1.5-2mL tubes

# Thermometers - Nonmercury

## Fisherbrand Partial-Immersion Thermometers With Standard Taper Joint



- 15-041-11A** - Ea.  
Thermometers; Fisherbrand Partial-Immersion; With Std. Taper Joint; High Range: +150 deg. C; Total Length: 225mm; Low Range: -20 deg. C
- 15-041-11C** - Ea.  
Total Length: 300mm; High Range: +150 deg. C; Low Range: -20 deg. C
- 15-041-11E** - Ea.  
Total Length: 325mm; High Range: +150 deg. C; Low Range: -20 deg. C
- 15-041-11B** - Ea.  
Total Length: 325mm; High Range: +250 deg. C; Low Range: -10 deg. C
- 15-041-11G** - Ea.  
Total Length: 350mm; High Range: +150 deg. C; Low Range: -20 deg. C
- 15-041-11D** - Ea.  
Total Length: 350mm; High Range: +250 deg. C; Low Range: -10 deg. C
- 15-041-11J** - Ea.  
Total Length: 375mm; High Range: +150 deg. C; Low Range: -20 deg. C
- 15-041-11F** - Ea.  
Total Length: 375mm; High Range: +250 deg. C; Low Range: -10 deg. C
- 15-041-11L** - Ea.  
Total Length: 400mm; High Range: +150 deg. C; Low Range: -20 deg. C
- 15-041-11H** - Ea.  
Total Length: 400mm; High Range: +250 deg. C; Low Range: -10 deg. C
- 15-041-11N** - Ea.  
Total Length: 425mm; High Range: +150 deg. C; Low Range: -20 deg. C
- 15-041-11K** - Ea.  
Total Length: 425mm; High Range: +250 deg. C; Low Range: -10 deg. C
- 15-041-11M** - Ea.  
Total Length: 450mm; High Range: +250 deg. C; Low Range: -10 deg. C
- 15-041-11W** - Ea.  
Total Length: 465mm; High Range: +150 deg. C; Low Range: -20 deg. C

- 15-041-11P** - Ea.  
Total Length: 475mm; High Range: +250 deg. C; Low Range: -10 deg. C
- 15-041-11S** - Ea.  
Total Length: 490mm; High Range: +250 deg. C; Low Range: -10 deg. C
- 15-041-11V** - Ea.  
Total Length: 500mm; High Range: +250 deg. C; Low Range: -10 deg. C
- 15-041-11Z** - Ea.  
Total Length: 525mm; High Range: +250 deg. C; Low Range: -10 deg. C
- 15-041-12E** - Ea.  
Total Length: 550mm; High Range: +150 deg. C; Low Range: -20 deg. C
- 15-041-12F** - Ea.  
Total Length: 600mm; High Range: +250 deg. C; Low Range: -10 deg. C

## Fisherbrand Red-Spirit Deep-Immersion Thermometers



- 15-041-6C** - Ea.  
Thermometer; Fisherbrand; Deep-Immersion; Red spirit-filled; High Range: +150 deg. C; Length: 380mm; Low Range: -20 deg. C
- 15-041-6A** - Ea.  
Length: 380mm; High Range: +50 deg. C; Length: 380mm; Low Range: -100 deg. C
- 15-041-6B** - Ea.  
Length: 380mm; High Range: +50 deg. C; Length: 380mm; Low Range: -50 deg. C
- 15-041-7C** - Ea.  
Length: 457mm; High Range: +150 deg. C; Length: 457mm; Low Range: -20 deg. C
- 15-041-6D** - Ea.  
Length: 457mm; High Range: +250 deg. C; Length: 457mm; Low Range: -10 deg. C
- 15-041-7A** - Ea.  
Length: 457mm; High Range: +50 deg. C; Length: 457mm; Low Range: -100 deg. C
- 15-041-7B** - Ea.  
Length: 475mm; High Range: +50 deg. C; Length: 457mm; Low Range: -50 deg. C
- 15-041-8C** - Ea.  
Length: 533mm; High Range: +150 deg. C; Length: 533mm; Low Range: -20 deg. C
- 15-041-7D** - Ea.  
Length: 533mm; High Range: +250 deg. C; Length: 533mm; Low Range: -10 deg. C
- 15-041-8A** - Ea.  
Length: 533mm; High Range: +50 deg. C; Length: 533mm; Low Range: -100 deg. C
- 15-041-8B** - Ea.  
Length: 533mm; High Range: +50 deg. C; Length: 533mm; Low Range: -50 deg. C
- 15-041-8D** - Ea.  
Length: 610mm; High Range: +250 deg. C; Length: 610mm; Low Range: -10 deg. C
- 15-041-9A** - Ea.  
Length: 685mm; High Range: +150 deg. C; Length: 685mm; Low Range: -20 deg. C
- 15-041-9B** - Ea.  
Length: 762mm; High Range: +250 deg. C; Length: 762mm; Low Range: -10 deg. C

# Timers

## Fisherbrand Countdown Timers



- 14-649-17** - Ea.  
Timer, Countdown; Traceable; Fisherbrand; Accuracy: 0.01%; Timer Range: 4 channels; 99 hr., 59 min., 59 sec. to 1 sec.; Display: 6-digit, 3/4 in. H LCD; Dimensions: 2-1/8L x 2-3/8W x 1/2 in. H
- 06-662-51** - Ea.  
Timer Range: 99 min., 59 sec. to 1 sec.; Dimensions: 1/2L x 2-1/2W x 2-1/2 in. H
- 06-662-44** - Ea.  
Timer Range: Two channels; 23 hours, 59 minutes, 59 seconds to 1 second; Display: Double-line, each line 3/8 in.H (1.0cm) LCD; Dimensions: 7/8L x 2 3/4W x 2 1/8 in.H



# Tubes

## Microcentrifuge Tubes

### Fisherbrand Locking-Lid Microcentrifuge Tubes with Polypropylene Snap-Cap



- 02-681-275** - Pk. of 1000, Cs. of 10  
Tube, Microcentrifuge; Fisherbrand; Polypropylene; w/ attached locking cap; Graduated; Capacity 0.6mL; Max. RCF 30,000xG; Blue
- 02-681-276** - Pk. of 1000, Cs. of 10  
Capacity 0.6mL; Max. RCF 30,000xG; Green
- 02-681-273** - Pk. of 1000, Cs. of 10  
Capacity 0.6mL; Max. RCF 30,000xG; Natural
- 02-681-278** - Pk. of 1000, Cs. of 10  
Capacity 0.6mL; Max. RCF 30,000xG; Orange
- 02-681-277** - Pk. of 1000, Cs. of 10  
Capacity 0.6mL; Max. RCF 30,000xG; Yellow
- 02-681-281** - Pk. of 1000, Cs. of 10  
Capacity: 0.6mL; Assorted; Max. RCF: 30,000xG
- 02-681-274** - Pk. of 1000, Cs. of 10  
Red
- 02-681-279** - Pk. of 1000, Cs. of 10  
Violet

### Fisherbrand Premium Microcentrifuge Tubes



- 05-408-125** - Pk. of 500, Cs. of 10  
Tube, Microcentrifuge; Fisherbrand; PP; w/Flat-top snap caps w/insertion spot; Graduated; Labeling area; Autoclavable; Nonpyrogenic; RNase/DNase-free; Max. RCF 30,000xG; 0.5mL capacity; 7.5 O.D. x 30.5mm L; Amber
- 05-408-123** - Pk. of 500, Cs. of 10  
0.5mL; Blue; O.D. x L: 7.5 x 30.5mm
- 05-408-124** - Pk. of 500, Cs. of 10  
0.5mL; Green; O.D. x L: 7.5 x 30.5mm
- 05-408-127** - Pk. of 500, Cs. of 10  
0.5mL; Lavender; O.D. x L: 7.5 x 30.5mm
- 05-408-120** - Pk. of 500, Cs. of 10  
0.5mL; Natural; O.D. x L: 7.5 x 30.5mm
- 05-408-126** - Pk. of 500, Cs. of 10  
0.5mL; Orange; O.D. x L: 7.5 x 30.5mm
- 05-408-121** - Pk. of 500, Cs. of 10  
0.5mL; Red; O.D. x L: 7.5 x 30.5mm
- 05-408-122** - Pk. of 500, Cs. of 10  
0.5mL; Yellow; O.D. x L: 7.5 x 30.5mm
- 05-408-128** - Pk. of 500, Cs. of 10  
0.5mL; Assorted; O.D. x L: 7.5 x 30.5mm
- 05-408-134** - Pk. of 500, Cs. of 10  
1.5mL; Amber; O.D. x L: 10.8 x 40.6mm
- 05-408-132** - Pk. of 500, Cs. of 10  
1.5mL; Blue; O.D. x L: 10.8 x 40.6mm
- 05-408-133** - Pk. of 500, Cs. of 10  
1.5mL; Green; O.D. x L: 10.8 x 40.6mm
- 05-408-136** - Pk. of 500, Cs. of 10  
1.5mL; Lavender; O.D. x L: 10.8 x 40.6mm
- 05-408-129** - Pk. of 500, Cs. of 10  
1.5mL; Natural; O.D. x L: 10.8 x 40.6mm
- 05-408-135** - Pk. of 500, Cs. of 10  
1.5mL; Orange; O.D. x L: 10.8 x 40.6mm
- 05-408-130** - Pk. of 500, Cs. of 10  
1.5mL; Red; O.D. x L: 10.8 x 40.6mm
- 05-408-131** - Pk. of 500, Cs. of 10  
1.5mL; Yellow; O.D. x L: 10.8 x 40.6mm
- 05-408-137** - Pk. of 500, Cs. of 10  
1.5mL; Assorted; O.D. x L: 10.8 x 40.6mm
- 05-408-145** - Pk. of 500, Cs. of 10  
2.0mL; Amber; O.D. x L: 11 x 40.6mm
- 05-408-141** - Pk. of 500, Cs. of 10  
2.0mL; Blue; O.D. x L: 11 x 40.6mm
- 05-408-138** - Pk. of 500, Cs. of 10  
2.0mL; Natural; O.D. x L: 11 x 40.6mm
- 05-408-139** - Pk. of 500, Cs. of 10  
2.0mL; Red; O.D. x L: 11 x 40.6mm
- 05-408-140** - Pk. of 500, Cs. of 10  
2.0mL; Yellow; O.D. x L: 11 x 40.6mm
- 05-408-146** - Pk. of 500, Cs. of 10  
2.0mL; Assorted; O.D. x L: 11 x 40.6mm

## Centrifuge Tubes

### Fisherbrand Easy Reader Centrifuge Tube



- 05-539-9** - Cs. of 500  
Tubes, Centrifuge; Easy Reader; Fisherbrand; Black printed graduation in 0.5mL subdivisions from 2 to 14.5mL; Conical bottom; NonSterile; Polypropylene; 50mL; 28 O.D. x 115mm Long; 6000 x G Max. RCF; Plug seal closure; Sleeve
- 05-539-1** - Cs. of 500  
Sterile; Polyethylene Terephthalate; Closure: Plug Seal; Type: Sterile; O.D. x L: 17 x 119mm
- 05-539-10** - Cs. of 500  
Sterile; Polyethylene; Closure: Plug Seal; Type: Sterile; O.D. x L: 28 x 115mm
- 05-539-5** - Cs. of 500  
Sterile; Polypropylene; 15mL; 17 O.D. x 119mm Long; 6000 x G Max. RCF; Plug seal closure; Foam Rack; Case of 500
- 05-539-8** - Cs. of 500  
Sterile; Polypropylene; 50mL; 28 O.D. x 115mm Long; 6000 x G Max. RCF; Flat top; Rack
- 05-539-6** - Cs. of 500  
Sterile; Polypropylene; 50mL; 28 O.D. x 115mm Long; 6000xG Max. RCF; Plug Seal Closure; Rack
- 05-539-7** - Cs. of 500  
Sterile; Polypropylene; 50mL; 28 O.D. x 115mm Long; 6000xG Max. RCF; Plug Seal Closure; Sleeve
- 05-539-12** - Cs. of 500  
Sterile; Polypropylene; Closure: Flat Top; Type: Sterile; O.D. x L: 17 x 119mm
- 05-539-13** - Cs. of 500  
Sterile; Polypropylene; Closure: Flat Top; Type: Sterile; O.D. x L: 28 x 115mm; 25/Sleeve

### Fisherbrand Graduated Conical-Bottom Tubes



- S50712** - Pk. of 50, Cs. of 10  
Tube, Centrifuge; Conical-Bottom; Blue Graduations; Sterile; Polypropylene; Racked; 15mL; Leak and Chemical resistant; Plug Seal; Withstand RFCs up to 6000xG

### Fisherbrand Higher-Speed Easy Reader Plastic Centrifuge Tubes



- 06-443-18** - Cs. of 500  
Tubes, Centrifuge; Easy Reader; Fisherbrand; Polypropylene; 50mL tubes; Graduated in 5mL increments from 5 to 50mL; Flat Top Closure; 28 O.D. x 115mm L; Max. RCF: 9400xG; Sterile; 25/Rack
- 06-443-19** - Cs. of 500  
Flat Top Closure; 25/Sleeve
- 06-443-20PK** - Pk. of 25  
Plug Seal Closure; 25/Pack
- 06-443-20** - Cs. of 500  
Plug Seal Closure; 25/Rack
- 06-443-21** - Cs. of 500  
Plug Seal Closure; 25/Sleeve

## Culture Tubes

### Fisherbrand Disposable Borosilicate Glass Tubes with Plain End



- 14-961-25** - Cs. of 1000  
Tube, Cell Culture; Disposable; Borosilicate Glass; Fisherbrand; Round Bottom; Plain End; O.D. x L: 10 x 75mm
- 14-961-26** - Cs. of 1000  
O.D. x L: 12 x 75mm
- 14-961-27** - Cs. of 1000  
O.D. x L: 13 x 100mm
- 14-961-28** - Cs. of 1000  
O.D. x L: 15 x 85mm
- 14-961-29** - Cs. of 1000  
O.D. x L: 16 x 100mm

- 14-961-31** - Cs. of 1000  
O.D. x L: 16 x 150mm
- 14-961-30** - Cs. of 1000  
O.D. x L: 16 x 25 mm
- 14-961-32** - Cs. of 500  
O.D. x L: 18 x 150mm
- 14-961-33** - Cs. of 500  
O.D. x L: 20 x 150mm
- 14-961-34** - Cs. of 500  
O.D. x L: 25 x 150mm

## Hybridization Tubes

### FisherBiotech Hybridization Incubator Tubes

- 13-247-100** - Pk. of 2  
Tubes, Incubator; FisherBiotech; Hybridization; For rotisseries that accept standard 35x70mm dia. tubes; Thick-walled borosilicate glass; With polypropylene cap, silicone O-ring; Size: 35 x 100mm
- 13-247-150** - Pk. of 2, Cs. of 3  
35mm ID x 150mmL
- 13-247-300** - Pk. of 2, Cs. of 3  
35mm ID x 300mmL
- 13-247-70** - Ea.  
70mm ID x 300mmL


## Test Tubes - Disposable Glass

### Fisherbrand Disposable Borosilicate Glass Tubes With Polypropylene Screw Cap



- 14-962-26D** - Cs. of 1000  
Test Tubes; Disposable; Borosilicate Glass; Round Bottom; Fisherbrand; With polypropylene screw cap; With marking spot; O.D. x L: 13 x 100mm; No. Per Tray: 500
- 14-962-26F** - Cs. of 1000  
O.D. x L: 16 x 100mm
- 14-962-26G** - Cs. of 1000  
O.D. x L: 16 x 125mm
- 14-962-26H** - Cs. of 1000  
O.D. x L: 16 x 150mm
- 14-962-26K** - Cs. of 500  
O.D. x L: 20 x 125mm
- 14-962-26M** - Cs. of 500  
O.D. x L: 20 x 150mm

### Fisherbrand Disposable Borosilicate Glass Tubes With Threaded End




- 14-959-35C** - Cs. of 1000  
Tubes, Test; Disposable; Borosilicate Glass; Fisherbrand; Threaded End; Round bottom; Closure: none; Style: With marking spot; O.D. x L: 13 x 100mm; GPI Thread: 13-415
- 14-959-35AA** - Cs. of 1000  
With marking spot; O.D. x L: 16 x 100mm
- 14-959-35A** - Cs. of 1000  
With marking spot; O.D. x L: 16 x 125mm
- 14-959-35B** - Cs. of 1000  
With marking spot; O.D. x L: 16 x 150mm
- 14-959-37A** - Cs. of 500  
With marking spot; O.D. x L: 20 x 125mm
- 14-959-37C** - Cs. of 500  
With marking spot; O.D. x L: 20 x 150mm
- 14-957-76A** - Cs. of 1000  
Without marking spot; O.D. x L: 13 x 100mm
- 14-957-76C** - Cs. of 1000  
Without marking spot; O.D. x L: 16 x 150mm
- 14-957-76E** - Cs. of 500  
Without marking spot; O.D. x L: 20 x 150mm

### Fisherbrand Disposable Flint Glass Tubes with Plain End




- 14-958A** - Cs. of 1000  
Tubes, Test; Disposable; Flint Glass; Plain End; Fisherbrand; O.D. x L: 6 x 50mm; 500/Tray
- 14-958B** - Cs. of 1000  
O.D. x L: 10 x 75mm; 250/tray
- 14-958J** - Cs. of 500  
O.D. x L: 12 x 75mm; 125/tray
- 14-958C** - Cs. of 1000  
O.D. x L: 12 x 75mm; 250/tray
- 14-958D** - Cs. of 1000  
O.D. x L: 13 x 100mm; 250/tray
- 14-958E** - Cs. of 1000  
O.D. x L: 15 x 85mm; 250/tray
- 14-958F** - Cs. of 1000  
O.D. x L: 16 x 100mm; 250/tray
- 14-958G** - Cs. of 1000  
O.D. x L: 16 x 125mm; 250/tray
- 14-958H** - Cs. of 1000  
O.D. x L: 16 x 150mm; 250/tray
- 14-958K** - Cs. of 500  
O.D. x L: 20 x 150mm; 250/tray

### Fisherbrand Disposable Heavy-Wall Borosilicate Glass Tubes



- 14-958-10A** - Cs. of 1000  
Test Tubes; Disposable; Borosilicate Glass; Heavy Wall; Round Bottom; Plain End; Fisherbrand; O.D. x L: 10 x 75mm; packaged in lint-free trays
- 14-958-10B** - Cs. of 1000  
O.D. x L: 12 x 75mm
- 14-958-10C** - Cs. of 1000  
O.D. x L: 13 x 100mm

### Fisherbrand Disposable Plain-End Glass Tubes with Color-Coded Label Area



- 14-957-20E** - Cs. of 1000  
Test Tube; Disposable Glass; Plain end; Round Bottom; Fisherbrand; With color-coded label area; O.D. x L: 10 x 75mm; Borosilicate Glass; Label Area: Blue
- 14-957-22E** - Cs. of 1000  
Borosilicate Glass; Label Color: Blue; O.D. x L: 12 x 75mm
- 14-957-24E** - Cs. of 1000  
Borosilicate Glass; Label Color: Blue; O.D. x L: 13 x 100mm
- 14-957-26E** - Cs. of 1000  
Borosilicate Glass; Label Color: Blue; O.D. x L: 16 x 100mm
- 14-957-20C** - Cs. of 1000  
Borosilicate Glass; Label Color: Green; O.D. x L: 10 x 75mm
- 14-957-22C** - Cs. of 1000  
Borosilicate Glass; Label Color: Green; O.D. x L: 12 x 75mm
- 14-957-24C** - Cs. of 1000  
Borosilicate Glass; Label Color: Green; O.D. x L: 13 x 100mm
- 14-957-26C** - Cs. of 1000  
Borosilicate Glass; Label Color: Green; O.D. x L: 16 x 100mm
- 14-957-20D** - Cs. of 1000  
Borosilicate Glass; Label Color: Orange; O.D. x L: 10 x 75mm
- 14-957-22D** - Cs. of 1000  
Borosilicate Glass; Label Color: Orange; O.D. x L: 12 x 75mm
- 14-957-24D** - Cs. of 1000  
Borosilicate Glass; Label Color: Orange; O.D. x L: 13 x 100mm
- 14-957-26D** - Cs. of 1000  
Borosilicate Glass; Label Color: Orange; O.D. x L: 16 x 100mm
- 14-957-20A** - Cs. of 1000  
Borosilicate Glass; Label Color: White; O.D. x L: 10 x 75mm
- 14-957-22A** - Cs. of 1000  
Borosilicate Glass; Label Color: White; O.D. x L: 12 x 75mm

**14-957-24A** - Cs. of 1000  
Borosilicate Glass; Label Color: White; O.D. x L: 13 x 100mm

**14-957-26A** - Cs. of 1000  
Borosilicate Glass; Label Color: White; O.D. x L: 16 x 100mm

**14-957-20B** - Cs. of 1000  
Borosilicate Glass; Label Color: Yellow; O.D. x L: 10 x 75mm

**14-957-22B** - Cs. of 1000  
Borosilicate Glass; Label Color: Yellow; O.D. x L: 12 x 75mm

**14-957-24B** - Cs. of 1000  
Borosilicate Glass; Label Color: Yellow; O.D. x L: 13 x 100mm

**14-957-26B** - Cs. of 1000  
Borosilicate Glass; Label Color: Yellow; O.D. x L: 16 x 100mm

**14-957-10E** - Cs. of 1000  
Flint Glass; Label Color: Blue; O.D. x L: 10 x 75mm

**14-957-12E** - Cs. of 1000  
Flint Glass; Label Color: Blue; O.D. x L: 12 x 75mm

**14-957-14E** - Cs. of 1000  
Flint Glass; Label Color: Blue; O.D. x L: 13 x 100mm

**14-957-16E** - Cs. of 1000  
Flint Glass; Label Color: Blue; O.D. x L: 16 x 100mm

**14-957-18E** - Cs. of 1000  
Flint Glass; Label Color: Blue; O.D. x L: 16 x 125mm

**14-957-10C** - Cs. of 1000  
Flint Glass; Label Color: Green; O.D. x L: 10 x 75mm

**14-957-12C** - Cs. of 1000  
Flint Glass; Label Color: Green; O.D. x L: 12 x 75mm

**14-957-14C** - Cs. of 1000  
Flint Glass; Label Color: Green; O.D. x L: 13 x 100mm

**14-957-16C** - Cs. of 1000  
Flint Glass; Label Color: Green; O.D. x L: 16 x 100mm

**14-957-10D** - Cs. of 1000  
Flint Glass; Label Color: Orange; O.D. x L: 10 x 75mm

**14-957-12D** - Cs. of 1000  
Flint Glass; Label Color: Orange; O.D. x L: 12 x 75mm

**14-957-14D** - Cs. of 1000  
Flint Glass; Label Color: Orange; O.D. x L: 13 x 100mm

**14-957-16D** - Cs. of 1000  
Flint Glass; Label Color: Orange; O.D. x L: 16 x 100mm

**14-957-18D** - Cs. of 1000  
Flint Glass; Label Color: Orange; O.D. x L: 16 x 125mm

**14-957-10A** - Cs. of 1000  
Flint Glass; Label Color: White; O.D. x L: 10 x 75mm

**14-957-12A** - Cs. of 1000  
Flint Glass; Label Color: White; O.D. x L: 12 x 75mm

**14-957-14A** - Cs. of 1000  
Flint Glass; Label Color: White; O.D. x L: 13 x 100mm

**14-957-16A** - Cs. of 1000  
Flint Glass; Label Color: White; O.D. x L: 16 x 100mm

**14-957-12B** - Cs. of 1000  
Flint Glass; Label Color: Yellow; O.D. x L: 12 x 75mm

**14-957-10B** - Cs. of 1000  
Flint Glass; Label Color: Yellow; O.D. x L: 10 x 75mm

**14-957-14B** - Cs. of 1000  
Flint Glass; Label Color: Yellow; O.D. x L: 13 x 100mm

**14-957-16B** - Cs. of 1000  
Flint Glass; Label Color: Yellow; O.D. x L: 16 x 100mm

**14-956-8J** - Cs. of 1000  
Clear polystyrene; 17 x 100mm; Cap: None

**14-959-16** - Pk. of 1000, Cs. of 4  
Translucent polypropylene; 12 x 75mm; Cap: None; 1000/Pk.

**14-961-11** - Cs. of 5000  
Translucent polypropylene; 12 x 75mm; Cap: None; 5000/Cs.

**14-959-10AA** - Cs. of 1000  
Translucent polypropylene; 12 x 75mm; Cap: Snap

**14-956-7A** - Cs. of 1000  
Translucent polypropylene; 13 x 100mm; Cap: None

**14-959-16A** - Cs. of 1000  
Translucent polypropylene; 13 x 100mm; Cap: None; Uniformly-oriented packaging

**14-956-7C** - Cs. of 1000  
Translucent polypropylene; 16 x 125mm; Cap: None

**14-956-7E** - Cs. of 1000  
Translucent polypropylene; 17 x 100mm; Cap: None

### Fisherbrand Sterile Plastic Tubes



**14-956-3B** - Cs. of 500  
Test Tubes; Disposable Fisherbrand; Sterile; Ideal for tissue culture procedures; No residual contamination; Clear Polystyrene; 12 x 75mm; Cap: Snap; Individually wrapped

**14-956-3C** - Cs. of 500  
Clear Polystyrene; 12 x 75mm; Cap: Snap; Bag of 25

**14-956-3D** - Cs. of 1000  
Clear Polystyrene; 12 x 75mm; Cap: Snap; Bag of 125

**14-956-4A** - Cs. of 1000  
Clear Polystyrene; 13 x 100mm; Cap: Screw; Blister pack of 125

**14-956-5A** - Cs. of 500  
Clear Polystyrene; 16 x 125mm; Cap: Screw; Individually wrapped

**14-956-5B** - Cs. of 1000  
Clear Polystyrene; 16 x 125mm; Cap: Screw; Bag of 125

**14-956-5C** - Cs. of 500  
Clear Polystyrene; 16 x 150mm; Cap: Screw; Individually wrapped

**14-956-6A** - Cs. of 500  
Clear Polystyrene; 17 x 100mm; Cap: Snap; Individually wrapped

**14-956-6D** - Cs. of 500  
Clear Polystyrene; 17 x 100mm; Cap: Snap; Blister pack of 25

**14-956-6B** - Cs. of 1000  
Clear Polystyrene; 17 x 100mm; Cap: Snap; Blister pack of 125

**14-956-6C** - Cs. of 1000  
Clear Polystyrene; 17 x 100mm; Cap: None; Blister pack of 125

**14-956-1C** - Cs. of 500  
Translucent Polypropylene; 12 x 75mm; Cap: Snap; Individually wrapped

**14-956-1D** - Cs. of 500  
Translucent Polypropylene; 12 x 75mm; Cap: None; Blister pack of 25

**14-956-1B** - Cs. of 1000  
Translucent Polypropylene; 12 x 75mm; Cap: None; Blister pack of 125

**14-956-1J** - Cs. of 500  
Translucent Polypropylene; 17 x 100mm; Cap: Snap; Individually wrapped

**14-956-1G** - Cs. of 1000  
Translucent Polypropylene; 17 x 100mm; Cap: None; Individually wrapped

**14-956-3A** - Pk. of 125, Cs. of 8  
Translucent Polypropylene; 17 x 100mm; Cap: None; Blister pack of 25

**14-956-1H** - Cs. of 500  
Translucent Polypropylene; 17 x 100mm; Cap: None; Blister pack of 125

## Test Tubes - Disposable Plastic

### Fisherbrand Nonsterile Plastic Tubes



**14-961-10** - Pk. of 1000, Cs. of 4  
Tubes, Test; Disposable Plastic; Fisherbrand; Nonsterile; Round Bottom; Extra-heavy wall construction; Clear polystyrene; 12 x 75mm; Cap: None; 1000/Pk.

**14-961-13** - Cs. of 5000  
Clear polystyrene; 12 x 75mm; Cap: None; 5000/Cs.

**14-961-10A** - Cs. of 1000  
Clear polystyrene; 12 x 75mm; Cap: None; Uniformly-oriented packaging

**14-961-20** - Cs. of 1000  
Clear polystyrene; 12 x 75mm; Cap: Snap

**14-956-8E** - Cs. of 1000  
Clear polystyrene; 13 x 100mm; Cap: None

**14-956-8G** - Cs. of 1000  
Clear polystyrene; 16 x 125mm; Cap: None



# Tubing

## Plastic Tubing


### Fisherbrand Clear PVC Tubing



- 14-169-7A** - Pk. of 50  
Tubing, Plastic; Fisherbrand; Clear PVC; I.D.: 1/8 in. (3.2mm); Wall Thickness: 1/16 in.; O.D.: 1/4 in.; In 50 ft. Lengths
- 14-169-7B** - Pk. of 50  
I.D.: 3/16 in.; Wall Thickness: 1/16 in.; O.D.: 5/16 in.
- 14-169-7C** - Pk. of 50  
I.D.: 1/4 in.; Wall Thickness: 1/16 in.; O.D.: 3/8 in.
- 14-169-7D** - Pk. of 50  
I.D.: 1/4 in.; Wall Thickness: 3/32 in.; O.D.: 7/16 in.
- 14-169-7E** - Pk. of 50  
I.D.: 1/4 in.; Wall Thickness: 1/8 in.; O.D.: 1/2 in.
- 14-169-7F** - Pk. of 50  
I.D.: 5/16 in.; Wall Thickness: 1/16 in.; O.D.: 7/16 in.
- 14-169-7G** - Pk. of 50  
I.D.: 3/8 in.; Wall Thickness: 1/16 in.; O.D.: 1/2 in. (12.7mm)
- 14-169-7H** - Pk. of 50  
I.D.: 3/8 in.; Wall Thickness: 1/8 in.; O.D.: 5/8 in. (15.9mm).
- 14-169-7J** - Pk. of 50  
I.D.: 1/2 in.; Wall Thickness: 1/16 in.; O.D.: 5/8 in. (15.9mm)
- 14-169-7K** - Pk. of 50  
I.D.: 1/2 in.; Wall Thickness: 3/32 in.; O.D.: 11/16 in.
- 14-169-7L** - Pk. of 50  
I.D.: 1/2 in.; Wall Thickness: 1/8 in.; O.D.: 3/4 in. (19mm)
- 14-169-7M** - Pk. of 50  
I.D.: 5/8 in.; Wall Thickness: 1/8 in.; O.D.: 7/8 in.
- 14-169-7N** - Pk. of 50  
I.D.: 3/4 in.; Wall Thickness: 1/8 in.; O.D.: 1 in.

## Rubber Tubing

### Fisherbrand Natural Rubber Latex 50 ft. Tubing




- 14-178-2AA** - Ea.  
Tubing, Pure Natural Rubber; Description: Amber, Thick Wall; Talc free; Nontacky; Autoclavable; I.D.: 1/8 in. (3.2mm); Length: 50 ft.
- 14-178-2BB** - Ea.  
Thick wall; I.D.: 1/8 in. (4.8mm)
- 14-178-2CC** - Ea.  
Thick wall; I.D.: 1/4 in. (6.4mm)
- 14-178-2B** - Ea.  
Thin wall; I.D.: 3/16 in. (4.8mm)
- 14-178-2C** - Ea.  
Thin wall; I.D.: 1/4 in. (6.4mm)
- 14-178-2D** - Ea.  
Thin wall; I.D.: 5/16 in. (7.9mm)

# Vials

## Autosampler Vials


### Fisherbrand Glass Shell Vials



- 03-339-26A** - Pk. of 144, Cs. of 16  
Vial, Shell; Borosilicate Glass, Type I, N-51A; Fisherbrand; Straight-sided; Ideal for dry contents; Closure: Titeseal; Package: Bulk; O.D. x H: 12 x 35mm; Capacity: 1/2 dr.
- 03-339-26B** - Pk. of 144, Cs. of 16  
Closure: Titeseal; Package: Bulk pack; O.D. x H: 15 x 45mm; Capacity: 1 dr.
- 03-339-26C** - Pk. of 144, Cs. of 12  
Closure: Titeseal; Package: Bulk pack; O.D. x H: 19 x 51mm; Capacity: 2 dr.
- 03-339-26D** - Pk. of 144, Cs. of 8  
Closure: Titeseal; Package: Bulk pack; O.D. x H: 19 x 65mm; Capacity: 3 dr.
- 03-339-26F** - Pk. of 144, Cs. of 3  
Closure: Titeseal; Package: Bulk; O.D. x H: 29 x 65mm; Capacity: 7 dr.
- 03-339-26G** - Pk. of 144, Cs. of 3  
Closure: Titeseal; Package: Bulk; O.D. x H: 29 x 80mm; Capacity: 10 dr.
- 03-339-26H** - Pk. of 144, Cs. of 3  
Closure: Titeseal; Package: Bulk; O.D. x H: 29 x 94mm; Capacity: 12 dr.
- 03-339-30A** - Pk. of 144, Cs. of 10  
Closure: None, accept corks; Package: Lab pack; O.D. x H: 9 x 30mm; Capacity: 1/4 dr.
- 03-339-30B** - Pk. of 144, Cs. of 10  
Closure: None, accept corks; Package: Lab pack; O.D. x H: 12 x 35mm; Capacity: 1/2 dr.
- 03-339-30C** - Pk. of 144, Cs. of 10  
Closure: None, accept corks; Package: Lab pack; O.D. x H: 15 x 45mm; Capacity: 1 dr.
- 03-339-30E** - Pk. of 144, Cs. of 6  
Closure: None, accept corks; Package: Lab pack; O.D. x H: 15 x 60mm; Capacity: 2 dr.
- 03-339-30F** - Pk. of 144, Cs. of 6  
Closure: None, accept corks; Package: Lab pack; O.D. x H: 19 x 65mm; Capacity: 3 dr.
- 03-339-30G** - Pk. of 144, Cs. of 4  
Closure: None, accept corks; Package: Lab Pack; O.D. x H: 21 x 70mm; Capacity: 4 dr.
- 03-339-30H** - Pk. of 144, Cs. of 4  
Closure: None, accept corks; Package: Lab Pack; O.D. x H: 23 x 75mm; Capacity: 6 dr.
- 03-339-30J** - Pk. of 144, Cs. of 4  
Closure: None, accept corks; Package: Lab Pack; O.D. x H: 25 x 95mm; Capacity: 8 dr.

## Cryogenic Storage Vials

### Fisherbrand Externally and Internally Threaded Cryogenic Storage Vials



- 10-500-25** - Pk. of 100, Cs. of 10  
Vial, Cryogenic Storage; Fisherbrand; Externally threaded; 1.2mL; Polypropylene; Self standing; HDPE closure; Graduations in 1/2mL increments; Large White Marking Area
- 10-500-26** - Pk. of 100, Cs. of 10  
Externally Threaded, 2.0mL
- 10-500-27** - Pk. of 50, Cs. of 20  
Externally Threaded, 5.0mL
- 12-567-500** - Pk. of 100, Cs. of 10  
Internally Threaded with Silicone Gasket; 1.2mL; Conical-Style Bottom with a Star-Style Foot
- 12-567-501** - Pk. of 100, Cs. of 10  
Internally Threaded with Silicone Gasket; 2.0mL; Round-Style Bottom with a Star-Style Foot
- 12-567-502** - Pk. of 100, Cs. of 10  
Internally Threaded with Silicone Gasket; 5.0mL; Round-Style Bottom with a Star-Style Foot

## EPA Vials

### Fisherbrand Economical Glass EPA Vials



**03-340-122** - Cs. of 72  
Vials, EPA; Fisherbrand; Capacity: 20mL; O.D. x H: 27 x 57mm; Color: Amber; Type I Borosilicate Glass; Unassembled; White polypropylene open-top closures with 0.125 in. TFE/silicone septa

**03-340-121** - Cs. of 72  
20mL; O.D. x H: 28 x 57mm; Color: Clear

**03-340-120** - Cs. of 72  
40mL; O.D. x H: 27 x 95mm; Color: Amber

**03-340-119** - Cs. of 72  
40mL; O.D. x H: 27 x 95mm; Color: Clear

### Fisherbrand EnviroWare Assembled Borosilicate Glass EPA Vials



**03-339-14C** - Pk. of 72, Cs. of 6  
Vial, EPA; Fisherbrand; Enviroware; Borosilicate Glass; GPI Thread: 24-400; Type: Clear; O.D. x H: 28 x 60mm; Capacity: 20mL

**03-339-14E** - Pk. of 72, Cs. of 6  
20mL; Amber

**03-339-14A** - Pk. of 72, Cs. of 6  
40mL; Clear

### Fisherbrand Glass EPA Vials



**02-912-382** - Cs. of 144  
Vial, EPA; Borosilicate Glass; Fisherbrand; Closure: Open-top Polypropylene screw caps, with 0.125 in. PTFE/silicone septa; Pre-cleaned;

**02-912-380** - Cs. of 144  
Type: Amber glass; Capacity: 10 dr. (40mL); O.D. x L: 28 x 95mm

**02-912-377** - Cs. of 144  
Type: Certified amber glass; Capacity: 10 dr. (40mL); O.D. x L: 28 x 95mm

**02-912-381** - Cs. of 144  
Type: Clear glass; Capacity: 5 dr. (20mL); O.D. x L: 28 x 57mm

**02-912-379** - Cs. of 144  
Type: Clear glass; Capacity: 10 dr. (40mL); O.D. x L: 28 x 95mm

**02-912-378** - Cs. of 144  
Type: Certified clear glass; Capacity: 5 dr. (20mL); O.D. x L: 28 x 57mm

**02-912-376** - Cs. of 144  
Type: Certified clear glass; Capacity: 10 dr. (40mL); O.D. x L: 28 x 95mm

## Sample Vials - Glass

### Fisherbrand Amber Glass Threaded Vials with Rubber-Lined Closures



**03-339-23A** - Pk. of 144, Cs. of 16  
Vial, General-purpose; Type I, Class B Amber Borosilicate Glass; Fisherbrand; GPI Thread Finish: 8-425; Closure: Black screw-top cap, White rubber liner; Autoclavable; Capacity: 1/2 dr. (1.8mL); O.D. x H: 12 x 35mm

**03-339-23B** - Pk. of 144, Cs. of 16  
Capacity: 1 dr. (3.7mL); O.D. x H: 15 x 45mm

**03-339-23C** - Pk. of 144, Cs. of 12  
Capacity: 2 dr. (7.4mL); O.D. x H: 17 x 60mm

**03-339-23D** - Pk. of 144, Cs. of 8  
Capacity: 4 dr. (14.8mL); O.D. x H: 21 x 70mm

**03-339-23E** - Pk. of 144, Cs. of 4  
Capacity: 8 dr. (29.6mL); O.D. x H: 25 x 95mm

### Fisherbrand Class A Clear Glass Threaded Vials with Attached Caps



**02-912-368** - Cs. of 144  
Vial, Sample; Type I, Class A Borosilicate glass; Fisherbrand; Closure: Black phenolic screw cap; Liner: Low-density polyethylene Poly-Seal cone; Capacity: 1 dr. (4mL); O.D. x L: 15 x 45mm; GPI Thread: 13-425

**02-912-369** - Cs. of 144  
Liner: Low-density polyethylene Poly-Seal cone; Capacity: 2 dr. (8mL); O.D. x L: 17 x 60mm; GPI Thread: 15-425

**02-912-374** - Cs. of 72  
Liner: Low-density polyethylene Poly-Seal cone; Capacity: 6-1/4 dr. (25mL); O.D. x L: 28 x 70mm; GPI Thread: 24-400

**02-912-375** - Cs. of 144  
Liner: Low-density polyethylene Poly-Seal cone; Capacity: 8 dr. (30mL); O.D. x L: 25 x 95mm; GPI Thread: 22-400

**02-912-351** - Cs. of 288  
Liner: PTFE-faced white 14B rubber; Capacity: 1/2 dr. (2mL); O.D. x L: 12 x 35mm; GPI Thread: 8-425

**02-912-352** - Cs. of 144  
Liner: PTFE-faced white 14B rubber; Capacity: 1 dr. (4mL); O.D. x L: 15 x 45mm; GPI Thread: 13-425

**02-912-353** - Cs. of 144  
Liner: PTFE-faced white 14B rubber; Capacity: 2 dr. (8mL); O.D. x L: 17 x 60mm; GPI Thread: 15-425

**02-912-354** - Cs. of 144  
Liner: PTFE-faced white 14B rubber; Capacity: 3 dr. (12mL); O.D. x L: 19 x 65mm; GPI Thread: 15-425

**02-912-355** - Cs. of 144  
Liner: PTFE-faced white 14B rubber; Capacity: 4 dr. (16mL); O.D. x L: 21 x 70mm; GPI Thread: 18-400

**02-912-357** - Cs. of 72  
Liner: PTFE-faced white 14B rubber; Capacity: 5 dr. (20mL); O.D. x L: 28 x 57mm; GPI Thread: 24-400

**02-912-356** - Cs. of 144  
Liner: PTFE-faced white 14B rubber; Capacity: 6 dr. (24mL); O.D. x L: 23 x 85mm; GPI Thread: 20-400

**02-912-358** - Cs. of 72  
Liner: PTFE-faced white 14B rubber; Capacity: 6-1/4 dr. (24mL); O.D. x L: 28 x 70mm; GPI Thread: 24-400

**02-912-359** - Cs. of 144  
Liner: PTFE-faced white 14B rubber; Capacity: 8 dr. (30mL); O.D. x L: 25 x 95mm; GPI Thread: 22-400

**02-912-370** - Cs. of 144  
PE Poly-seal liner; Capacity: 12mL; O.D. x H: 19 x 65mm

**02-912-371** - Cs. of 144  
PE Poly-seal liner; Capacity: 16mL; O.D. x H: 21 x 70mm

**02-912-373** - Cs. of 72  
PE Poly-seal liner; Capacity: 20mL; O.D. x H: 28 x 27mm

**02-912-372** - Cs. of 144  
PE Poly-seal liner; Capacity: 24mL; O.D. x H: 23 x 85mm

**02-912-365** - Cs. of 200  
PTFE liner; Capacity: 2mL; O.D. x H: 15 x 28mm

**02-912-366** - Cs. of 200  
PTFE liner; Capacity: 4mL; O.D. x H: 17 x 38mm

**02-912-367** - Cs. of 200  
PTFE liner; Capacity: 6mL; O.D. x H: 19 x 40mm

Fisherbrand Class A Clear Glass Threaded Vials with Rubber-Lined Closures



- 03-339-25A** - Pk. of 144, Cs. of 16  
Vial, Sample; Type I, Class A Borosilicate Glass; Fisherbrand; Low alkali content; Closure: Black phenolic screw cap; Liner: White rubber; GPI Thread Finish: 8-425; Capacity: 1/2 dr. (1.8mL); O.D. x L: 12 x 35mm
- 03-339-25B** - Pk. of 144, Cs. of 16  
GPI Thread Finish: 13-425; Capacity: 1 dr.
- 03-339-25C** - Pk. of 144, Cs. of 12  
GPI Thread Finish: 15-425; Capacity: 2 dr.
- 03-339-25D** - Pk. of 144, Cs. of 8  
GPI Thread Finish: 15-425; Capacity: 3 dr.
- 03-339-25E** - Pk. of 144, Cs. of 8  
GPI Thread Finish: 18-400; Capacity: 4 dr.
- 03-339-25F** - Pk. of 144, Cs. of 6  
GPI Thread Finish: 20-400; Capacity: 6 dr.

Fisherbrand Class B Amber Glass Threaded Vials with Attached Caps



- 02-912-360** - Cs. of 288  
General-purpose, Amber glass, Class B, Threaded vials, 1/2 dr. (2mL), 12 dia. x 35mm high, 8-425 GPI thread finish, With attached black phenolic cap with PTFE liner
- 02-912-361** - Cs. of 144  
1 dr. (4mL)
- 02-912-362** - Cs. of 144  
2 dr. (8mL)
- 02-912-363** - Cs. of 72  
5 dr. (20mL)
- 02-912-364** - Cs. of 72  
10 dr. (40mL)

Fisherbrand Class B Clear Glass Threaded Vials with Rubber-Lined Closures



- 03-339-22A** - Pk. of 144, Cs. of 16  
Vial, Sample, Fisherbrand; Class B Clear Borosilicate Glass Threaded Vial; Autoclavable black phenolic cap, white rubber liner; 1.8mL; 12 x 35mm
- 03-339-22B** - Pk. of 144, Cs. of 16  
3.7mL; 15 x 45mm
- 03-339-22C** - Pk. of 144, Cs. of 16  
5.5mL; 16 x 50mm
- 03-339-22D** - Pk. of 144, Cs. of 12  
7.4mL; 17 x 60mm
- 03-339-22E** - Pk. of 144, Cs. of 8  
11.1mL; 19 x 65mm
- 03-339-22F** - Pk. of 144, Cs. of 8  
14.8mL; 21 x 70mm
- 03-339-22G** - Pk. of 144, Cs. of 6  
22.2mL; 23 x 85mm
- 03-339-22H** - Pk. of 144, Cs. of 4  
29.6mL; 25 x 95mm
- 03-339-22J** - Pk. of 144, Cs. of 3  
14.8mL; 28 x 57mm
- 03-339-22K** - Pk. of 144, Cs. of 3  
22.2mL; 28 x 70mm
- 03-339-22M** - Pk. of 144, Cs. of 3  
35.1mL; 28 x 95mm
- 03-339-22N** - Pk. of 144, Cs. of 3  
40.7mL; 28 x 108mm

Fisherbrand Clear Glass Threaded Vials with Closures Attached



- 03-338AA** - Pk. of 144, Cs. of 36  
Vial, Sample; Borosilicate Glass, Type I, Class B; Fisherbrand; With attached caps; Closure: Black phenolic screw cap; Liner: Polyvinyl-faced pulp; GPI Thread Finish: 8-425; Capacity: 1/2 dr. (1.8mL); O.D. x L: 12 x 35mm
- 03-338A** - Pk. of 144, Cs. of 16  
Capacity: 1 dr.; O.D. x L: 15mm x 45mm;
- 03-338B** - Pk. of 144, Cs. of 18  
Capacity: 1-1/2 dr. (5.5mL); O.D. x L: 16 x 50mm
- 03-338C** - Pk. of 144, Cs. of 12  
Capacity: 2 dr. (7.4mL); O.D. x L: 17 x 60mm
- 03-338D** - Pk. of 144, Cs. of 12  
Capacity: 3 dr. (11.1mL); O.D. x L: 19 x 65mm
- 03-338E** - Pk. of 72, Cs. of 10  
Capacity: 4 dr. (14.8mL); O.D. x L: 28 x 57mm
- 03-338F** - Pk. of 144, Cs. of 8  
Capacity: 4 dr. (14.8mL); O.D. x L: 21 x 70mm
- 03-338G** - Pk. of 72, Cs. of 8  
Capacity: 6 dr. (22.2mL); O.D. x L: 23 x 85mm
- 03-338H** - Pk. of 144, Cs. of 6  
Capacity: 6 dr. (22.2mL); O.D. x L: 23mm x 85mm
- 03-338J** - Pk. of 72, Cs. of 8  
Capacity: 8 dr. (29.6mL); O.D. x L: 25 x 95mm
- 03-338K** - Pk. of 72, Cs. of 6  
Capacity: 9-1/2 dr. (35.1mL); O.D. x L: 28 x 95mm
- 03-338L** - Pk. of 72, Cs. of 6  
Capacity: 11 dr. (40.7mL); O.D. x L: 28 x 108mm

Fisherbrand Clear Glass Threaded Vials With Closures Packaged Separately



- 03-339-21A** - Pk. of 144, Cs. of 16  
Vial, Sample; Borosilicate Glass, Type I, Class B; Fisherbrand; Closure: Black phenolic screw cap; Liner: Polyvinyl-faced pulp; GPI Thread Finish: 8-425; Capacity: 1/2 dr. (1.8mL); O.D. x L: 12 x 35mm
- 03-339-21B** - Pk. of 144, Cs. of 16  
Capacity: 1 dr. (3.7mL); O.D. x L: 15 x 45mm
- 03-339-21C** - Pk. of 144, Cs. of 16  
Capacity: 1-1/2 dr. (5.5mL); O.D. x L: 16 x 50mm
- 03-339-21D** - Pk. of 144, Cs. of 12  
Capacity: 2 dr. (7.4mL); 17 x 60mm
- 03-339-21E** - Pk. of 144, Cs. of 8  
Capacity: 3 dr.; GPI Thread Finish: 15-425; O.D. x L: 19mm x 65mm
- 03-339-21F** - Pk. of 144, Cs. of 8  
Capacity: 4 dr. (14.8mL); O.D. x L: 21 x 70mm
- 03-339-21G** - Pk. of 144, Cs. of 6  
GPI Thread Finish: 20-400; Capacity: 6 dr.
- 03-339-21H** - Pk. of 144, Cs. of 4  
GPI Thread Finish: 22-400; Capacity: 8 dr. (29.6mL); O.D. x L: 25 x 95mm
- 03-339-21J** - Pk. of 144, Cs. of 3  
Capacity: 4 dr. (14.8mL); O.D. x L: 28 x 57mm
- 03-339-21K** - Pk. of 144, Cs. of 3  
Capacity: 6 dr. (22.2mL); O.D. x L: 28 x 70mm
- 03-339-21M** - Pk. of 144, Cs. of 3  
GPI Thread Finish: 24-400; Capacity: 9 dr.
- 03-339-21N** - Pk. of 144, Cs. of 3  
Capacity: 11 dr. (40.7mL); O.D. x L: 28 x 108mm



## Fisherbrand Glass Tooled-Neck Vials with Polyethylene Closures



- 03-339-27A** - Pk. of 144, Cs. of 16  
Vial, Sample; Borosilicate, Type 1, N-51A Glass; Fisherbrand; Tooled neck for accurate fit w/closures; Closure: Open-bottom; O.D. x H: 15 x 45mm; Capacity: 1 dr. (3.7mL)
- 03-339-27B** - Pk. of 144, Cs. of 10  
Closure: Open-bottom; O.D. x H: 19 x 48mm; Capacity: 2 dr. (7.4mL)
- 03-339-27C** - Pk. of 144, Cs. of 8  
Closure: Closed-bottom; O.D. x H: 21 x 50mm; Case of: 8; Capacity: 3 dr. (11.1mL)
- 03-339-27D** - Pk. of 144, Cs. of 6  
Closure: Closed-bottom; O.D. x H: 25 x 52mm; Capacity: 4 dr. (14.8mL)
- 03-339-27E** - Pk. of 144, Cs. of 4  
Closure: Closed-bottom; O.D. x H: 27 x 55mm; Case of: 4; Capacity: 5 dr. (18.5mL)
- 03-339-27F** - Pk. of 144, Cs. of 3  
Type 1, N-51A borosilicate glass; Tooled-Neck; PE Closure; 6 dr. (25.9mL)

## Sample Vials - Plastic

### Fisherbrand Polyethylene Sample Vials with Hinged Cap



- 03-338-1A** - Cs. of 100  
Vial, Sample, Fisherbrand; Polyethylene Vial; Withstand temperatures to 80 deg. C (176 deg. F); Snap-top hinged lids; 1/4 dr.
- 03-338-1B** - Cs. of 100  
Polyethylene; Hinged snap cap; 1-1/16 dr. (2.5mL); 14 x 32mm
- 03-338-1C** - Cs. of 100  
Polyethylene; Hinged snap cap; 1-3/8 dr. (5mL); 15 x 50mm
- 03-338-1D** - Cs. of 100  
Polyethylene; Hinged snap cap; 2-1/4 dr. (8mL); 17 x 56mm
- 03-338-1E** - Cs. of 100  
Polyethylene; Hinged snap cap; 2-3/4 dr. (10mL); 22 x 32mm
- 03-338-1F** - Cs. of 100  
Polyethylene; Hinged snap cap; 6-7/8 dr. (25mL); 24 x 72mm
- 03-338-1G** - Cs. of 100  
Polyethylene; Hinged snap cap; 13-3/4 dr. (50mL); 30 x 74mm

### Fisherbrand Polypropylene Sample Vials with Hinged Cap



- 03-341-75A** - Cs. of 100  
Vials, Sample; Fisherbrand; With hinged cap; High-density polypropylene; Flat surface on cap top; Withstand temperatures up to 100 deg. C (212 deg. F); Color: Opaque; Capacity: 1-1/2 oz. (45mL); 100/Cs.
- 03-341-75C** - Cs. of 600  
Color: Opaque; Capacity: 1-1/2 oz. (45mL); 600/Cs.
- 03-341-75E** - Cs. of 100  
Color: Opaque; Capacity: 4 oz. (120mL); 100/Cs.
- 03-341-75G** - Cs. of 300  
Color: Opaque; Capacity: 4 oz. (120mL); 300/Cs.
- 03-341-76A** - Cs. of 50  
Color: Opaque; Capacity: 8 oz. (240mL)

### Fisherbrand Polystyrene Containers



- 03-338-3A** - Cs. of 144  
Vial, Sample, Fisherbrand; Polystyrene Container; Snap-on cap; 3 dr.
- 03-338-3B** - Cs. of 144  
5 dr. (18.5mL)
- 03-338-3C** - Cs. of 144  
7 dr. (25.9mL)
- 03-338-3D** - Cs. of 144  
10 dr. (37.0mL)

- 03-338-3E** - Cs. of 144  
13 dr. (48.1mL)
- 03-338-3F** - Cs. of 144  
15 dr. (55.5mL)
- 03-338-3G** - Cs. of 144  
20 dr. (74mL)
- 03-338-3H** - Cs. of 72  
30 dr. (111mL)
- 03-338-3J** - Cs. of 72  
40 dr. (148mL)

## Scintillation Vials - Glass

### Fisherbrand 7mL Borosilicate Glass Scintillation Vial



- 03-337-26** - Cs. of 1000  
Vials, Scintillation; 7mL Borosilicate Glass; Fisherbrand; White plastic screw cap; GPI Thread: 15-425; O.D. x H: 17 x 54mm

### Fisherbrand 20mL Borosilicate Glass Scintillation Vials



- 03-337-14** - Cs. of 500  
Vials, Scintillation; Fisherbrand; Borosilicate glass; 20mL; Potassium-free; Fit all common counters; White polyethylene; Cap Style: No liner; Packaging: Caps separate
- 03-337-4** - Cs. of 500  
White polypropylene; Cap Style: Pulp-backed metal foil liner
- 03-337-15** - Cs. of 500  
White polypropylene; Cap style: Pulp-backed metal foil liner; Caps separate; 500/CS
- 03-337-7** - Cs. of 500  
White urea; Cap Style: Cone-shaped plastic liner
- 03-337-5** - Cs. of 500  
White urea; Cap Style: Cork-backed metal foil liner; Case of 500
- 03-340-4E** - Cs. of 500  
White urea; Cap Style: Cork-backed metal liner

## Scintillation Vials - Plastic

### Fisherbrand 4, 6, and 8mL HDPE ScintiVials



- 03-337-39** - Cs. of 1000  
Vial, Scintillation; Fisherbrand ScintiVials; HDPE; White polypropylene caps; Cap Style: Flush cap; O.D. x H (with cap): 16 x 56.6mm; Capacity: 6mL
- 03-337-42** - Cs. of 1000  
Cap Style: Flush cap; 8mL
- 03-337-40** - Cs. of 1000  
Cap Style: String of 6 push caps; 4mL
- 03-337-41** - Cs. of 1000  
Cap Style: Twist-off cap with shoulder; 6mL

**Fisherbrand Bulk packed HDPE 7mL Scintillation Vial**



**03-337-20** - Cs. of 2000  
Vial, Scintillation; Fisherbrand; HDPE; Wide mouth;  
Resists swelling; Closure: Screw Cap; Liner: none;  
Capacity: 7mL; Bulk Pack

**Fisherbrand HDPE 7mL Scintillation Vial**



**03-337-1** - Cs. of 1000  
Vial, Scintillation; Fisherbrand; HDPE; Capacity: 7mL;  
Wide mouth; Closure: White polyethylene; Liner: none;  
Resists swelling; GPI Thread: 15-425; O.D. x H: 17 x  
54mm

**Fisherbrand 20mL HDPE Scintillation Vials**



**03-337-23A** - Cs. of 500  
Vials, Scintillation; Fisherbrand; HDPE; 20mL; Resist  
leaking and swelling; Fit all common counters; Provide  
low background counts and excellent light transmission;  
White Polypropylene; Cap Style: Pulp-backed metal foil  
liner; Attached caps  
**03-337-23** - Cs. of 1000  
White Polypropylene; Cap Style: Pulp-backed metal foil  
liner; Bulk pack  
**03-337-23B** - Cs. of 500  
White Polypropylene; Cap Style: Pulp-backed metal foil  
liner; Separate caps  
**03-337-2** - Cs. of 500  
White Urea; Cap Style: Cork-backed metal foil liner  
**03-337-23C** - Cs. of 500  
White Urea; Cone-shaped plastic liner

**Fisherbrand Polystyrene Snap-Cap Vial**



**03-341-11** - Cs. of 144  
Vial, Scintillation; Fisherbrand; 26mL; High-Density  
Polyethylene; Threaded Cap

**Fisherbrand Rack-Mounted Scinti-6 Vials**

**03-337-46** - Cs. of 1000  
Vials, Scintillation; Fisherbrand; 6mL capacity; HDPE; White PP cap; O.D. x  
H (with cap), 16 x 56.6mm. Each tray holds 200 vials in a 10 x 20 array. Case  
contains five rack-mounted trays with five polybags of 200 flush caps each.


**Fisherbrand Scintillation Vials**



**03-337-12B** - Cs. of 500  
Vial, Scintillation; Fisherbrand; 18mL; High-Density  
Polyethylene; Threaded Cap  
**03-337-11B** - Cs. of 500  
Capacity: 20mL

# Weighing Dishes

**Fisherbrand Hexagonal Polystyrene Weighing Dishes**




**02-202B** - Pk. of 500, Cs. of 6  
Dish, Weighing; Fisherbrand; Hexagonal; 2-1/2 in.  
(64mm) Top Opening; 2 in. (51mm) Base; Polystyrene;  
Disposable; 0.05 oz. (1.4g)  
**02-202C** - Pk. of 500, Cs. of 4  
4 in. (102mm) Top Opening; 3 in. (76mm) Base  
**02-202-100** - Pk. of 500, Cs. of 8  
I.D. Top/Base: 1-1/2 in./1 in.; Appx. Weight: 0.02 oz.  
**02-202-103** - Pk. of 500, Cs. of 4  
I.D. Top/Base: 5 in./ 3-1/2 in. Depth: 1-1/4 in.; Appx.  
weight: 0.25 oz.

**Fisherbrand Low Nitrogen Weighing Paper**



**09-898-12A** - Pk. of 500, Cs. of 12  
Paper, Weighing; Fisherbrand; Low nitrogen; Sufficient  
body to support sample; Does not absorb moisture; 3 x 3  
in. (76 x 76mm); 500 sheets per pack  
**09-898-12B** - Pk. of 500, Cs. of 12  
4 x 4 in.  
**09-898-12C** - Pk. of 500, Cs. of 12  
6 x 6 in.

**Fisherbrand Polystyrene Antistatic Weighing Dishes**



**08-732-116** - Pk. of 250, Cs. of 8  
Dish; Weighing; Fisherbrand; Anti-Static; 12mil thick  
Polystyrene; Pour Boat; 1-5/8 L x 1-1/4 W x 5/16 in. H;  
Natural; disposable; flat bottom  
**08-732-117** - Cs. of 250  
Pour Boat; 3-1/2 L x 5-1/4 W x 1 in. H; Natural  
**08-732-118** - Cs. of 250  
Pour Boat; 4-3/4 L x 7-3/4 W x 1 in. H; Natural  
**08-732-112** - Pk. of 500, Cs. of 8  
Standard Style; 1-5/8 L x 1-5/8 W x 5/16 in. H; Natural  
**08-732-114** - Cs. of 500  
Standard style; 3-1/2 L x 3-1/2 W x 1 in. H; Blue  
**08-732-113** - Cs. of 500  
Standard style; 3-1/2 L x 3-1/2 W x 1 in. H; Natural  
**08-732-115** - Cs. of 500  
Standard style; 5-1/2 L x 5-1/2 W x 1 in. H; Blue

# Wipes

**Fisherbrand Premoistened Wipes**



**06-664-14** - Pk. of 100  
Wipes, Disinfectant; Fisherbrand; Dry; Cleanroom-grade  
cellulose/polyester blend cloth; Safe for cleaning optics;  
Will not scratch surfaces or fall apart; Recloseable pop-up  
canister; Size: 6 x 9 in.  
**06-665-23** - Pk. of 100  
Premoistened with Deionized Water  
**06-665-24** - Pk. of 100  
Premoistened with isopropyl alcohol and deionized water





Cat. No.	Page	Cat. No.	Page	Cat. No.	Page	Cat. No.	Page	Cat. No.	Page
01-213	15	02-681-155	25	02-707-436	24	02-911-401	10	02-911-759	5
01-213-1	15	02-681-156	25	02-707-437	24	02-911-402	10	02-911-760	5
01-213-100	15	02-681-157	25	02-707-438	24	02-911-403	10	02-911-761	3
01-213-101	15	02-681-158	25	02-707-440	24	02-911-404	10	02-911-762	3
01-213-102	15	02-681-163	25	02-707-441	24	02-911-405	10	02-911-763	3
01-213-103	15	02-681-164	25	02-707-446	24	02-911-432	6	02-911-764	3
01-213-104	15	02-681-165	25	02-707-447	24	02-911-433	6	02-911-765	3
01-213-105	15	02-681-166	25	02-707-448	24	02-911-434	6	02-911-766	7
01-369A	1	02-681-168	25	02-707-449	24	02-911-435	6	02-911-767	7
01-369B	1	02-681-169	25	02-707-450	24	02-911-436	7	02-911-768	7
01-369D	1	02-681-170	25	02-707-451	24	02-911-437	7	02-911-769	7
01-814-1	1	02-681-171	25	02-707-452	24	02-911-438	7	02-911-770	7
01-814-2	1	02-681-273	33	02-707-453	24	02-911-439	7	02-911-772	3
01-814-3	1	02-681-274	33	02-707-454	24	02-911-440	7	02-911-773	3
01-816-1A	2	02-681-275	33	02-707-455	24	02-911-441	7	02-911-774	3
01-816-1B	2	02-681-276	33	02-707-500	25	02-911-442	7	02-911-775	3
01-816-1C	2	02-681-277	33	02-707-501	25	02-911-443	7	02-911-776	3
01-816-1D	2	02-681-278	33	02-707-502	25	02-911-444	6	02-911-777	5
01-816-1E	2	02-681-279	33	02-707-503	25	02-911-445	6	02-911-778	5
01-816A	2	02-681-281	33	02-707-504	25	02-911-446	6	02-911-779	5
01-816B	2	02-681-436	25	02-707-505	25	02-911-447	6	02-911-780	5
01-816C	2	02-681-437	23	02-707-506	25	02-911-448	6	02-911-781	5
01-816D	2	02-681-438	25	02-707-507	25	02-911-449	6	02-911-782	5
01-816E	2	02-681-440	23	02-707-508	25	02-911-450	7	02-911-783	5
01-816F	2	02-681-454	25	02-707-509	25	02-911-451	7	02-911-784	5
01-826-1	2	02-681-455	25	02-707-510	25	02-911-452	7	02-911-785	5
01-826-2	2	02-681-457	25	02-707-511	25	02-911-453	8	02-911-790	5
01-826-3	2	02-681159	25	02-707-81	23	02-911-454	8	02-911-791	3
01-826-4	2	02-707-100	25	02-707-82	23	02-911-455	8	02-911-793	5
01-826-5	2	02-707-101	25	02-707-83	23	02-911-456	8	02-911-794	5
01-826-6	2	02-707-105	25	02-707-84	23	02-911-457	8	02-911-796	7
01-826-7	2	02-707-106	25	02-707-86	23	02-911-458	8	02-911-798	7
01-826A	2	02-707-107	25	02-707-88	23	02-911-459	8	02-911-799	6
01-826B	2	02-707-108	25	02-707-89	23	02-911-460	8	02-911-812	6
01-826C	2	02-707-110	25	02-895-1A	9	02-911-461	8	02-911-813	6
01-826D	2	02-707-111	25	02-895-1B	9	02-911-462	8	02-911-814	6
01-826E	2	02-707-112	25	02-895-1C	9	02-911-463	8	02-911-815	6
01-827A	1	02-707-113	25	02-895-1D	9	02-911-464	8	02-911-816	6
01-827B	1	02-707-114	25	02-895A	9	02-911-465	8	02-911-817	6
01-827C	1	02-707-115	25	02-895B	9	02-911-466	8	02-911-818	6
01-827D	1	02-707-119	25	02-895C	9	02-911-467	8	02-911-819	6
01-828A	2	02-707-135	23	02-895D	9	02-911-468	8	02-911-820	6
01-828B	2	02-707-136	23	02-896-1A	9	02-911-469	8	02-911-821	6
01-828C	2	02-707-137	23	02-896-1B	9	02-911-470	8	02-911-822	6
01-828D	2	02-707-138	23	02-896-1C	9	02-911-471	8	02-911-823	6
01-828E	2	02-707-139	23	02-896-1D	9	02-911-472	4	02-911-824	3
01-829A	2	02-707-146	23	02-896-1E	9	02-911-473	8	02-911-825	3
01-829B	2	02-707-147	23	02-896-1F	9	02-911-474	8	02-911-827	3
01-829C	2	02-707-149	23	02-896-2A	9	02-911-475	8	02-911-828	8
01-829D	2	02-707-181	23	02-896-2B	9	02-911-476	8	02-911-829	8
01-829E	2	02-707-182	23	02-896-2C	9	02-911-721	3	02-911-830	8
01-829F	2	02-707-400	24	02-896-2D	9	02-911-723	5	02-911-831	8
01-830A	2	02-707-401	24	02-896-2E	9	02-911-724	5	02-911-832	8
01-830B	2	02-707-402	24	02-896-2F	9	02-911-725	5	02-911-833	5
01-830C	2	02-707-403	24	02-896A	9	02-911-726	5	02-911-852	6
01-830D	2	02-707-405	24	02-896B	9	02-911-729	5	02-911-854	7
01-830E	2	02-707-406	24	02-896C	9	02-911-730	5	02-911-857	6
01-831A	2	02-707-407	24	02-896D	9	02-911-731	7	02-911-859	3
01-831B	2	02-707-408	24	02-896E	9	02-911-732	5	02-911-860	6
01-832A	1	02-707-409	24	02-896F	9	02-911-734	4	02-911-861	4
02-202-100	40	02-707-410	24	02-911-379	10	02-911-735	4	02-911-862	4
02-202-103	40	02-707-411	24	02-911-380	10	02-911-736	4	02-911-863	8
02-202B	40	02-707-412	24	02-911-381	10	02-911-737	5	02-911-865	7
02-202C	40	02-707-413	24	02-911-382	10	02-911-738	5	02-911-866	7
02-591-44	19	02-707-414	24	02-911-383	10	02-911-739	5	02-911-867	3
02-591-45	19	02-707-415	24	02-911-384	10	02-911-740	5	02-911-868	3
02-591-46	19	02-707-416	24	02-911-385	10	02-911-741	5	02-911-869	7
02-617-158	18	02-707-417	24	02-911-386	10	02-911-742	5	02-911-871	7
02-617-159	18	02-707-418	24	02-911-387	10	02-911-743	4	02-911-872	6
02-617-160	18	02-707-419	24	02-911-388	10	02-911-744	4	02-911-875	10
02-617-161	18	02-707-420	24	02-911-389	10	02-911-745	5	02-911-876	3
02-617-162	18	02-707-421	24	02-911-390	10	02-911-746	5	02-911-877	10
02-617-163	18	02-707-422	24	02-911-391	10	02-911-747	5	02-911-878	4
02-681-142	25	02-707-423	24	02-911-392	10	02-911-748	5	02-911-879	5
02-681-145	25	02-707-424	24	02-911-393	10	02-911-749	5	02-911-880	7
02-681-146	25	02-707-425	24	02-911-394	10	02-911-750	5	02-911-881	10
02-681-147	25	02-707-426	24	02-911-395	10	02-911-751	5	02-911-883	3
02-681-150	25	02-707-427	24	02-911-396	10	02-911-752	5	02-911-884	10
02-681-151	25	02-707-428	24	02-911-397	10	02-911-755	5	02-911-885	9
02-681-152	25	02-707-429	25	02-911-398	10	02-911-756	5	02-911-887	10
02-681-153	25	02-707-433	24	02-911-399	10	02-911-757	5	02-911-888	7
02-681-154	25	02-707-434	24	02-911-400	10	02-911-758	5	02-911-889	5

Cat. No.	Page
02-911-890	9
02-911-891	7
02-911-892	5
02-911-893	9
02-911-894	6
02-911-895	4
02-911-896	6
02-911-897	6
02-911-898	4
02-911-899	4
02-911-900	4
02-911-901	6
02-911-902	5
02-911-903	5
02-911-904	5
02-911-905	5
02-911-906	8
02-911-907	8
02-911-908	5
02-911-910	6
02-911-911	8
02-911-912	3
02-911-913	7
02-911-914	5
02-911-916	7
02-911-917	7
02-911-918	7
02-911-919	7
02-911-920	3
02-911-921	6
02-911-922	6
02-911-923	8
02-911-924	8
02-911-925	8
02-911-926	8
02-911-927	6
02-911-928	5
02-911-929	7
02-911-930	7
02-911-931	8
02-911-932	7
02-911-933	7
02-911-934	7
02-911-935	7
02-911-936	7
02-911-937	7
02-911-938	7
02-911-939	6
02-911-940	6
02-911-941	8
02-911-942	7
02-911-943	7
02-911-944	6
02-911-945	6
02-911-946	10
02-911-947	10
02-911-964	3
02-911-969	6
02-911-970	6
02-911-971	6
02-911-972	6
02-911-973	6
02-911-993	10
02-911-994	10
02-911-995	10
02-911-996	3
02-911-997	3
02-911-998	6
02-911-999	6
02-912-000	6
02-912-001	6
02-912-002	6
02-912-003	7
02-912-004	7
02-912-005	7
02-912-006	7
02-912-007	7
02-912-008	7
02-912-009	6
02-912-010	6
02-912-011	6
02-912-012	6

Cat. No.	Page
02-912-013	6
02-912-014	6
02-912-016	10
02-912-017	10
02-912-039	7
02-912-040	7
02-912-045	3
02-912-114	7
02-912-116	3
02-912-117	7
02-912-118	3
02-912-119	8
02-912-255	8
02-912-256	8
02-912-257	8
02-912-270	4
02-912-271	4
02-912-272	4
02-912-273	4
02-912-274	4
02-912-275	4
02-912-276	4
02-912-277	4
02-912-278	4
02-912-279	4
02-912-280	4
02-912-296	10
02-912-297	10
02-912-298	3
02-912-299	7
02-912-300	10
02-912-301	10
02-912-302	10
02-912-303	8
02-912-304	3
02-912-305	3
02-912-306	3
02-912-307	4
02-912-308	10
02-912-309	3
02-912-310	3
02-912-311	3
02-912-312	3
02-912-313	3
02-912-314	7
02-912-315	7
02-912-316	7
02-912-317	8
02-912-318	6
02-912-337	5
02-912-338	5
02-912-339	4
02-912-340	5
02-912-341	5
02-912-342	5
02-912-343	5
02-912-344	5
02-912-345	3
02-912-351	37
02-912-352	37
02-912-353	37
02-912-354	37
02-912-355	37
02-912-356	37
02-912-357	37
02-912-358	37
02-912-359	37
02-912-360	38
02-912-361	38
02-912-362	38
02-912-363	38
02-912-364	38
02-912-365	37
02-912-366	37
02-912-367	37
02-912-368	37
02-912-369	37
02-912-370	37
02-912-371	37
02-912-372	37
02-912-373	37
02-912-374	37

Cat. No.	Page
02-912-375	37
02-912-376	37
02-912-377	37
02-912-378	37
02-912-379	37
02-912-380	37
02-912-381	37
02-912-382	37
03-007-30	14
03-007-31	14
03-007-32	14
03-007-33	14
03-007-34	14
03-007-35	14
03-007-36	14
03-007-37	14
03-007-38	14
03-007-39	14
03-007-40	14
03-007-41	14
03-007-42	14
03-007-43	14
03-007-44	14
03-007-45	14
03-007-46	9
03-007-47	9
03-007-48	9
03-007-49	9
03-007-52	9
03-007-53	9
03-007-54	9
03-007-55	9
03-337-1	40
03-337-11B	40
03-337-12B	40
03-337-14	39
03-337-15	39
03-337-2	40
03-337-20	40
03-337-23	40
03-337-23A	40
03-337-23B	40
03-337-23C	40
03-337-26	39
03-337-39	39
03-337-4	39
03-337-40	39
03-337-41	39
03-337-42	39
03-337-46	40
03-337-5	39
03-337-7	39
03-338-1A	39
03-338-1B	39
03-338-1C	39
03-338-1D	39
03-338-1E	39
03-338-1F	39
03-338-1G	39
03-338-3A	39
03-338-3B	39
03-338-3C	39
03-338-3D	39
03-338-3E	39
03-338-3F	39
03-338-3G	39
03-338-3H	39
03-338-3J	39
03-338A	38
03-338AA	38
03-338B	38
03-338C	38
03-338D	38
03-338E	38
03-338F	38
03-338G	38
03-338H	38
03-338J	38
03-338K	38
03-338L	38
03-339-14A	37
03-339-14C	37

Cat. No.	Page
03-339-14E	37
03-339-21A	38
03-339-21B	38
03-339-21C	38
03-339-21D	38
03-339-21E	38
03-339-21F	38
03-339-21G	38
03-339-21H	38
03-339-21J	38
03-339-21K	38
03-339-21M	38
03-339-21N	38
03-339-22A	38
03-339-22B	38
03-339-22C	38
03-339-22D	38
03-339-22E	38
03-339-22F	38
03-339-22G	38
03-339-22H	38
03-339-22J	38
03-339-22K	38
03-339-22M	38
03-339-22N	38
03-339-23A	37
03-339-23B	37
03-339-23C	37
03-339-23D	37
03-339-23E	37
03-339-25A	38
03-339-25B	38
03-339-25C	38
03-339-25D	38
03-339-25E	38
03-339-25F	38
03-339-26A	36
03-339-26B	36
03-339-26C	36
03-339-26D	36
03-339-26F	36
03-339-26G	36
03-339-26H	36
03-339-27A	39
03-339-27B	39
03-339-27C	39
03-339-27D	39
03-339-27E	39
03-339-27F	39
03-339-30A	36
03-339-30B	36
03-339-30C	36
03-339-30E	36
03-339-30F	36
03-339-30G	36
03-339-30H	36
03-339-30J	36
03-340-119	37
03-340-120	37
03-340-121	37
03-340-122	37
03-340-4E	39
03-340-7L	12
03-341-11	40
03-341-75A	39
03-341-75C	39
03-341-75E	39
03-341-75G	39
03-341-76A	39
03-391-100	32
03-391-101	31
03-391-102	31
03-391-103	31
03-391-104	31
03-391-105	31
03-391-106	32
03-391-107	32
03-391-108	32
03-391-109	31
03-391-110	32
03-391-111	31
03-391-112	32

Cat. No.	Page
03-391-113	32
03-391-114	32
03-391-115	32
03-391-116	32
03-391-117	32
03-391-118	32
03-391-119	32
03-391-187	31
03-391-188	31
03-391-189	31
03-391-200	11
03-391-201	11
03-391-202	11
03-391-203	11
03-391-204	11
03-391-205	11
03-395-100	11
03-395-101	11
03-395-102	11
03-395-103	11
03-395-104	11
03-395-105	11
03-395-106	11
03-395-107	11
03-395-108	11
03-395-109	11
03-395-110	11
03-395-111	11
03-395-112	11
03-395-113	11
03-395-114	11
03-395-400	10
03-395-401	11
03-395-402	10
03-395-403	10
03-395-404	10
03-395-600	29
03-395-601	29
03-395-602	29
03-395-603	29
03-395-604	29
03-395-605	29
03-395-606	29
03-395-607	29
03-395-608	29
03-395-609	29
03-395-610	29
03-395-611	29
03-395-612	29
03-395-613	29
03-395-614	29
03-395-615	29
03-395-616	29
03-395-617	29
03-395-618	29
03-395-619	29
03-395-620	29
03-395-621	29
03-395-622	29
03-395-623	29
03-395-624	29
03-395-625	29
03-395-626	29
03-395-627	29
03-395-628	29
03-395-629	29
03-395-630	29
03-395-631	29
03-395-632	29
03-395-634	29
03-395-638	29
03-396-49	21
03-396-50	21
03-396-51	21
03-396-52	21
03-396-53	21
03-396-54	21
03-396-55	21
03-396-56	22
03-409-23A	4
03-409-23B	4
03-409-23C	4

Cat. No.	Page	Cat. No.	Page	Cat. No.	Page	Cat. No.	Page	Cat. No.	Page
03-409-23D	4	05-541-43	28	08-773-1	12	09-795K	16	10-500-2	18
03-409-23E	4	05-541-44	29	08-773-2	12	09-800	16	10-500-25	36
03-409-23G	4	05-541-45	28	08-773-3	12	09-801A	16	10-500-26	36
03-409-23H	4	05-541-46	28	08-887	17	09-801AA	16	10-500-27	36
03-409-23J	4	05-541-47	28	08-890	17	09-801B	16	10-500-28	29
03-409-23K	4	05-541-48	28	08-895	17	09-801BB	16	10-500-29	29
03-409-23L	4	05-541-49	29	08-905	18	09-801C	16	10-500-3	18
03-409-23M	4	05-714-1	12	08-907	18	09-801D	16	10-500-5	18
03-409-23N	4	05-714-4	12	08-914-5	19	09-801E	16	10-500-6	18
03-409-23P	4	06-443-18	33	08-916-5A	19	09-801F	16	10-500-7	18
03-409-23Q	4	06-443-19	33	08-916-5B	19	09-801G	16	10-500-8	18
03-409-23R	4	06-443-20	33	08-916-5C	19	09-801H	16	10-500-9	18
03-409-23S	4	06-443-20PK	33	08-916-5D	19	09-801J	16	12-544-1	23
03-448-1	22	06-443-21	33	08-917-5	19	09-801K	16	12-544-4	23
03-448-10	22	06-662-44	32	08-920A	20	09-802	16	12-545-89	22
03-448-11	28	06-662-51	32	08-920B	20	09-802-1A	16	12-550-15	22
03-448-12	28	06-664-14	40	08-935	30	09-802-1B	16	12-550-16	22
03-448-13	28	06-665-23	40	08-940	30	09-803-5A	16	12-550-16A	22
03-448-14	28	06-665-24	40	08-951-10	29	09-803-5B	16	12-550-17	22
03-448-15	28	07-781Y	12	08-951-20	30	09-803-5C	16	12-550-18	22
03-448-16	28	07-782A	12	08-951-25	30	09-803-5D	16	12-550-19	22
03-448-17	28	07-782B	12	08-951-5	29	09-803-5E	16	12-550-20	22
03-448-2	22	07-782C	12	09-719-5	17	09-803-5F	16	12-550-33	22
03-448-20	28	07-782D	12	09-719-7	17	09-803-6A	16	12-550-343	22
03-448-3	22	07-782E	12	09-719A	17	09-803-6B	16	12-550-A3	22
03-448-4	22	07-782F	12	09-719B	17	09-803-6C	16	12-565-500	20
03-448-5	22	07-782G	12	09-719C	17	09-803-6D	16	12-565-501	20
03-448-6	22	07-782H	12	09-719D	17	09-803-6E	16	12-565-502	31
03-448-7	22	07-782J	12	09-719E	17	09-803-6F	16	12-565-503	31
03-448-8	22	07-782K	12	09-719F	17	09-803-6G	16	12-565-504	31
03-448-9	22	07-791A	12	09-719G	17	09-803-6H	16	12-565-505	31
03-565	11	07-791B	12	09-719H	17	09-803-6J	16	12-565-506	20
03-572-5	11	07-791C	12	09-720-3	17	09-803-6K	16	12-565-507	31
03-578	11	07-791D	12	09-720-4	17	09-804	16	12-565-508	31
04-320-5	14	07-791E	12	09-720-5	17	09-804-110A	16	12-565-509	31
04-342	15	07-791F	12	09-720-6	17	09-804-150A	16	12-565-511	20
04-343	15	07-791G	12	09-720-7	17	09-804-21A	15	12-565-633	11
05-408-120	33	07-791H	12	09-720-8	17	09-804-21B	15	12-565-634	11
05-408-121	33	07-791J	12	09-730-16	17	09-804-21C	15	12-565-638	11
05-408-122	33	07-791K	12	09-730-17	17	09-804-21D	15	12-565-639	11
05-408-123	33	07-791L	12	09-730-19	17	09-804-24A	15	12-565-640	11
05-408-124	33	07-791M	12	09-730-21	17	09-804-24C	15	12-566-120	31
05-408-125	33	07-905-2	13	09-730-25	17	09-804-25D	15	12-566-121	31
05-408-126	33	07-905-6	13	09-730-27	17	09-804-42A	15	12-567-500	36
05-408-127	33	08-100-10	12	09-790-12C	16	09-804-42B	15	12-567-501	36
05-408-128	33	08-100-11	12	09-790-12D	16	09-804-42C	15	12-567-502	36
05-408-129	33	08-100-12	12	09-790-12E	16	09-804-42D	15	13-070-2	19
05-408-130	33	08-100-13	12	09-790-12F	16	09-804-55A	15	13-247-100	34
05-408-131	33	08-408-219	21	09-790-12G	17	09-804-55B	15	13-247-150	34
05-408-132	33	08-555A	14	09-790-14D	16	09-804-55C	15	13-247-300	34
05-408-133	33	08-555B	14	09-790-14E	16	09-804-55D	15	13-247-70	34
05-408-134	33	08-555C	14	09-790-14F	16	09-804-70A	16	13-379	20
05-408-135	33	08-555D	14	09-790-14G	16	09-804-70C	16	13-379-1	20
05-408-136	33	08-555E	14	09-790-2A	17	09-804-70D	16	13-379-2	20
05-408-137	33	08-555F	14	09-790-2B	17	09-804-90A	16	13-379-3	20
05-408-138	33	08-555G	14	09-790-2C	17	09-804-90C	16	13-379-4	20
05-408-139	33	08-557-1A	14	09-790-2D	17	09-898-12A	40	13-379-5	20
05-408-140	33	08-557-1B	14	09-790-2E	17	09-898-12B	40	13-379-6	20
05-408-141	33	08-557-1C	14	09-790-2F	17	09-898-12C	40	13-379-7	20
05-408-145	33	08-557-1D	14	09-790-2G	17	10-270	18	13-650-2A	27
05-408-146	33	08-557-1E	14	09-790-4C	17	10-275	17	13-650-2B	27
05-408-150	24	08-557-1F	14	09-790-4D	17	10-280	18	13-650-2C	27
05-408-151	24	08-557-1G	14	09-790-4E	17	10-295	18	13-650-2D	27
05-408-152	24	08-557-1H	14	09-790-4F	17	10-300	17	13-650-2E	27
05-408-153	24	08-732-112	40	09-790A	17	10-316A	17	13-650-2F	27
05-408-154	24	08-732-113	40	09-790B	17	10-316B	17	13-650-2G	27
05-539-1	33	08-732-114	40	09-790C	17	10-316C	17	13-650-2H	27
05-539-10	33	08-732-115	40	09-790D	17	10-353	15	13-650-2J	27
05-539-12	33	08-732-116	40	09-790E	17	10-354	15	13-650-2K	27
05-539-13	33	08-732-117	40	09-790F	17	10-355	15	13-650-2L	27
05-539-5	33	08-732-118	40	09-790G	17	10-356	15	13-650-2M	27
05-539-6	33	08-757-11C	23	09-795A	16	10-358	15	13-650-2N	27
05-539-7	33	08-757-12	23	09-795AA	16	10-500	18	13-650-2P	27
05-539-8	33	08-757-13	23	09-795B	16	10-500-1	18	13-650-2Q	27
05-539-9	33	08-757-14	23	09-795BB	16	10-500-10	18	13-650-2R	27
05-541-29	28	08-757-150	23	09-795C	16	10-500-11	18	13-650-2S	27
05-541-36	28	08-757-90	23	09-795D	16	10-500-12	18	13-650-2T	27
05-541-37	28	08-757-91	23	09-795E	16	10-500-13	18	13-650-2U	27
05-541-38	29	08-757-9A	23	09-795F	16	10-500-14	18	13-650-2W	27
05-541-39	28	08-757-9B	23	09-795G	16	10-500-15	18	13-650-2X	27
05-541-40	28	08-758-2	23	09-795H	16	10-500-16	18	13-650-2Y	27
05-541-41	28	08-758-3	23	09-795J	16	10-500-17	18	13-650-3A	27

Cat. No.	Page
13-650-3B	27
13-650-3C	27
13-650-3D	27
13-650-3E	28
13-650-3F	28
13-650-3G	28
13-650-3H	28
13-650-3J	28
13-650-3K	28
13-650-3L	28
13-650-3M	28
13-650-3N	28
13-650-3P	28
13-650-3Q	28
13-650-3R	28
13-650-3S	28
13-650-3T	28
13-650-3V	28
13-650-3W	28
13-650-3X	28
13-650-3Y	28
13-665A	27
13-665F	27
13-665G	27
13-665H	27
13-665K	27
13-665M	27
13-665N	27
13-671-108A	27
13-671-108B	27
13-671-108C	27
13-671-108D	27
13-671-108E	27
13-674-26H	27
13-674-26M	27
13-675-2C	26
13-675-3C	26
13-675D	27
13-675F	27
13-675G	27
13-676-10B	26
13-676-10C	26
13-676-10F	26
13-676-10G	26
13-676-10H	26
13-676-10J	26
13-676-10K	26
13-676-10M	26
13-676-10Q	26
13-676-10R	26
13-678-11	26
13-678-11B	26
13-678-11C	26
13-678-11D	26
13-678-11E	26
13-678-11F	26
13-678-12B	27
13-678-12C	27
13-678-12D	27
13-678-12E	27
13-678-14B	27
13-678-14C	27
13-678-20A	26
13-678-20B	26
13-678-20C	26
13-678-20D	26
13-678-25A	26
13-678-25B	26
13-678-25C	26
13-678-25D	26
13-678-25E	26
13-678-26A	26
13-678-27B	26
13-678-27C	26
13-678-27D	26
13-678-27E	26
13-678-27F	26
13-678-31E	26
13-678-31F	26
13-678-31H	26
13-678-31J	26
13-678-36A	26

Cat. No.	Page
13-678-36B	26
13-678-36C	26
13-678-36D	26
13-678-36E	26
13-681-24	28
13-681-50	28
13-681-500	28
13-681-501	28
13-681-502	28
13-681-503	28
13-681-504	28
13-681-51	28
13-681-800	28
13-681-801	28
13-681-802	28
13-683-61A	23
13-683-61B	23
13-683-61C	23
13-683-61D	23
13-683-61E	23
13-683-61F	23
13-683-61G	23
13-683-61H	23
13-683-62A	23
13-683-62C	23
13-683-62D	23
13-683-62E	23
13-683-62F	23
13-683-62G	23
13-711-7M	27
13-711-9D	27
13-804-12	30
13-804-2	29
13-804-4	29
13-804-5	29
13-804-6	30
13-804-8	30
13-806	30
13-806-2	30
13-808-2	30
13-808-4	30
13-810-2	30
13-810-4	30
13-812-10	18
13-812-14	18
13-812-16	18
13-812-18	18
13-812-24	18
13-812-26	18
13-812-27	18
13-812-4	18
13-812-45	18
13-812-8	18
14-127-46	2
14-127-47	2
14-127-48	3
14-127-49	3
14-127-50	3
14-127-52	3
14-127-53	3
14-135A	13
14-135B	13
14-135C	13
14-135D	13
14-135E	13
14-135F	13
14-135G	13
14-135H	13
14-135J	13
14-135K	13
14-135L	13
14-135M	13
14-135N	13
14-135P	13
14-135Q	13
14-135R	13
14-135S	13
14-135T	13
14-135U	13
14-140A	13
14-140B	13
14-140C	13

Cat. No.	Page
14-140D	13
14-140E	13
14-140F	13
14-140G	13
14-140H	13
14-140J	13
14-140K	13
14-140L	13
14-140M	13
14-140N	13
14-140P	13
14-140Q	13
14-140R	13
14-140S	13
14-141A	13
14-141B	13
14-141C	13
14-141D	13
14-141E	13
14-141F	13
14-141G	13
14-141H	13
14-141J	13
14-141K	13
14-141L	13
14-141M	13
14-141N	13
14-141P	13
14-141Q	13
14-141R	13
14-141S	13
14-141T	13
14-141V	13
14-141W	13
14-141X	13
14-169-7A	36
14-169-7B	36
14-169-7C	36
14-169-7D	36
14-169-7E	36
14-169-7F	36
14-169-7G	36
14-169-7H	36
14-169-7J	36
14-169-7K	36
14-169-7L	36
14-169-7M	36
14-169-7N	36
14-178-2AA	36
14-178-2B	36
14-178-2BB	36
14-178-2C	36
14-178-2CC	36
14-178-2D	36
14-206-62	2
14-206-63	2
14-206-64	2
14-206-65	2
14-230-200	21
14-230-201	21
14-230-202	21
14-230-203	21
14-230-204	21
14-230-205	21
14-230-206	21
14-230-207	21
14-230-208	21
14-230-209	21
14-230-210	21
14-230-211	21
14-230-212	21
14-230-213	21
14-230-214	21
14-230-215	21
14-230-216	21
14-230-217	21
14-230-218	21
14-230-219	21
14-230-220	21
14-230-221	21
14-230-222	21
14-230-223	21

Cat. No.	Page
14-230-224	21
14-230-225	21
14-230-226	21
14-230-227	21
14-230-228	21
14-230-229	21
14-230-230	21
14-230-231	21
14-230-232	20
14-230-233	20
14-230-234	20
14-230-235	20
14-230-236	20
14-230-237	20
14-230-238	20
14-230-239	20
14-230-240	20
14-230-241	20
14-230-242	21
14-230-243	21
14-230-244	20
14-230-245	21
14-230-246	21
14-230-247	21
14-230-248	21
14-357Q	30
14-373	30
14-373-25A	30
14-373-25B	30
14-374	30
14-375-10	30
14-375-258	29
14-375-259	29
14-375-260	29
14-512-120	31
14-512-121	31
14-512-122	30
14-512-123	31
14-512-129	31
14-512-145	30
14-512-1451	30
14-512-1452	30
14-512-146	30
14-512-1461	30
14-512-1462	30
14-512-147	30
14-512-1471	30
14-512-1472	30
14-512-148	30
14-512-1481	30
14-512-1482	30
14-512-149	30
14-512-1491	30
14-512-1492	30
14-512-150	30
14-512-151	30
14-512-152	30
14-512-153	30
14-512-154	30
14-512-155	30
14-512-156	30
14-512-158	30
14-512-159	30
14-512-160	30
14-513-51	31
14-513-51SIX	31
14-513-52	31
14-513-53	30
14-513-54	31
14-513-55	31
14-513-56	31
14-513-57	31
14-513-57SIX	31
14-513-58	31
14-513-58SIX	31
14-513-59	31
14-513-59SIX	31
14-513-60	31
14-513-60SIX	31
14-513-61SIX	31
14-513-62	31
14-513-63SIX	31

Cat. No.	Page
14-513-64SIX	31
14-513-65	31
14-513-65SIX	31
14-513-66	31
14-513-66SIX	31
14-513-67	31
14-513-68	31
14-649-17	32
14-826-340	10
14-905-30	20
14-905-35	20
14-956-1B	35
14-956-1C	35
14-956-1D	35
14-956-1G	35
14-956-1H	35
14-956-1J	35
14-956-3A	35
14-956-3B	35
14-956-3C	35
14-956-3D	35
14-956-4A	35
14-956-5A	35
14-956-5B	35
14-956-5C	35
14-956-6A	35
14-956-6B	35
14-956-6C	35
14-956-6D	35
14-956-7A	35
14-956-7C	35
14-956-7E	35
14-956-8E	35
14-956-8G	35
14-956-8J	35
14-957-10A	35
14-957-10B	35
14-957-10C	35
14-957-10D	35
14-957-10E	35
14-957-12A	35
14-957-12B	35
14-957-12C	35
14-957-12D	35
14-957-12E	35
14-957-14A	35
14-957-14B	35
14-957-14C	35
14-957-14D	35
14-957-14E	35
14-957-16A	35
14-957-16B	35
14-957-16C	35
14-957-16D	35
14-957-16E	35
14-957-18D	35
14-957-18E	35
14-957-20A	34
14-957-20B	35
14-957-20C	34
14-957-20D	34
14-957-20E	34
14-957-22A	34
14-957-22B	35
14-957-22C	34
14-957-22D	34
14-957-22E	34
14-957-24A	35
14-957-24B	35
14-957-24C	34
14-957-24D	34
14-957-24E	34
14-957-26A	35
14-957-26B	35
14-957-26C	34
14-957-26D	34
14-957-26E	34
14-957-76A	34
14-957-76C	34
14-957-76E	34
14-957-82A	12
14-957-82E	12



Cat. No.	Page	Cat. No.	Page	Cat. No.	Page
14-958-10A	34	15-904	31	22-305-654	15
14-958-10B	34	15-905	31	22-363-250	11
14-958-10C	34	18-960A	1	22-363-251	11
14-958A	34	18-960B	1	22-363-252	12
14-958B	34	18-960C	1	22-363-253	11
14-958C	34	19-041-189A	18	22-363-254	12
14-958D	34	19-041-189B	18	22-363-255	11
14-958E	34	19-041-189C	18	22-363-547	12
14-958F	34	19-041-189D	18	22-363-548	12
14-958G	34	19-041-189E	18	22-363-549	12
14-958H	34	19-041-190A	19	22-363-595	19
14-958J	34	19-041-190B	19	22-363-596	19
14-958K	34	19-041-190C	19	22-363-597	19
14-959-10AA	35	19-041-190D	19	22-363-598	19
14-959-16	35	19-041-190E	19	22-363-599	19
14-959-16A	35	19-050-548A	19	22-363-600	19
14-959-35A	34	19-050-548B	19	22-363-601	19
14-959-35AA	34	19-050-548C	19	22-363-608	19
14-959-35B	34	19-050-548D	19	23-022-524	1
14-959-35C	34	19-050-550A	19	FB101	14
14-959-36A	12	19-050-550B	19	FB102	14
14-959-36B	12	19-050-550C	19	FB104	14
14-959-37A	34	19-050-550D	19	S30688A	2
14-959-37C	34	19-098-910	1	S50712	33
14-959-96A	1	19-098-912	1	SB105	11
14-959-96C	1	19-098-913	1		
14-961-10	35	19-098-914	1		
14-961-10A	35	19-098-915	1		
14-961-11	35	19-098-916	1		
14-961-13	35	19-098-917	1		
14-961-20	35	19-098-918	1		
14-961-25	33	19-098-921	1		
14-961-26	33	19-098-922	1		
14-961-27	33	19-098-923	1		
14-961-28	33	19-098-924	1		
14-961-29	33	19-098-925	1		
14-961-30	34	19-098-926	1		
14-961-31	34	19-130-1597B	19		
14-961-32	34	19-130-1597C	19		
14-961-33	34	19-130-1597D	19		
14-961-34	34	19-130-1597F	19		
14-962-26D	34	19-130-2088	15		
14-962-26F	34	19-130-2089	15		
14-962-26G	34	19-181-503	15		
14-962-26H	34	19-181-504	15		
14-962-26K	34	19-181-516	20		
14-962-26M	34	19-181-517	20		
15-041-11A	32	19-181-518	20		
15-041-11B	32	19-181-519	20		
15-041-11C	32	19-181-520	20		
15-041-11D	32	21-197-2E	24		
15-041-11E	32	21-197-2F	24		
15-041-11F	32	21-277-2A	23		
15-041-11G	32	21-277-2B	23		
15-041-11H	32	21-377-203	20		
15-041-11J	32	21-377-204	20		
15-041-11K	32	21-377-205	20		
15-041-11L	32	22-031-145	22		
15-041-11M	32	22-034-985	22		
15-041-11N	32	22-037-160	22		
15-041-11P	32	22-037-161	22		
15-041-11S	32	22-037-162	22		
15-041-11V	32	22-037-163	22		
15-041-11W	32	22-037-169	22		
15-041-11Z	32	22-037-170	22		
15-041-12E	32	22-037-171	22		
15-041-12F	32	22-044-560	2		
15-041-6A	32	22-230-482	26		
15-041-6B	32	22-230-890	22		
15-041-6C	32	22-230-891	22		
15-041-6D	32	22-230-892	22		
15-041-7A	32	22-265-446	22		
15-041-7B	32	22-266-858	22		
15-041-7C	32	22-266-882	22		
15-041-7D	32	22-266-908	22		
15-041-8A	32	22-266-940	22		
15-041-8B	32	22-266-973	22		
15-041-8C	32	22-267-005	22		
15-041-8D	32	22-267-013	22		
15-041-9A	32	22-267-062	22		
15-041-9B	32	22-267-096	22		
15-903	31	22-293-232	22		

# What's on Your Shelf?

## Fisher Scientific Lab Essentials Catalog

Your single source for high-quality Fisherbrand consumables and Fisher Scientific lab equipment. The Fisher Scientific Lab Essentials Catalog makes it easier than ever to choose top-quality lab products, equipment, instruments, and more!

- More than 7500 Fisher Scientific products you use every day
- Completely revised, easy-to-read product descriptions
- Intuitive navigational features



# BLUE CREW

WHERE THE SCIENTIFIC COMMUNITY CONNECTS!

## Get Connected!

*The Fisher Scientific Blue Crew is just another exciting way to make BLUE work for you.*

## Registration is fast, easy and free!

### Members Get More

- Official Blue Crew welcome package\*
- Members-only promotions and online offers
- Members-only access to our game archive
- Exclusive participation in focus groups, surveys, and new product testing
- Free wall posters, downloadable screen savers and wallpapers
- More chances to win free stuff

## Sign up online today – You're in good company!

[www.fishersci.com/BlueCrew/Register](http://www.fishersci.com/BlueCrew/Register)

*\*While supplies last*



HOW BLUE IS YOUR LAB? • HOW BLUE IS YOUR LAB? • HOW BLUE IS YOUR LAB?

**Fisherbrand**<sup>®</sup>