

Centrifuges

Fisherbrand micro

Centrifuge, micro, Fisherbrand®



For PCR*, microfilter separation, clinical applications and HPLC protocols.

- For 1.5mL/2mL and 0.5mL and 0.2mL tubes
- Variable speed setting 500rpm to 12,500rpm, in 100rpm increments
- Variable timer 1s to 30min, in 1s increments
- Selectable digital display of speed in RCF or rpm
- Low noise, vibration and a small footprint
- Imbalance and tilt safety cut out
- Quick spin button
- Supplied complete with two rotors and a universal power adapter for worldwide voltages

Compact centrifuge, ideal for most protocols requiring fast spins and quick spin downs.

Operation is extremely simple. For duplicate spins, simply reload rotor and press start. Designed with a Quick-Spin option for commonly used 20s spins. Includes a tool free, quick change rotor system. Safety features ensure the unit will not operate with the lid open, a jammed rotor, or when out of balance. With a 12 place standard rotor and lid for 1.5mL/2mL tubes (lid does not have to be in place to use this rotor), an 8 place strip rotor for 4 x 0.2mL strips or 32 single tubes, 12 x 0.2mL and 12 x 0.5mL adapters for customised use of the standard rotor.

Technical Specification

Speed range, rpm	500 to 12,000
RCF, x g	15 to 9,100
Dimensions, internal [w x d x h], mm	203 x 171 x 114
Electrical supply	100-240V, 50/60Hz

Catalogue No	Legacy No	Description
12972041	-	High speed mini centrifuge

Accessories and spare parts

Catalogue No	Legacy No	Description
12982041	-	Replacement power adapter
12992041	-	Replacement universal cord set
12902051	-	Replacement standard rotor with cover
12912051	-	Replacement strip tube rotor with knob
12922051	-	Replacement rotor knob
12932051	-	Standard rotor cover

*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffmann-La Roche

Have you thought about.....

Trays, utility, Fisherbrand®



- Lightweight solution for transporting tubes and smaller items efficiently



Catalogue No	Legacy No	Description
11740634	-	Utility tray with 13mm OneRack™
11750634	-	Utility tray with 16mm OneRack™

Accessories

Catalogue No	Legacy No	Description	Pack qty
11760634	-	Disposable inserts for 13mm OneRack™	4
11770634	-	Disposable inserts for 16mm OneRack™	4

Look for the Fisherbrand® logo throughout the catalogue. It's your guide to great value products.