

Be Safe



Safe Hands



*Introducing*

# LabServ Nitrile Gloves

The LabServ range of Nitrile Powder free gloves have been specifically designed for laboratory, pharmacy, food processing, food preparation and general industry uses.

#### Product Features and Benefits:

- Made of soft 100% nitrile stronger than latex.
- More puncture resistant than latex, Rubber & vinyl
- Ambidexterous
- Light weight and disposable
- Protection against short term chemical contact
- Odorless & easy donning
- Secure grip – in wet and dry conditions
- Ultra low protein content reduced risk of reactions

#### Best Suited for

Medical Personnel | Lab Technicians | Automotive Workers  
Food Handlers | Product Inspectors



## Product Specification

Product Name	Non-sterile Powder Free Nitrile Examination Gloves		
Thickness	3 MIL		
Average Weight of Medium Size	4.0gm		
Usage	Nitrile (Carboxylated Butadiene Acrylonitrile)		
Cuff Finishing	Beaded		
Color	Cobalt Blue		
Type	Ambidextrous		
Surface finish	Finger Textured		
Product Characteristic	Size	Palm Width (mm)	Length (mm)
Dimension	X-Small	80 ± 5mm	240min
	Small	85 ± 5mm	240min
	Medium	95 ± 5mm	240min
	Large	106 ± 5mm	240min
	Extra Large	115 ± 5mm	240min
Thickness	Single wall thickness (mm)		
	Cuff	0.05 - 0.08	
	Palm	0.07 - 0.10	
	Finger	0.09 - 0.13	
Physical Properties (ASTM Standard)		Un-aged	Aged
	Tensile Strength (MPa)	16-24	18-26
	Ultimate Elongation (%)	550-650	500-600
Physical Properties (EN Standard)		Un-aged	Aged
	Force at break (N)	Min 6	Min 6
Residual Powder	2.0 mg/glove maximum		

## Ordering Information

Product Code	Product Description	Pack
LSG781637	LabServ Nonsterile Powderfree Nitrile Exam <b>Gloves-L</b>	Pack of 1000 Pcs
LSG781635	LabServ Nonsterile Powderfree Nitrile Exam <b>Gloves-M</b>	Pack of 1000 Pcs
LSG781635	LabServ Nonsterile Powderfree Nitrile Exam <b>Gloves-S</b>	Pack of 1000 Pcs



Part of Thermo Fisher Scientific

For customer service, call India toll-free 1-800-209-7001.  
 To fax an order, use 022 6680 3001 or 3002.  
 Email: [qfc.customer@thermofisher.com](mailto:qfc.customer@thermofisher.com) • web: [www.fishersci.in](http://www.fishersci.in)  
 © 2012 Thermo Fisher Scientific Inc. All rights reserved